

Jun 22, 2018 Version 3

Adult mouse thymus dissociation (on ice) V.3

DOI

dx.doi.org/10.17504/protocols.io.q7edzje



Andrew Potter¹

¹Cincinnati Children's Hospital Medical Center

Human Cell Atlas Metho...



Andrew Potter

CCHMC

OPEN  ACCESS



DOI: dx.doi.org/10.17504/protocols.io.q7edzje

Protocol Citation: Andrew Potter 2018. Adult mouse thymus dissociation (on ice). **protocols.io**

<https://dx.doi.org/10.17504/protocols.io.q7edzje>

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: June 22, 2018

Last Modified: June 22, 2018

Protocol Integer ID: 13254

Keywords: thymus, CAP, single cell, dissociation

Abstract

Protocol for thymus dissociation (10-week old CD-1 female).

Attachments



thymus_final_1_31_18...

50KB

Guidelines

Storage Conditions of Reagents

Reagent	Storage Condition
DPBS (Thermofisher, 14190144)	4°C
0.5 M EDTA (Ambion, AM9260G)	room temp.
BSA (Sigma, A8806)	4°C
Protease from <i>Bacillus Licheniformis</i> (Sigma, P5380)	Store 100 µL aliquots (100 mg/mL) in DPBS at -80°C
DNase 1 (Applichem, A3778)	Store 10 µL aliquots (250 U/10 µL) in DPBS at -80°C

***Bacillus Licheniformis* enzyme mix (1 mg/mL enzyme):**

492 µL DPBS (No added Ca, Mg)
 0.5 mM EDTA (0.5 µL of 0.5 M EDTA/mL)
 125 U DNase 1 / mL (2.5 µL)
 5 µL of 100 mg/mL enzyme (final conc. 1 mg/mL)


***Bacillus Licheniformis* enzyme mix (2 mg/mL enzyme):**

487 µL DPBS (No added Ca, Mg)
 0.5 mM EDTA (0.5 µL of 0.5 M EDTA/mL)
 125 U DNase 1 / mL (2.5 µL)
 10 µL of 100 mg/mL enzyme (final conc. 2 mg/mL)








+12.5 mg of tissue

Materials

MATERIALS

 Please see Guidelines for required materials.



- 1 Quickly isolate thymus and immerse in ice-cold PBS.
- 2 Place thymus on petri dish on ice using sterile forceps.
- 3 Remove red regions rich in red blood cells using razorblade.
- 4 Using razor blade, mince whole thymus on petri dish, on ice 2 min until fine paste.
 00:02:00 Mincing
- 5 Weigh out 12.5 mg tissue on petri dish.
- 6 Using razor blade, place tissue in 1.5 mL tube containing 0.5 mL digest mix (1 mg/mL) on ice.
 0.5 mL Digest mix (1 mg/mL)
- 7 Shake every 30 seconds to re-suspend tissue for 2 minutes.
 00:02:00 Re-suspending
 00:00:30 Shake
- 8 At 2 min, triturate gently 10X using 1 mL pipet set to 400 μ L.
- 9 For 3 additional minutes (5 min total time), every minute remove tube and triturate gently 10X using 1 mL pipet set to 400 μ L.
 00:03:00 Digest on ice
 00:01:00 Remove tube and triturate gently 10X
- 10 Let tissue chunks settle for 1 min on ice.
 00:01:00 On ice
- 11 At 6 mins total time, remove 80% (400 μ L) of supernatant consisting of dissociated cells (leaving undissociated tissue chunks at the bottom of the tube) and apply to 30 μ M filter on sterile 50 mL conical- rinse filter with 6 mL ice-cold PBS/BSA 0.04%. Save 50 mL conical and filter for next steps.



🧪 6 mL ice-cold PBS/BSA 0.04%

- 12 Add additional 0.5 mL enzyme mix (2 mg/mL) to residual tissue chunks in 1.5 mL tube.

🧪 0.5 mL Enzyme mix (2 mg/mL)

- 13 For 6 additional min (12 min. total), continue triturating gently (10x) every minute on ice.

⌚ 00:06:00 Digest

⌚ 00:01:00 Triturate gently (10x)

- 14 After 12 min. total digest time, triturate digest mix 10X and transfer to 30 μ M filter (the same tube/filter as used previously).

- 15 Rinse filter with 6 mL ice-cold PBS/BSA 0.04%.

🧪 6 mL ice-cold PBS/BSA 0.04%

- 16 Transfer flow-through to 15 mL conical and spin down 650 G for 5 minutes at 4° C.

🌡️ 4 °C Spin down

⌚ 00:05:00 Spin down

- 17 Remove supernatant and re-suspend in 1 mL total volume PBS/BSA 0.04% in a 1.5 mL tube.

🧪 1 mL PBS/BSA 0.04%

- 18 Spin 610 G for 5 minutes at 4 °C.

⌚ 00:05:00 Spinning

🌡️ 4 °C

- 19 Remove supernatant and re-suspend in 1 mL ice-cold PBS/BSA 0.04%.

🧪 1 mL ice-cold PBS/BSA 0.04%

- 20 Determine cell yield and viability using hemocytometer with trypan blue. Adjust concentration to 1,000 cells / μ L for 10x Chromium or 100 cells / μ L for DropSeq.