

Sep 17, 2025

Agarose–gelatin double embedding prior to paraffin processing

DOI

<https://dx.doi.org/10.17504/protocols.io.bp2l6z7e5gqe/v1>

Hyoyoung Maeng¹, Hyuk Song¹

¹School of Advanced Biotechnology, Konkuk University, 120, Neungdong-ro, Gwangjin-gu, Seoul, Republic of Korea, 05029

BioTechniques



songh

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.bp2l6z7e5gqe/v1>

Protocol Citation: Hyoyoung Maeng, Hyuk Song 2025. Agarose–gelatin double embedding prior to paraffin processing. protocols.io <https://dx.doi.org/10.17504/protocols.io.bp2l6z7e5gqe/v1>

License: This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: September 16, 2025

Last Modified: September 17, 2025

Protocol Integer ID: 227372

Keywords: Agarose, Gelatin, Double embedding, Neutral buffered formalin, Paraffin, Spheroid, small tissue fragment, Histology, spheroids for routine paraffin, gelatin, paraffin processing, routine paraffin, agarose, paraffin processing this protocol detail, small tissue fragment, handling of small specimen, spheroid

Funders Acknowledgements:

National Research Foundation of Korea

Grant ID: NRF-2022R1A2C1004503

Abstract

This protocol details an agarose–gelatin double embedding workflow to immobilize small tissue fragments or spheroids for routine paraffin embedding. Samples are embedded in a 2% agarose–5% gelatin mixture, then fixed in neutral buffered formalin (NBF) and processed by standard FFPE methods. Compared with conventional methods, it facilitates handling of small specimens, preserves orientation, and minimizes sample damage.

Attachments



[Agarose–gelatin doub...](#)

179KB

Guidelines

Troubleshooting

- Ensure that the agarose–gelatin solution is well mixed and not solidified before embedding.
- Prevent the sample from drying out.
- Process samples in small batches to avoid solidification of the agarose–gelatin gel.
- If the agarose–gelatin gel does not fully cover the sample, add additional solution to embed completely.

Materials

Reagents

- 
 Agarose **Invitrogen Catalog #75510-019**
- 
 Gelatin, from porcine skin **Merck MilliporeSigma (Sigma-Aldrich) Catalog #G1890**
- 
 DPBS **SolBio Catalog #DPB 001**
- 
 Formalin solution, neutral buffered, 10% **Merck MilliporeSigma (Sigma-Aldrich) Catalog #HT501128**

Materials

- 
 200 µL pipet tip **Corning Catalog #T-200-Y**
- 
 1000 µL pipet tip **Thermo Fisher Scientific Catalog #112NXL-Q**
- 
 50mL Conical Tube **SPL life sciences Catalog #50050**
- 
 Petri Dish **SPL life sciences Catalog #10090**
- 
 10mL Serological Pipet **SPL life sciences Catalog #91010**
- 
 Microtubes **Corning Catalog #MCT-175-C**
- 
 Single Edge Blade **Dorco Catalog #DN-52**
- 
 Kimtech Science Wipers **Yuhan-Kimberly Catalog #41112**

Equipment

Equipment		
AX1502KR	NAME	
ADVENTURER®	BRAND	
-	SKU	

Equipment

MWO-18M1	NAME
SK Magic	BRAND
-	SKU

Equipment

BW-10E	NAME
Jeiotech	BRAND
-	SKU

Safety warnings

- ! -Exercise caution to avoid burns when handling heated gels.
- During fixation, use appropriate equipment to prevent exposure to fixatives.

Before start

- Ensure all reagents and materials are ready before you begin.
- Verify that the samples are properly prepared.
- Confirm the water bath and microwave oven are working properly.



Procedure

- 1 Set the water bath to $60\text{ }^{\circ}\text{C}$.
- 2 Prepare a $4\text{ Mass} / \% \text{ volume}$ agarose solution($200\text{ }\mu\text{L}$ per sample) : add agarose powder to distilled water (DW) in a 50 mL conical tube, heat in a microwave oven until fully dissolved.
- 3 Keep the $4\text{ Mass} / \% \text{ volume}$ agarose solution in the $60\text{ }^{\circ}\text{C}$ water bath.
- 4 Prepare a $10\text{ Mass} / \% \text{ volume}$ gelatin solution: add gelatin powder to a volume of DW equal to 1.2x the volume of the $4\text{ Mass} / \% \text{ volume}$ agarose solution in a 50 mL tube.

Note

A larger volume of gelatin solution is required because some loss may occur during transfer of the viscous solution.

- 5 Place the $10\text{ Mass} / \% \text{ volume}$ gelatin solution in the $60\text{ }^{\circ}\text{C}$ water bath for $00:20:00$ to dissolve completely.
- 6 While maintaining $60\text{ }^{\circ}\text{C}$, mix the $10\text{ Mass} / \% \text{ volume}$ gelatin solution with the $4\text{ Mass} / \% \text{ volume}$ agarose solution at a 1:1 volume ratio using a serological pipette.
- 7 Adjust the water bath temperature to $40\text{ }^{\circ}\text{C}$ to cool the final $2\text{ Mass} / \% \text{ volume}$ agarose– $5\text{ Mass} / \% \text{ volume}$ gelatin solution.

Sample Preparation

- 8 Transfer small tissue fragments or spheroids into Microtubes and wash twice with DPBS.



- 9 Add > 20× sample volume of neutral buffered formalin (NBF) solution to the tube and fix at Room temperature for 00:10:00 - 00:30:00 , depending on sample size.

Note

If the sample is not fixed before double embedding, it may be distorted by water imbibition.

- 10 Wash the sample twice with DW, each for 00:05:00 .

First Embedding

- 11 Place the sample onto a Petri dish.

Note

Transfer the samples with a small amount of water using a pipette with the tip cut off to minimize damage.

- 12 Remove excess water around the sample with tissue paper.

Note

Proceed immediately to the next step before the sample dries.

- 13 Add 200 µL of 2 Mass / % volume agarose- 5 Mass / % volume gelatin solution onto the sample.

Note

Keep the area around the sample free of air bubbles.

- 14 Using pipet tips or forceps, position the sample at the center and adjust it so that the desired cutting direction faces downward.

- 15 Allow the gel to solidify at Room temperature for 00:40:00 .



- 16 Trim the solidified gel into a pyramid shape (~25 mm³ base) using a blade.

Note

The pyramid shape helps maintain the orientation of the sample during subsequent processing.

Citation

McClelland KS, Ng ET, Bowles J (1970)
. Agarose/gelatin immobilisation of tissues or embryo segments for orientated paraffin embedding and sectioning.
Differentiation; research in biological diversity.

<https://doi.org/10.1016/j.diff.2015.12.001>

LINK

- 17 Transfer the pyramid-shaped gel into a Microtube containing  1 mL of NBF solution.

Note

Unfixed gels tend to stick together; therefore, fixation must be performed in separate Microtubes for each sample.

- 18 Fix at  4 °C for  24:00:00 .

- 19 Proceed with routine paraffin embedding according to standard procedures.

Protocol references

[1] McClelland KS, Ng ET, Bowles J. Agarose/gelatin immobilisation of tissues or embryo segments for orientated paraffin embedding and sectioning. Differentiation. 2016;91(4-5):68-71 doi:10.1016/j.diff.2015.12.001.



Citations

Step 16

McClelland KS, Ng ET, Bowles J. Agarose/gelatin immobilisation of tissues or embryo segments for orientated paraffin embedding and sectioning.

<https://doi.org/10.1016/j.diff.2015.12.001>