


Nov 26, 2019

## Aichivirus 3C3D RT-PCR

DOI

[dx.doi.org/10.17504/protocols.io.2i7gchn](https://dx.doi.org/10.17504/protocols.io.2i7gchn)



Judy A Northill<sup>1</sup>

<sup>1</sup>Public Health Virology, Forensic and Scientific Services



Judy A Northill

Public Health Virology, Forensic and Scientific Services

### Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.2i7gchn>

**Protocol Citation:** Judy A Northill 2019. Aichivirus 3C3D RT-PCR. **protocols.io**  
<https://dx.doi.org/10.17504/protocols.io.2i7gchn>

**Manuscript citation:**

**License:** This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

**Protocol status:** Working

**We use this protocol and it's working**

**Created:** May 02, 2019


**Last Modified:** November 26, 2019

**Protocol Integer ID:** 22847

**Keywords:** AiV, aichivirus, aichi, RT-PCR, gastro, faecal, faeces, aichivirus 3c3d rt, pcr this rt, pcr, 3d, junction region of 3c, rt, 3c, human sample,

## Abstract

This RT-PCR will detected Aichivirus A from human samples. It spans the junction region of 3C and 3D and is used for genotyping.



Schematic of Aichivirus A with the 3C3D RT-PCR target region.

## Guidelines

Method assumes the user is familiar with the thermocycler and software used to run the protocol. Mix preparation should occur in a different laboratory or room to the amplification or post-PCR area.

## Materials

### STEP MATERIALS

- ⊗ 5ml Ethidium Bromide Solution [0.625mg/ml] **G-Biosciences Catalog #R041**
- ⊗ Agarose low EEO (Agarose Standard) **AppliChem Catalog #A21140100**
- ⊗ 100bp DNA Ladder, 250ul (50 lanes) **Promega Catalog #G2101**
- ⊗ SuperScript™ III One-Step RT-PCR System with Platinum™ Taq High Fidelity DNA Polymerase **Thermo Fisher Scientific Catalog #12574035**

## Protocol materials

⊗ 5ml Ethidium Bromide Solution [0.625mg/ml] **G-Biosciences Catalog #R041**

⊗ Agarose low EEO (Agarose Standard) **AppliChem Catalog #A21140100**

⊗ 100bp DNA Ladder, 250ul (50 lanes) **Promega Catalog #G2101**

⊗ SuperScript™ III One-Step RT-PCR System with Platinum™ Taq High Fidelity DNA Polymerase **Thermo Fisher Scientific Catalog #12574035**

⊗ SuperScript™ III One-Step RT-PCR System with Platinum™ Taq High Fidelity DNA Polymerase **Thermo Fisher Scientific Catalog #12574035**

⊗ 5ml Ethidium Bromide Solution [0.625mg/ml] **G-Biosciences Catalog #R041**

⊗ Agarose low EEO (Agarose Standard) **AppliChem Catalog #A21140100**

⊗ 100bp DNA Ladder, 250ul (50 lanes) **Promega Catalog #G2101**

## Troubleshooting

## Oligonucleotide sequences

1

| Name      | 5'-3'                  |
|-----------|------------------------|
| AiV-6213F | ACTGGGCCACCCTCCAGACG   |
| AiV-7044R | GGTTGATTTCAGCTTGGAGTTC |

## Reaction set-up

2

- Prepare sufficient for number of reaction plus a 'dead volume', usually 2 extra. Adjust as necessary if using a robotic dispenser.



SuperScript™ III One-Step RT-PCR System with Platinum™ Taq High Fidelity DNA Polymerase **Thermo Fisher Scientific Catalog #12574035**

| Reagent   | Vol (μL)<br>x1 | Final reaction concentration |
|---|----------------|------------------------------|
| Nuclease free water                             | 3.6            |                              |
| Primer AiV-6213F (20pmol/μl)                    | 0.5            | 500nM                        |
| Primer AiV-7044R (20pmol/μl)                    | 0.5            | 500nM                        |
| 2 X Reaction mix                                | 10             | 1X                           |
| Superscript III RT/Platinum Taq HiFi enzyme mix | 0.4            |                              |
| <b>TOTAL VOLUME</b>                             | <b>15</b>      |                              |

Dispense 15μL to each reaction well.

Add 5μL of template, extracted RNA, controls or NTC (nuclease-free water).  
Total reaction volume is 20μL

## Amplification

3 The assay has been used with Eppendorf thermocyclers.

PCR cycling times

|  |                            |                            |                           |      |
|--|----------------------------|----------------------------|---------------------------|------|
|  | 1<br>cycl<br>e             | 40<br>cycl<br>es           | 1<br>cycl<br>e            | Hold |
|  | 50°<br>C 30<br>minu<br>tes | 94°<br>C 15<br>seco<br>nds | 68°<br>C 5<br>minu<br>tes | 15°C |
|  | 94°<br>C 2<br>minu<br>tes  | 50°<br>C 30<br>seco<br>nds |                           |      |
|  |                            | 68°<br>C 60<br>seco<br>nds |                           |      |

4

Amplified products are analysed by gel electrophoresis or equipment such as a QIAxcel.

For gel electrophoresis, a 1.5% agarose gel with ethidium bromide was made using the following recipe.

|  | <b>Reagent</b>   | <b>Volu<br/>me</b> |
|--|------------------|--------------------|
|  | Agarose          | 1.5g               |
|  | 0.5 X TBE buffer | 100<br>ml          |

Boil in microwave for 1.5 -2 minutes until agarose powder is dissolved.

Add 2 drops of ethidium bromide (0.625mg/ml) and mix before pouring warm into the gel form with a comb.

Leave for approx 30minutes at room temperature to set.

#### Note

For the novice user to gel electrophoresis there a detailed protocol at <https://dx.doi.org/10.17504/protocols.io.s38egrw> that will give you further guidance.

Experienced users may run their favorite gel recipe with appropriate dye. The aim is to visualise the amplified products. Gel running times and voltage can be adjusted depending on the equipment used.



⌘ 5ml Ethidium Bromide Solution [0.625mg/ml] **G-Biosciences Catalog #R041**

⌘ Agarose low EEO (Agarose Standard) **AppliChem Catalog #A21140100**


⌘ 100bp DNA Ladder, 250ul (50 lanes) **Promega Catalog #G2101**

Place gel (on a tray) into the tank and cover with 0.5X TBE with ethidium bromide.

Mix 2µl of loading dye with 10µl of amplified product and add to a well in the gel.

Mix 2µl of loading dye with 5µl of 100bp marker and add to a well in the gel. Ideally, the first and last lane but this is not essential.

Run the gel for 60-90 minutes at 80 volts or until the bands on the marker have separated adequately without run off the end of the gel.



Example of 3C-3D products analysed using a QIAxcel.

#### Expected result

Amplified products will produce a band of 832bp in size. No template controls (NTC) should be not detected.