

Sep 17, 2024

Akoya Phenocycler-Fusion Tissue Staining and Imaging Protocol for Formalin-Fixed Paraffin-Embedded Samples adapted from Indiana University

 Forked from [Indiana University adapted Akoya Phenocycler-Fusion Tissue Staining and Imaging Protocol for Fresh Frozen Kidney Samples](#)

DOI

<https://dx.doi.org/10.17504/protocols.io.4r3l229xpl1y/v1>

Archibald Enniful¹, negin.farzad¹, yao.lu.yl¹, Rong Fan¹

¹Yale University



Archibald Enniful

Yale University

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.4r3l229xpl1y/v1>

Protocol Citation: Archibald Enniful, negin.farzad , yao.lu.yl , Rong Fan 2024. Akoya Phenocycler-Fusion Tissue Staining and Imaging Protocol for Formalin-Fixed Paraffin-Embedded Samples adapted from Indiana University. **protocols.io**

<https://dx.doi.org/10.17504/protocols.io.4r3l229xpl1y/v1>



License: This is an open access protocol distributed under the terms of the **[Creative Commons Attribution License](#)**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: September 21, 2023

Last Modified: September 17, 2024

Protocol Integer ID: 88119

Keywords: embedded human lymph node sample, ffpe human lymph node sample, human lymph node sample, imaging protocol for fresh frozen kidney sample, fusion tissue staining, human biomolecular atlas program consortia, cell type marker, cellular senescence network, fresh frozen kidney sample, imaging protocol for formalin, other immune cell type, various cell type, cell type, yale university tmc akoya phenocycler

Funders Acknowledgements:

Yale TMC for Cellular Senescence in Lymphoid Organs

Grant ID: 1U54AG076043-01

Abstract

This protocol presents the Yale University TMC Akoya Phenocycler-Fusion Tissue Staining & Imaging Protocol for Formalin-Fixed Paraffin-Embedded Human Lymph Node Samples. This protocol was adapted from Indiana University adapted Akoya Phenocycler-Fusion Tissue Staining and Imaging Protocol for Fresh Frozen Kidney Samples.

This protocol has been used for FFPE human lymph node samples for the Cellular Senescence Network (SenNet) Program and the Human BioMolecular Atlas Program consortia. The size of the marker panels has ranged between 35 to 47 targets and includes cell type markers for B cells, T cells, macrophages, and other immune cell types. that will label various cell types Multiple sections have been placed on a single slide and imaged simultaneously.

Materials

Akoya Biosciences Phenocycler-Fusion and Staining Kit

Antibodies of choice

Coplin jars

Methanol



Tissue Baking - Day 1

12h

- 1 Bake tissue slides in a 60°C drying oven overnight for optimal tissue adherence to the glass slide. For some tissue samples, baking at 60°C for 1 hour will be sufficient to stick to the glass slide.
 - For optimal use of reagents, 2-3 tissue slides are prepared in one batch.

Tissue deparaffinization and rehydration

1h

- 2 Locate and/or label 10 plastic Coplin jars and fill them with the following reagents:
 - Two Coplin jars filled with Xylene (HistoChoice can be used as an alternative)
 - Two Coplin jars filled with 100% Ethanol
 - One Coplin jar filled with 90% Ethanol
 - One Coplin jar filled with 70% Ethanol
 - One Coplin jar filled with 50% Ethanol
 - One Coplin jar filled with 30% Ethanol
 - Two Coplin jars filled with double-distilled water
- 3 Remove tissue slides from the 60°C drying oven and proceed with the tissue deparaffinization steps. 1m
- 4 Place the tissue slide into the first Coplin jar filled with Xylene for 5 minutes. 5m
- 5 Remove the tissue slide from the first Coplin jar filled with Xylene and place it into the second Coplin jar filled with Xylene for 5 minutes. 5m
- 6 Remove the tissue slide from the second Coplin jar filled with Xylene and place it into the first Coplin jar filled with 100% ethanol for 5 minutes. 5m
- 7 Remove the tissue slide from the first Coplin jar filled with 100% ethanol and place it into the second Coplin jar filled with 100% ethanol for 5 minutes. 5m
- 8 Remove the tissue slide from the second Coplin jar filled with 100% ethanol and place it into the Coplin jar filled with 90% ethanol for 5 minutes. 5m
- 8.1 Remove the tissue slide from the Coplin jar filled with 90% ethanol and place it into the Coplin jar filled with 70% ethanol for 5 minutes. 5m



9 Remove the tissue slide from the Coplin jar filled with 70% ethanol and place it into the Coplin jar filled with 50% ethanol for 5 minutes.

5m

10 Remove the tissue slide from the Coplin jar filled with 50% ethanol and place it into the Coplin jar filled with 30% ethanol for 5 minutes.

5m

11 Remove the tissue slide from the Coplin jar filled with 30% ethanol and place it into the first Coplin jar filled with water for 5 minutes.

5m

12 Remove the tissue slide from the first Coplin jar filled with water and place it into the second Coplin jar filled with water for 5 minutes.

5m

Antigen retrieval and pre-staining steps

2h

13 Akoya Biosciences provides two antigen retrieval buffers: AR6 (Citrate buffer, pH 6) and AR9 (Tris EDTA buffer, pH 9). The choice of antigen retrieval buffer depends on the target antigen. In our experience, AR9 buffer works optimally for most antigens.

A pressure cooker is used for antigen retrieval. Alternatives could be a microwave or a hot plate.

14 Prepare 50 ml 1X AR9 buffer from the stock 10X AR9 buffer. (Dilute with water)

15 Place the slide into the 1X AR9 buffer in and seal the Coplin jar.
Fill the pressure cooker with water to approximately a quarter or a third full.

16 Place the coplin jar into the pressure cooker and set the timer for 20 mins at high pressure.

20m

17 After 20 minutes, turn the pressure valve on the pressure cooker to release pressure. Apply caution since high temperatures could potentially lead to burns

18 Remove the Coplin jar with the slide from the pressure cooker and allow it to cool to room temperature for 30 minutes to 1 hour.

1h

19 While the coplin jar cools to room temperature, locate and/or label 5 plastic coplin jars and fill them with the following reagents:

Two Coplin jars filled with double-distilled water

Two Coplin jars filled with hydration buffer



One Coplin jar filled with staining buffer

- 20 Move the slide into the first coplin jar filled with water and dip the slide 2-3 times to remove the AR9 buffer solution.
- this is not an incubation step, just a quick rinse 1m
- 21 Move the slide into the next coplin jar filled with water and incubate the slide for 2 minutes to thoroughly remove any remaining AR9 buffer solution. 2m
- 22 Move the slide into the first Hydration Buffer coplin jar and incubate the slide for 2 minutes. 2m
- 23 Move the slide into the second Hydration Buffer coplin jar and incubate the slide for another 2 mins. 2m
- 24 Place slide in Staining Buffer and incubate for 20-30 minutes
- timing depends on how quickly the antibody cocktail can be made 30m
- 25 During the Staining Buffer incubation, prepare the Antibody Cocktail

Preparing Antibody Cocktail

- 26 Remove the selected antibodies, spin down if necessary, and place on ice
- 27 Prepare the appropriate volume of Blocking Buffer for the number of slides being stained
*See attached file for Blocking Buffer components
 CODEX_Master_Mix.xlsx 18.9KB
- 28 Label one tube for each **unique** Antibody Cocktail being prepared
- 28.1 The FINAL volume for each tube should be 200 uL - counting antibodies. Remove 1 uL of Blocking Buffer per 1 uL of antibody being added.
- Amounts will vary depending on experiment and desired antibody concentration
- 29 Add the desired amount of antibodies to the appropriate tubes



- Our concentrations are listed in the attached table. Antibody clones and vendors available upon request



CODEX_Master_Mix.xlsx 18.9KB

30 Pipette gently to mix the solution

Tissue Staining

3h

31 Optional step: Cut a rectangular piece of parafilm that is approximately the same size as the sample slide

32 Pre-load a pipette with 190 uL of the prepared Antibody Cocktail

33 Remove slide from Staining Buffer and use a kimwipe to gently dry the slide, without touching the tissue sample(s)

34 Place slide on the humidity chamber, tissue side up, and dispense the Antibody Cocktail
- make sure the tissue sample is completely covered, and there are no air bubbles present

35 Optional step: Gently place parafilm over tissue sample

36 Incubate Sample Slide at room temperature for 3 hours.

3h

FFPE Tissue post-staining

37 Locate and/or label the following plastic coplin jars:
2 x Staining Buffer
1 x Post-Stain Fixative
1 x 100% Methanol
4 x PBS (1x concentration)
1 x CODEX Buffer + Buffer Additive (1x concentration)

38 Fill Methanol coplin jar and place in freezer until needed



Post Staining

51m

- 39 If parafilm steps were included, gently remove parafilm from sample slide
- 40 Place sample slide in Staining Buffer coplin jar and incubate for 2 minutes
 - Dip slide 2-3 times to ensure Antibody Cocktail is fully washed away
- 41 Place sample slide in the second Staining Buffer coplin jar and incubate for 2 minutes
- 41.1 During this incubation, prepare the Post-Staining fixative solution
 - For one coplin jar, use 36 mL of storage buffer and 4 mL of fresh 16% PFA
- 42 Move slide into the Post-Stain coplin jar and incubate for 10 minutes
- 43 Remove slide from the Post-Stain fix and wash in the first PBS coplin jar
 - this is not an incubation step, just a few quick dips to rinse the tissue
- 44 Repeat step 43 in the second and third PBS coplin jars
- 45 With slide still in the third PBS coplin jar, remove Methanol coplin jar from the freezer
- 46 Place slide in the ice cold Methanol for 5 minutes
- 47 After Methanol incubation, immediately place slide in the first PBS coplin jar to rinse
 - As in step 43, this is not an incubation, just a quick rinse
- 48 Repeat step 43 with the second and third PBS coplin jars
- 49 With the slide still in the third PBS coplin jar, prepare the final fixative solution
 - for up to 5 slides, use 1000 uL of 1xPBS and 20 uL of Fixative Reagent
- 50 Preload a pipette with 200 uL of final fixative solution



51 Remove sample slide from the third PBS coplin jar and place in the humidity chamber, tissue side up

52 Dispense 200 uL of Final Fixative solution and incubate for 20 minutes

20m

53 After the 20 minute incubation, rinse the slides in each of the three PBS coplin jars, as was done previously.

1m

54 There are two options for how to proceed next:

a) if imaging is not taking place immediately, sample slide(s) can be stored in Storage Buffer at 4C for up to 5 days.

b) if imaging is taking place the same day staining is finished, proceed to step 45

Cover-slipping the Slide for Imaging

10m 30s

55 Place slide to be cover-slipped in a fresh jar of 1x PBS and incubate for 10 minutes

10m

56 After the 10 minute incubation, carefully dry the slide with a kimwipe, so the coverslip can adhere properly

57 Place coverslip - sticky side up - on the stage. Make sure the coverslip is straight, as overhang will create an inadequate seal or cause the slide to not fit in the flow cell

58 Gently place the slide - tissue side facing the sticky side of the coverslip - on the stage, making sure the slide is lined up with the coverslip

59 Push the stage under the pressure arm, lower the arm, and let rest with for 30 seconds

30s

60 Place the now cover-slipped slide into 1x CODEX Buffer for 10 minutes to allow for equilibration and proper adherence of coverslip to take place.

10m

Preparing the Reporter Plate

61 Prepare the reporter plate stock solution as described in the attached spreadsheet.

62 There will be one unique well per cycle, including the blanks. Each well gets a different volume of stock solution based on the amount of reporter to be added. Final volume of each well should total 250 uL.

In general, the volumes are as follows, per number of reporters being added

3 reporters → 235 uL of stock solution

2 reporters → 240 uL of stock solution

1 reporter → 245 uL of stock solution

0 reporters → 250 uL of stock solution (typically the two blank cycles)

62.1 The stock solution for the first and last cycles go into wells H1 and H2, respectively. Stock solution + reporters will start in A1 and continue in order until all cycles are accounted for.

*An example of a reporter plate set up is shown in the attached spreadsheet

 CODEX_Master_Mix.xlsx 18.9KB

63 Once all the necessary wells are filled, cover wells with foil plate seal to protect from debris and prevent the wells from drying out

*we have prepared reporter plates up to a week in advance, but generally try to prepare them no sooner than the day before an imaging session

Imaging the Slides

64 When reporter plates and slides are ready to use, follow the instructions on the Phenocycler-Fusion to begin the imaging experiment. An example of our experimental settings is shown in the attached spreadsheet.

 Imaging Parameters.xlsx

Analysis Options Available

65 Segmentation and analysis for each sample are conducted using Qupath software (<https://qupath.github.io/>).

Downstream analysis following cell segmentation and phenotyping is done using Seurat in a custom R Studio script that can be made available.