

Sep 08, 2025

Alkaline Lysis Plasmid Purification

DOI

<https://dx.doi.org/10.17504/protocols.io.yxmvmby9og3p/v1>

Agustín Sanchez-Temprano¹

¹CICA - Centro Interdisciplinar de Química e Biología.



Agustín Sanchez-Temprano

CICA - Centro Interdisciplinar de Química e Biología.

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.yxmvmby9og3p/v1>

Protocol Citation: Agustín Sanchez-Temprano 2025. Alkaline Lysis Plasmid Purification . **protocols.io**

<https://dx.doi.org/10.17504/protocols.io.yxmvmby9og3p/v1>

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: September 05, 2025



Last Modified: September 08, 2025

Protocol Integer ID: 226560

Keywords: plasmid dna purification by alkaline lysis, plasmid dna purification, alkaline lysis plasmid purification, plasmid dna, unique structural properties of plasmid dna, quality plasmid suitable for downstream application, quality plasmid, chromosomal dna, dna, alkaline lysis, cloning, such as cloning, molecular biology technique, used molecular biology technique

Abstract

Plasmid DNA purification by alkaline lysis is a widely used molecular biology technique for isolating plasmid DNA from *Escherichia coli* cultures. This method leverages the unique structural properties of plasmid DNA to separate it effectively from chromosomal DNA and cellular debris, yielding high-quality plasmid suitable for downstream applications such as cloning, sequencing, PCR, and transfection.

Guidelines

- Use sterile techniques to avoid contamination that may degrade plasmid DNA.
- Adjust reagent volumes proportionally according to the culture volume (mini, midi, or maxi prep).
- Mix buffers and cell suspensions gently by inversion rather than vortexing to prevent shearing of genomic DNA.
- Keep all buffers and samples cold when indicated to maintain plasmid stability and reduce nuclease activity.
- Always include RNase A in the resuspension buffer to remove RNA contamination.
- Perform centrifugation steps at recommended speeds and temperatures to ensure efficient separation of plasmid DNA.
- Validate plasmid DNA quality and concentration post-prep by spectrophotometry or gel electrophoresis.



Materials

Resuspension Buffer

- 50 mM Tris base
- 10 mM EDTA
- 100 µg/ml of RNase (optional)
- pH 8.0

Lysis Buffer

- 200 mM NaOH
- 1% SDS

Neutralization Buffer

- 3M Potassium acetate
- pH 5.0

Elution Buffer

- 10 mM Tris base
- 1 mM EDTA
- pH 8.0

Isopropanol

Ethanol

Sterile Water

Shaking incubator

Laboratory centrifuge

Safety warnings

- ! ▪ Sodium hydroxide (NaOH) and sodium dodecyl sulfate (SDS) are caustic reagents: wear gloves, protective eyewear, and lab coat when handling.
- Avoid vigorous vortexing or shaking after lysis buffer addition to reduce shearing of chromosomal DNA that contaminates plasmid preps.
- Do not exceed recommended incubation times during lysis to prevent irreversible denaturation of plasmid DNA.
- Handle pellets gently during washing and resuspension to avoid loss of plasmid DNA.
- Dispose of bacterial waste and chemical reagents according to institutional biosafety and chemical safety guidelines.



Before start

- Prepare all buffers fresh or ensure they are properly stored and equilibrated to appropriate temperatures before starting.
- Pre-cool centrifuges, tubes, and neutralization buffer to 4°C or on ice as required.
- Grow *E. coli* cultures overnight in LB medium with appropriate antibiotics to reach stationary phase (14–16 hours, 37°C, shaking).
- Confirm culture density to avoid overloading the lysis step, which can reduce purity and yield.
- Have all materials including pipettes, filter tips, tubes, and reagents ready and labeled for efficient workflow.



Bacterial culture

1 Inoculate *E. coli* harboring the plasmid into LB broth with appropriate antibiotic.

- Mini prep: 2-5 mL
- Midi prep: 25-100 mL
- Maxi prep: 250-1000 mL

2 Incubate Overnight at 150 rpm, 37°C



Bacterial harvesting

2m

3 Centrifuge

- Mini prep: 12000 x g, 4°C, 00:01:00
- Midi prep: 10000 x g, 4°C, 00:10:00
- Maxi prep: 10000 x g, 4°C, 00:10:00

2m



4 Discard supernatant completely.

Resuspension

5 Resuspend cell pellet in resuspension buffer:

- Mini prep: 250 µL
- Midi prep: 4 mL
- Maxi prep: 10 mL

6 Mix gently by pipetting until cells are completely resuspended.

Alkaline lysis

4m

7 Add lysis buffer to resuspended cells:

- Mini prep: 250 µL
- Midi prep: 4 mL
- Maxi prep: 10 mL



8 Mix gently by inverting the tube 4-6 times. Avoid vortexing to prevent shearing genomic DNA.

9 Incubate at Room temperature for 00:04:00 . Do not exceed 5 minutes to prevent irreversible denaturation.

4m



Neutralization

10m

10 Add ice-cold neutralization buffer

▪ Mini prep: 350 μ L

▪ Midi prep: 6 mL

▪ Maxi prep: 15 mL

11 Mix gently by inverting until a thick white precipitate forms.

12 Incubate on ice for 00:10:00 to enhance precipitate formation.

10m



Lysate clearing.

13 Centrifuge



▪ Mini prep: 12000 x g, 4°C, 00:10:00

▪ Midi prep: 10000 x g, 4°C, 00:15:00

▪ Maxi prep: 10000 x g, 4°C, 00:15:00

14 ▪ Carefully transfer the clear supernatant to a fresh tube without disturbing the pellet.

Precipitation of plasmid DNA

30m

15 Add ice-cold isopropanol

▪ Mini prep: 250 μ L

▪ Midi prep: 4 mL

▪ Maxi prep: 10 mL

16 Mix by gentle inversion.



17 Incubate at -20 °C for 00:30:00

30m



18 Centrifuge

45m

- Mini prep: 12000 x g, 4°C, 00:15:00
- Midi prep: 10000 x g, 4°C, 00:15:00
- Maxi prep: 10000 x g, 4°C, 00:15:00



19 ▪ Discard supernatant gently without disturbing pellet.

Wash DNA pellet

5m

20 Add 1 mL of [M] 70 % (v/v) ice-cold ethanol.

21 Centrifuge

- Mini prep: 12000 x g, 4°C, 00:15:00
- Midi prep: 10000 x g, 4°C, 00:15:00
- Maxi prep: 10000 x g, 4°C, 00:15:00



22 ▪ Remove ethanol carefully and air-dry pellet for 00:05:00 . Avoid over-drying.

5m



Resuspension

5m

23 Resuspend DNA pellet TE buffer or nuclease-free water.

- Mini prep: 50-100 µL
- Midi prep: 500-2000 µL
- Maxi prep: 5-20 mL

24 ▪ Incubate at 37 °C or Room temperature for 00:05:00 until fully dissolved.

5m





Protocol references

Birnboim, H. C., & Doly, J. (1979). A rapid alkaline extraction procedure for screening recombinant plasmid DNA. *Nucleic acids research*, 7(6), 1513–1523. <https://doi.org/10.1093/nar/7.6.1513>.

Green, M. R., & Sambrook, J. (2018). Preparation of Plasmid DNA by Alkaline Lysis with Sodium Dodecyl Sulfate: Maxipreps. *Cold Spring Harbor protocols*, 2018(1). <https://doi.org/10.1101/pdb.prot093351>.