

Mar 14, 2025

Alpha-synuclein Seed Amplification Assay

DOI

<https://dx.doi.org/10.17504/protocols.io.kxygxwj4v8j/v1>

Sunil Srivastav¹, Davin Henderson², Jessika Suescun Ocampo¹, Divya Onkar Mondhe¹, * Richard Gordon¹

¹School of Biomedical Science, Translational Research Institute, Queensland University of Technology, Brisbane, 4102, Australia;

²CWD Evolution LLC, Fort Collins, CO, United States



Sunil Srivastav

Queensland University of Technology

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.kxygxwj4v8j/v1>

Protocol Citation: Sunil Srivastav, Davin Henderson, Jessika Suescun Ocampo, Divya Onkar Mondhe, * Richard Gordon 2025. Alpha-synuclein Seed Amplification Assay. **protocols.io** <https://dx.doi.org/10.17504/protocols.io.kxygxwj4v8j/v1>

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited



Protocol status: Working

We have been using this protocol successfully, and it is working well.

Created: March 04, 2025

Last Modified: March 14, 2025

Protocol Integer ID: 123808

Keywords: seed amplification, alpha synculein, synuclein seed amplification assay preparation of fecal sample, synuclein seed amplification assay preparation, synuclein seed amplification, fecal sample, synuclein

Abstract

Preparation of fecal samples for alpha-synuclein seed amplification assay.

Materials

1. RT-QulC buffer
2. alpha-synuclein substrate
3. optic bottom 96 well plates
4. Fluorometer

Protocol

- 1 Homogenize fecal samples in phosphate buffered saline
- 2 Centrifuge the homogenate and carefully remove the supernatant.
- 3 To the supernatant add metal-ion nanoparticles and resuspend them.
- 4 Next, perform magnetic separation of α -synuclein seeds.
- 5 Add these seeds to the 96 well plate and now resuspend them with Real-time quaking-induced conversion (RT-QuIC) buffer containing α -synuclein substrate.
- 6 Using a standard fluorometer, take fluorescence readings with an excitation of 450 nm and emission of 480 nm.
- 7 After a few hours are elapsed, an increase in fluorescence intensity could be observed, indicating alpha-synuclein seed amplification.

Protocol references

1. L. Concha-Marambio, S. Pritzkow, M. Shahnawaz, C. M. Farris, C. Soto, Seed amplification assay for the detection of pathologic alpha-synuclein aggregates in cerebrospinal fluid. *Nat Protoc* **18**, 1179-1196 (2023).
2. D. M. Henderson *et al.*, Detection of chronic wasting disease prion seeding activity in deer and elk feces by real-time quaking-induced conversion. *J Gen Virol* **98**, 1953-1962 (2017).
3. N. D. Denkers, D. M. Henderson, C. K. Mathiason, E. A. Hoover, Enhanced prion detection in biological samples by magnetic particle extraction and real-time quaking-induced conversion. *J Gen Virol* **97**, 2023-2029 (2016).