


Dec 11, 2020

# Amplification and Pooling

In 5 collections

DOI

[dx.doi.org/10.17504/protocols.io.beqkjduw](https://dx.doi.org/10.17504/protocols.io.beqkjduw)



Franziska Aron<sup>1</sup>, Guido Brandt<sup>2</sup>

<sup>1</sup>Friedrich-Schiller Universität Jena; <sup>2</sup>Max Planck Institute for the Science of Human History

WarinnerGroup

MPI EVA Archaeogenetics



Franziska Aron

Friedrich-Schiller Universität Jena, Max Planck Institute fo...

OPEN ACCESS



DOI: [dx.doi.org/10.17504/protocols.io.beqkjduw](https://dx.doi.org/10.17504/protocols.io.beqkjduw)

**Protocol Citation:** Franziska Aron, Guido Brandt 2020. Amplification and Pooling . **protocols.io**

<https://dx.doi.org/10.17504/protocols.io.beqkjduw>

**License:** This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

**Protocol status:** Working

**We use this protocol and it's working**

**Created:** April 07, 2020

**Last Modified:** December 11, 2020

**Protocol Integer ID:** 35308

**Keywords:** DNA library, NGS, dual-index, ancient DNA, sequencing, nonUDG, double-stranded, DNA, genomic DNA, genomics, palaeogenetics, archaeogenetics, paleogenetics, archeogenetics, aDNA, Illumina, library preparation, nucleic acids, Amplification, PCR, Index Amplification



## Abstract

This protocol describes the amplification procedure of dual-indexed double-stranded DNA libraries, for shotgun Illumina sequencing. It is typically used for libraries indexed using the following protocol:

(<https://dx.doi.org/10.17504/protocols.io.bakticwn>)

## Image Attribution

Franziska Aron

## Guidelines

### **Working in an Molecular Biology Laboratory**

This protocol can place in a typical DNA-based molecular biology lab.

Please keep in mind the safety guidelines of your specific country and institution.

Recommendations include wearing of:

- lab coats
- closed shoes and trousers
- safety glasses
- nitril or latex gloves

## Materials

### MATERIALS

✕ 0.2 ml PCR Tube strips **Eppendorf Catalog #0030124359**

✕ DNA LoBind Tube 1.5ml **Eppendorf Catalog #022431021**

✕ 2 ml LoBind Tubes **Eppendorf Catalog #0030108078**

✕ Eppendorf Tubes® 5.0 mL with snap cap **Eppendorf Catalog #30119460**

✕ dNTP Mix (25 mM each) **Thermo Fisher Scientific Catalog #R1121**

✕ Sodium Acetate buffer solution 3M pH 5.2 for molecular biology **Merck MilliporeSigma (Sigma-Aldrich) Catalog #S7899-500ML**

✕ Tween 20 **Merck MilliporeSigma (Sigma-Aldrich) Catalog #P9416-50ML**

✕ Water HPLC Plus **Merck MilliporeSigma (Sigma-Aldrich) Catalog #34877-2.5L-M**

✕ D1000 Ladder **Agilent Technologies Catalog #5067-5586**

✕ D1000 ScreenTape **Agilent Technologies Catalog #5067-5582**

✕ D1000 Reagents **Agilent Technologies Catalog #5067-5583**

✕ Herculase II Fusion DNA Polymerase **Agilent Technologies Catalog #600679**

✕ High Sensitivity D1000 Ladder **Agilent Technologies Catalog #5067-5587**

✕ High Sensitivity D1000 Reagents **Agilent Technologies Catalog #5067-5585**

✕ High Sensitivity D1000 ScreenTape **Agilent Technologies Catalog #5067-5584**

✕ MinElute PCR Purification Kit **Qiagen Catalog #28004**

### Primers

	Oligo_ID	Sequence (5'-3')	Concentration
	IS5	AATGATACGGCGACCACCGA	10 µM
	IS6	CAAGCAGAAGACGGCATACGA	10 µM

### Lab equipment

PCR Thermocycler (e.g. Eppendorf Thermomaster Nexus)

Centrifuge 1.5/2.0 ml (e.g. Eppendorf 5424)

Rotor 1.5/2.0ml (e.g. Eppendorf F-45-24-11)

Mini table centrifuge

TapeStation (e.g. Agilent Technologies, 4200 TapeStation System, SKU: G2991AA)

Vortex mixer (e.g. Scientific Industries Vortex-Genie® 2)

## Safety warnings

### ! Reagents

*Sodium Acetate*

- H139: Causes serious eye irritation



*Ethanol*

- H225 Highly flammable liquid and vapour.
- H319 Causes serious eye irritation.



*Guanidinium hydrochloride (GuHCl)* (in PB buffer of Qiagen MinElute kit)

- H302 Harmful if swallowed.
- H332 Harmful if inhaled.
- H315 Causes skin irritation.
- H319 Causes serious eye irritation.



### Kits

Check manufacturer's safety information for the TapeStation Kits used in this protocol.

Check manufacturer's safety information for the MinElute PCR Purification kit used in this protocol.

- Note that PBI must be stored at room temperature in the dark. PBI is light sensitive.






## Before start

### Planning

This protocol takes 1 day.

Check all waste disposal guidance for all reagents in this protocol against your corresponding laboratory regulations.

### Preparation of buffers (Qiagen MinElute kit):

- Add ethanol to PE wash buffer according to manufacturer's instructions.
- Add  200  $\mu\text{L}$  pH-Indicator and  300  $\mu\text{L}$  Sodium Acetate to  48.5 mL of PB binding buffer. This solution is referred to as **PBI** throughout the protocol. Must be stored at room temperature in the dark. PBI is light sensitive.
- Add Tween-20 to EB elution buffer to a final concentration of 0.05% Tween-20 in EB. This solution is referred to as **EBT** throughout the protocol.

### Equipment

Make sure all necessary equipment is available (see Materials).

### Abbreviations



EBT = modified EB-Buffer (MinElute Kit), see Preparation of buffers

HPLC = High Performance Liquid Chromatography (-Grade Water)


PBI = modified PB-Buffer (MinElute Kit), see Preparation of buffers

PE = PE-Buffer from Qiagen MinElute Kit

### Samples

This protocol is designed for the amplification of indexed libraries as prepared by the protocol described in (<https://dx.doi.org/10.17504/protocols.io.bakticwn>). The indexing protocol generates  50  $\mu\text{L}$  of indexed library, of which  20  $\mu\text{L}$  will be used for this protocol. Ensure sufficient indexed library is available before starting this protocol.

## Calculations

1 Prepare amplification assay [  100 µL per reaction ]




Based on the quantification results of the indexed libraries (<https://dx.doi.org/10.17504/protocols.io.bakticwn>) calculate the number of PCR cycles (amplification factor) needed to reach  $10^{13}$  copies of DNA per indexed Library.

### Note

Formula in Excel to get the Cycles needed

=LOG((1\* $10^{13}$ /Copies per rxn),2)  
(log base 2)

**Example:** The following calculation is for  5 µL per reaction, with two indexed library samples (A and B) having different concentrations of DNA copies.

Optional Changes: 1.If the Calculation shows up less then 3 Cycles, you also have the Option to add less then 5 µL.

2. Instead of 4 reactions of 5 µL each you can also split in 8 reactions of 2 µL each

Instead of  $4 \times 5 \mu\text{L}$  reactions you can also split in  $8 \times 2 \mu\text{L}$  reactions

Sample Name	Copies per µL	µL per rxn	Copies per rxn	Cycles needed	Real Cycles	Amplification Factor	Output per rxn [Copies]
A	7.32E+10	5	3.66E+11	4.772998	5	32	1.17E+13
B	5.79E+06	5	2.32E+07	18.72018251	19	524288	1.21E+13

### Note

Do not calculate the amount of cycles for a higher amount of copies than  $1.4 \times 10^{13}$  to avoid heteroduplexes.

## Preparation

- 2 Prepare cleaned workspace with all necessary reagents and equipment.

### Note




Label all  0.2 mL PCR strips for the PCR reactions.

## PCR


- 3 Set up four amplification reactions of  100 µL each per library

Reagent	Stock concentration	Final concentration	1x Volume [µl]
Herculase II Reaction buffer	5x	1x	20
IS5 primer	10 µM	0.4 µM	4
IS6 primer	10 µM	0.4 µM	4
dNTP's	25 mM	0.25 mM	1
Herculase II Fusion	1 U	0.01 U	1
DNA			5
HPLC-Water			65
<b>Total</b>			<b>100</b>

- 3.1 Vortex master mix before adding the enzyme. After adding the enzyme, mix by pipetting or inverting the tube.

- 3.2 Pipette  95 µL mastermix and  5 µL indexed library into each tube (use  0.2 mL PCR strips).

### Safety information

Keep the remaining library at  -20 °C until further use.

4 Amplify in a thermocycler with the following program:



Tempreature	Time	
95°C	2 min	Initial denaturation
95°C	30 sec	Cycles (see Step 1)
60°C	30 sec	
72°C	30 sec	
72°C	5 min	Final elongation
Finally hold the reactions at 10 °C.		


### Note




Adjust the number of cycles according to the amplification factor as calculated in step 1.

During this incubation take **MinElute columns** out of the fridge so they warm up to **room temperature** before use in the next step.

This is an ideal point to prepare downstream steps, including labelling of final elution tubes, MinElute Columns etc.

## MinElute Purification

5 Purify with MinElute kit with the following modifications to the manufacturer's protocol: Use one column for all four reactions [=  400 µL *PCR product*] of a sample.

5.1 Add  2400 µL PB or PBI\* buffer to a  5 mL tube for each sample (this is  600 µL buffer for each PCR reaction). Add all 4 PCR reactions per sample to the










same tube with PB buffer and vortex briefly.

#### Safety information

After the PCR product is mixed with the PBI, the PBI should keep its yellow colour. If it turns purple the pH is too high and the efficiency of the MinElute columns is not guaranteed.



- 5.2 Load  700  $\mu$ L of the mixture onto one MinElute column, incubate for  00:02:00 , spin  15800 x g, 00:01:00 , and discard flow-through.

#### Note

Pour off the liquid into a waste tube, and pat the rim of the collection tube dry on a paper tissue or towel. Use just one spot on the paper tissue per sample. Be careful not to touch the rim of the tube on the waste container. After you are finished with all samples, discard the paper and wipe clean the surface underneath with water and soap.




- 5.3 Repeat loading until the complete mixture was run through the column.

 [go to step #5.2](#)

- 5.4 Add  700  $\mu$ L PE (wash) buffer, spin  15800 x g, 00:01:00 , and discard flow-through.

- 5.5 Dry spin  15800 x g, 00:01:00 ,

- 5.6 Put column into new  1.5 mL LoBind tubes.

- 5.7 Add  50  $\mu$ L EBT buffer to the center of the filter, incubate for  00:02:00 , and spin  15800 x g, 00:01:00 to elute the amplified indexed library.

#### Note

Carefully pipette EBT directly onto the center of the membrane without touching the membrane.



## Measurement and Dilution

- 6 Dilute amplified index library 1:10 with HPLC- water and check for fragment size, concentration, and heteroduplexes. (for example with the D1000 Kit's Tape, Reagent and Buffer - following the manufacturer's protocol on the TapeStation)

#### Note

if you see heteroduplexes you need to perform a reconditoning PCR.



Reconditoning PCR: one cycle PCR using 100 ng library template in a 100 µl Herculase PCR reaction (same set up as in 3) and amplified with 1 cycles of 95°C for 2 min, 58°C for 2 min, and 72°C for 5 min. Purify with MinElute kit following the instructions from Step 5, but elute in 20µl EBT.

- 7 Dilute each amplified indexed library to  10 nM with EBT buffer or HPLC-water for shotgun sequencing. Then pool the  10 nM amplified indexed libraries in equimolar amounts (take the same volume for each sample).

#### Note

The final concentration of a pool of several amplified indexed libraries should be

 10 nM .

- 8 Check the  10 nM library or the  10 nM library pool for the correct concentration, (for example with the HighSensitivity D1000 Kit's Tape, Reagent and Buffer following the manufacturer's protocol on the TapeStation.)