

🌐 Amplification Free Paired End Library Construction Protocol




[dx.doi.org/10.17504/protocols.io.bd3ti8nn](https://doi.org/10.17504/protocols.io.bd3ti8nn)

¹The Earlham Institute, Norwich Research Park, Norwich, NR4 7UZ, United Kingdom

 **Graham Etherington**
The Earlham Institute



Created: March 23, 2020

























Last Modified: March 26, 2020

Protocol Integer ID: 34643



Keywords: polecat, vertebrate, non-model organism, Illumina, chromium, Bionano, assembly, sequencing,

Abstract

Amplification Free Paired End Library Construction Protocol.

- 1 A total of  600 ng of DNA was sheared in a  60 μL volume on a Covaris S2 (Covaris, Massachusetts, USA) for 1 cycle of  00:00:40 with a duty cycle of 5%, cycles per burst of 200 and intensity of 3.
- 2 The fragmented molecules were then end repaired in  100 μL volume using the NEB End Repair Module (NEB, Hitchin, UK) incubating the reaction at  22 $^{\circ}\text{C}$ for  00:30:00 .
- 3 Post incubation  58 μL beads of CleanPCR beads (GC Biotech, Alphen aan den Rijn, The Netherlands) were added using a positive displacement pipette to ensure accuracy and the DNA precipitated onto the beads.
- 4 This is then washed twice with 70% ethanol and the end repaired molecules eluted in  25 μL Nuclease free water (Qiagen, Manchester, UK).
- 5 End repaired molecules were then A tailed in  30 μL volume using in the NEB A tailing module (NEB) incubating the reaction at  37 $^{\circ}\text{C}$ for  00:30:00 .
- 6 To the A tailed library molecules  1 μL of an appropriate Illumina TruSeq Index adapter (Illumina, San Diego, USA) is added and mixed, then  31 μL of Blunt/ TA ligase (NEB) is added and incubated at  22 $^{\circ}\text{C}$ for  00:10:00 .
- 7 Post incubation  5 μL of stop ligation is added and the reaction incubated at  Room temperature for  00:05:00 .
- 8 Following this incubation  67 μL beads of CleanPCR beads (GC Biotech, Alphen aan den Rijn, The Netherlands) were added and the DNA precipitated onto the beads.
- 9 The samples are then washed twice with 70% ethanol and the end repaired molecules eluted in  100 μL nuclease free water.



- 10 Two further CleanPCR bead based purifications were undertaken to remove any adapter dimer molecules that may have formed during the adapter ligation step. The first with 0.9x volume beads, the second with 0.6x and the final library eluted in  25 μL Resuspension Buffer (Illumina).
- 11 Library QC was performed by running a  1 μL aliquot on a High Sensitivity BioAnalyser chip (Agilent, Stockport, UK) and the DNA concentration measured using the High Sensitivity Qubit (Thermo Fisher, Cambridge, UK).
- 12 To determine the number of viable library molecules the library was subjected to quantification by the Kappa qPCR Illumina quantification kit (Kapa Biosystems, London, UK) and a test lane run at 10pM on a MiSeq (Illumina) with 2×300bp reads to allow the library to be characterised prior to generation of the 60x coverage required on the HiSeq2500s (Illumina) with a 2×250bp read metric.