

Aug 07, 2025

Amplifying Bacteriophage

DOI

<https://dx.doi.org/10.17504/protocols.io.36wggpjdvk5/v1>

Madison Altieri¹

¹Bowling Green State University



Madison Altieri

Bowling Green State University

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.36wggpjdvk5/v1>

Protocol Citation: Madison Altieri 2025. Amplifying Bacteriophage. protocols.io
<https://dx.doi.org/10.17504/protocols.io.36wggpjdvk5/v1>

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: August 07, 2025



Last Modified: August 07, 2025

Protocol Integer ID: 224266

Keywords: amplifying bacteriophage, amplification phage, phage, liquid culture, spot test

Abstract

This is for the amplification phage by making webbed plates (Path A; ~3days) or by liquid culture (~2days). This protocol also details how to titer phage through a spot test.

Materials

Phage Buffer

- 100mL ddH₂O
- 2mL 100mM MgSO₄ (or 200μL 1M)
- 1mL 500mM CaCl₂
- 1mL 1M Tris-Cl
- filter sterilize 0.22μm

T-top

- 10g Tryptone
- 8g NaCl
- 1L ddH₂O
- Mix and aliquot 100mL into bottles
- Into bottles, put 0.75g agar
- autoclave

T-broth

- 10g Tryptone
- 8g NaCl
- 1L ddH₂O

T-plate

- 10g Tryptone
- 8g NaCl
- 15g Agar
- 1L ddH₂O

Day 1 Path A: Plating Bacteria and Phage to make a webbed plate

15m

1 Incubate Bacteria and Phage

100µL of overnight culture

100µL of bacteriophage culture

You will need to adjust this accordingly based on the phage titer.

Time of incubation will also be dependent on the phage: ~00:10:00-00:30:00

15m



2 Making Top agar (With nutrient deficient T-plates)

Make melted agar from T-top, T-broth, and 1M CaCl₂

Each plate should have a total of 3mL top agar

- 1:1 ratio of T-top and T-broth and for every 10mL 100µL of CaCl₂
- Mix 1M CaCl₂ and T-broth
- Add melted T-top
- Ex.) if you were making 4 plates (4 * 3 = 12mL). Meaning you would add 6mL T-top, 6mL T-broth, and ~100µL CaCl₂


1d

Take the 200µL phage and bacterial mix and mix with 3mL top agar - pour immediately onto T-plate

Sit at room temperature to solidify

Typically 2-3 plates will give a good titer

Note: make sure that the top agar is not too hot that it will immediately kill the bacteria around 40-60°C depending on bacteria



Webbed Plate



Day 2 Path A: Scrape Method

15m

3 Prepare items

- Beaker with ethanol and metal scraping device
- Conical tube with 5mL phage buffer
- 0.22µm filter and syringe
- Sterile conical tube for final filtrate

4 Scraping

15m

Sterilize the scraper with ethanol and bunsen burner - wait to cool

Physically scrape off the top bacteria and phage from the T-plate and place in the 5mL of phage buffer. If you have more than one plate combine all scrapings into the same 5mL phage buffer

Centrifuge 4°C 00:10:00 5,000 RPM

- Filter supernatant with 0.22µm
- Store phage at 4°C
- To freeze phage in -80°C: 1mL phage + 1mL 50% Glycerol

Day 1 Path B: Liquid Amplification of phage

4h 45m

5 Liquid Culture

45m

In a test tube:

5mL of desired broth (usually LB)

100µL Bacteria

10µL 1M CaCl₂

Incubate until log phase (*E. coli* ~00:45:00)

6 Bacteriophage

4h

Add Bacteriophage: the amount you add will depend on your titer. For *E. coli* lytic phage, 10µL will result in amplified phage.

Note: incubation times after adding bacteriophage vary depending on host bacteria. For E. coli ~04:00:00 will give you titers of 10⁸⁺ phage.

7 Filter

Filter the amplified phage with a 0.22µm filter into a sterile conical tube

- Store phage at 4°C
- To freeze phage in -80°C: 1mL phage + 1mL 50% Glycerol

Checking Phage Titer

1d

8 Making Top agar (With nutrient deficient T-plates)

Make melted agar from T-top, T-broth, and 1M CaCl₂

Each plate should have a total of 3mL top agar



- 1:1 ratio of T-top and T-broth and for every 10mL 100μL of CaCl₂
- Mix 1M CaCl₂ and T-broth
- Add melted T-top
- Ex.) if you were making 4 plates (4 * 3 = 12mL). Meaning you would add 6mL T-top, 6mL T-broth, and ~100μL CaCl₂
- Take the 100μL host bacteria and mix with 3mL top agar - pour immediately onto T-plate

Sit at room temperature to solidify

9 **Serial Dilution**

Perform 10-fold serial dilution of filtered phage out to -8 to -12

Spot 5μL of each serial dilution onto the solidified T-plate

Incubate overnight

Note: if you want to have a more accurate titer, serial dilute and mix onto an entire plate with 100μL host bacteria and 3mL top agar

1d