



Jan 31, 2024

Amyloid beta (A β) aggregates N-terminal labeling

DOI

<https://dx.doi.org/10.17504/protocols.io.kxygx3dk4g8j/v1>

Patricia Yuste Checa¹, F Ulrich Hartl¹

¹Department of Cellular Biochemistry, Max Planck Institute of Biochemistry



Patricia Yuste-Checa

MPI Biochemistry

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.kxygx3dk4g8j/v1>

Protocol Citation: Patricia Yuste Checa, F Ulrich Hartl 2024. Amyloid beta (A β) aggregates N-terminal labeling. **protocols.io** <https://dx.doi.org/10.17504/protocols.io.kxygx3dk4g8j/v1>

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: January 25, 2024

Last Modified: May 31, 2024

Protocol Integer ID: 94480

Keywords: ASAPCRN, using amyloid beta, amyloid beta, amyloid beta as example, protein aggregate, a β , labeling, terminal labeling, label, terminal labeling this protocol detail, protein

Funders Acknowledgements:

Aligning Science Across Parkinson's

Grant ID: ASAP-000282

Abstract

This protocol details how to efficiently label protein aggregates at the N-terminal using Amyloid beta as example.

Attachments



[974-2531.docx](#)

19KB

Materials

Buffers:

- Labeling buffer: 1M 0.1 Mass Percent sodium bicarbonate buffer pH 8.3
- 1x PBS pH 7.2
- Dimethyl sulfoxide (DMSO)

 Alexa 488 NHS ester Thermo Fisher Scientific Catalog #A20000

 pHrodo Red Succinimidylester Thermo Fisher Scientific Catalog #P36600



Amyloid beta (A β) aggregates N-terminal labeling

30m

- 1 Centrifuge A β aggregates (e.g.,  1 mL at  10 micromolar (μ M)) at  20000 x g, 4°C, 00:10:00 . 
- 2 Wash the pellet containing the aggregates with 1x PBS  7.2 . 
- 3 Centrifuge at  20000 x g, 4°C, 00:10:00 . 
- 4 Wash the pellet with Labeling buffer (refer materials section). 
- 5 Centrifuge at  20000 x g, 4°C, 00:10:00 . 
- 6 Resuspend the pellet in  200 μ L Labeling buffer.
- 7 Dissolve the dye (Alexa488 (A488) NHS ester, Thermo Fisher Scientific, A20000; pHrodo Red Succinimidylester, Thermo Fisher Scientific, P36600) in DMSO.
- 7.1 With a pipette tip gently touch the dye powder which will stick to the tip. 
- 7.2 Immerse the tip in some DMSO previously dispensed in a tube.
- 7.3 Repeat the procedure several times until the solution reaches the desired color.

Note

If labeling a protein for uptake assays analyzed by flow cytometry, it is recommended to use A488 because the A488 signal outside the cell can be easily quenched by adding Trypan blue right before measurement. pHrodo red dye is a pH sensitive dye which fluoresces brightly only in acidic environments and therefore can be used to specifically monitor phagocytosis and endocytosis.



- 8 Quantify diluted dye concentration by nanodrop. Dilute the sample in water to reach a $\lambda < 1$ for an accurate measurement.

Note

Physical characteristics of the dyes to be set in the nanodrop:

Alexa488: Absorbance maximum (λ_{max}): 495 nm; Extinction coefficient (ϵ): 71,000 $\text{cm}^{-1}\text{M}^{-1}$; Correction factor at 280 nm (CF_{280}): 0.11; Correction factor at 260 nm (CF_{260}): 0.3.

pHrodo Red: Absorbance maximum (λ_{max}): 560 nm; Extinction coefficient (ϵ): 65,000 $\text{cm}^{-1}\text{M}^{-1}$; Correction factor at 280 nm (CF_{280}): 0.12; Correction factor at 260 nm (CF_{260}): 0.36.

- 9 Add the corresponding amount of dye to a final protein:dye ratio of 1:10 (e.g.,  50 μL of dye at  2.5 millimolar (mM) for  200 μL protein at  50 micromolar (μM)).



- 10 Incubate  01:00:00 at  Room temperature in the dark.

1h



- 11 Centrifuge at  20000 x g, 4°C for  00:10:00 .

10m



- 12 Wash the pellet with  1 mL methanol to remove excess dye (just when labeling with pHrodo dye).

**Note**

Skip the methanol washing step if using Alexa dye for labeling.

- 13 Centrifuge at  20000 x g, 4°C for  00:10:00 .

10m



- 14 Wash the pellet 2 to 4 times with 1x PBS  7.2 .





- 15 Resuspend the pellet with 1x PBS $\text{pH } 7.2$ or desired final buffer (e.g., in $200 \mu\text{L}$ buffer resulting in $50 \mu\text{M}$ labeled aggregates).
- 16 Sonicate labeled aggregates in a Bioruptor sonication bath (Diagenode) (5 cycles of 5 seconds on – 5 seconds off), or similar.
- 17 Aliquot, flash-freeze in liquid nitrogen and store at -80°C .

Note

This protocol can be used for other aggregates like Tau or α -Synuclein aggregates.