



Mar 27, 2025

Amyloid beta (M1-42) Aggregation Monitored by Thioflavin-T (ThT) Fluorescence in a Plate Reader

DOI

<https://dx.doi.org/10.17504/protocols.io.bp2l68ky5gqe/v1>

Patricia Yuste-Checa¹, F Ulrich Hartl¹

¹Department of Cellular Biochemistry, Max Planck Institute of Biochemistry



Patricia Yuste-Checa

MPI Biochemistry

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.bp2l68ky5gqe/v1>

Protocol Citation: Patricia Yuste-Checa, F Ulrich Hartl 2025. Amyloid beta (M1-42) Aggregation Monitored by Thioflavin-T (ThT) Fluorescence in a Plate Reader. **protocols.io** <https://dx.doi.org/10.17504/protocols.io.bp2l68ky5gqe/v1>

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: March 20, 2025

Last Modified: March 27, 2025

Protocol Integer ID: 124796

Keywords: ASAPCRN, amyloid beta, fluorescence, plate reader, thioflavin, a β

Funders Acknowledgements:

Aligning Science Across Parkinson's

Grant ID: ASAP-000282

Abstract

This protocol is based on Linse S. *Methods Mol Biol.* (2020), section "3.10 Preparation of Samples for Kinetics" and details how to efficiently monitor A β (M1–42) aggregation by thioflavin T fluorescence in a plate reader.

Materials

Buffers:

- **Gu6:** pH pH 8.5

	A	B
	Guanidinium hydrochloride	6 M
	Sodium phosphate	20 mM

- **PE 8.5:** pH pH 8.5

	A	B
	Sodium phosphate	20 mM
	EDTA	0.2 mM

- [M] 1 millimolar (mM) Thioflavin T (ThT) in water. This can be stored at 🌡️ -20 °C .

Equipment

- CLARIOstar plate reader (BMG Labtech)



Corning® 96-well Half Area Black/Clear Flat Bottom Polystyrene NBS Microplate, 25 per Bag, without L Corning Catalog #3881



Size Exclusion Chromatography

1

Note

To discard any possible aggregate that could have formed during the freezing lyophilization process, size exclusion chromatography of the purified peptide should be performed before the kinetic experiment.

Dissolve  0.5 mg of lyophilized A β (M1-42) in  700 μ L Gu6 buffer.

Note


See protocol "Purification of recombinant Amyloid beta peptide (M1-42) from *Escherichia coli*".

2 Load the sample onto a Cytiva Superdex 75 Increase 10/300 GL ( 24 mL) SEC column previously equilibrated with PE 8.5 buffer.

3 Run the column with PE 8.5 at  0.5 μ L .

Note

NOTE: Collect the fractions preferably in low binding tubes.



Size exclusion chromatogram

- 4 $A\beta$ (M1-42) monomers elute at around 14-15 mL (fractions around 27-31).
- 5 Measure absorbance of fractions at 205 nm wavelength using a UV spectrophotometer using PE 8.5 buffer as blank.

Note

NOTE: $A\beta$ (M1-42) contains just one tyrosine and no tryptophan. Therefore, measurement at 280 nm wavelength of low concentrated peptide may not be accurate.

- 6 Pool fractions containing $A\beta$ (M1-42) monomers and measure the concentration at A_{205nm} .
- 7 Prepare aggregation reactions:
 - Mix reagents to a final concentration of 4 micromolar (μ M) $A\beta$ (M1-42),
 10 micromolar (μ M) ThT in PE 8.5 buffer.
 - Prepare a mix for 4.5 reactions per condition (quadruplicates in each plate as technical replicates), final volume 80 μ L per well.



**Note**


Mix of  360 μL =  80 μL reaction x 4.5.

- Concentration of A β (M1–42) and ThT should be tested empirically in order to achieve a final concentration in a range that provides a linear response of the fluorescence intensity versus fibril mass concentration, see Linse S. Methods Mol Biol. (2020).
- Molecular or chemical chaperones can be included in the reaction in order to study their effect on A β (M1–42) aggregation.
- Keep A β (M1–42) peptide as well as the 96 well plate during preparation always on ice.

- 8 Dispense  80 μL of the mix into a well of a 96-well half-area plate of black polystyrene with a clear bottom (Corning 3881).
- 9 If possible, the outer wells should not be used and should be filled with water.
- 10 Seal the plate with parafilm to avoid excessive evaporation.
- 11 In a CLARIOstar plate reader (BMG Labtech) set the following parameters and start the reaction:
 - Fluorescence measurement: ThT signal, excitation 440 nm, emission 480 nm, measured every  00:03:00 .
 - Temperature  37 $^{\circ}\text{C}$.

Note

NOTE: Under these conditions A β (M1–42) rapidly aggregates, reaching a fluorescence plateau after approximately  02:00:00 aggregation. The data can be fitted using Sigma plot software (Sigmoidal, Sigmoid, 3 Parameter function) to obtain the half time for reaching the aggregation plateau.



Protocol references

Linse S. Expression and Purification of Intrinsically Disordered A β Peptide and Setup of Reproducible Aggregation Kinetics Experiment. *Methods Mol Biol.* 2020;2141:731 754. doi: 10.1007/978-1-0716-0524-0_38. PMID: 32696387.