

Jun 20, 2022

Version 5

# An X-HTDC method for estimating particulate phosphorus from microalgae V.5

DOI

[dx.doi.org/10.17504/protocols.io.kqdg35dq7v25/v5](https://dx.doi.org/10.17504/protocols.io.kqdg35dq7v25/v5)



Ying-Yu Hu<sup>1</sup>, Zoe V. Finkel<sup>1</sup>

<sup>1</sup>Dalhousie University

Marine Microbial Macroecology Lab

Tech. support email: [ruby.hu@dal.ca](mailto:ruby.hu@dal.ca)



Ying-Yu Hu

Dalhousie University

## Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.kqdg35dq7v25/v5>

**Protocol Citation:** Ying-Yu Hu, Zoe V. Finkel 2022. An X-HTDC method for estimating particulate phosphorus from microalgae. **protocols.io** <https://dx.doi.org/10.17504/protocols.io.kqdg35dq7v25/v5> Version created by **Ying-Yu Hu**



**License:** This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

**Protocol status:** Working

**The detailed study of this protocol has been submitted for peer review.**

**Created:** April 07, 2022

**Last Modified:** June 20, 2022

**Protocol Integer ID:** 60447

**Keywords:** particulate phosphorus, intracellular phosphorus, phosphomolybdenum-ascorbic reduction, orthophosphate, oxalate reagent, adsorbed phosphorus, X-HTDC, High temperature dry combustion, estimating particulate phosphorus, measuring particulate phosphorus, particulate phosphorus from organic compound, total particulate phosphorus, particulate phosphorus from microalgae, particulate phosphorus, phosphorus recovery from several organic phosphorus compound, phosphorus recovery, tpp recovery from laboratory phytoplankton culture, intracellular phosphorus in microalgae, intracellular phosphorus, laboratory phytoplankton culture, dry combustion ash, marine phytoplankton culture, adsorbed phosphorus, hydrolysis of the ash, several organic phosphorus compound, used combustion temperature, particulate sample, combustion temperature, polyphosphate, um orthophosphate, high temperature dry combustion, particulate samples from the field, minimum sampling biomass



## Abstract

Total particulate phosphorus (TPP) is often determined using the High Temperature Dry Combustion (HTDC) method followed by hydrolysis of the ash and then molybdenum colorimetry. Here we show that a higher than traditionally-used combustion temperature, 800 °C vs. 450 - 550 °C, improves phosphorus recovery from several organic phosphorus compounds, marine phytoplankton cultures and particulate samples from the field. In aggregate these improvements to the method double the P recovery from phospholipids to 97%. TPP recovery from laboratory phytoplankton cultures and field samples increased an average of 13%, primarily due to improvements in P recovery from phospholipids, polyphosphates, and nucleic acids. We refer to this new method as the eXtra high temperature dry combustion ash/hydrol method (X-HTDC) and recommend its application for measuring particulate phosphorus from organic compounds in aquatic systems.

The working range of this assay is 1.22 to 500 uM orthophosphate. Minimum sampling biomass is 0.19 ug P/filter. In order to assess the intracellular phosphorus in microalgae, we recommend an oxalate reagent (Tovar-Sanchez 2003) to wash the microalgae collected on the filter to remove surface adsorbed phosphorus.

### Citation

P.S. Chen, T.Y. Toribara and Huber Warner. Microdetermination of Phosphorus. Anal. Chem..

<https://doi.org/10.1021/ac60119a033>

LINK

### Citation

AntonioTovar-Sanchez, Sergio A Sañudo-Wilhelmy, Manuel Garcia-Vargas, Richard S Weaver, Linda C Popels, David A Hutchins

. A trace metal clean reagent to remove surface-bound iron from marine phytoplankton. Marine Chemistry.

[https://doi.org/10.1016/S0304-4203\(03\)00054-9](https://doi.org/10.1016/S0304-4203(03)00054-9)

LINK

## Protocol materials

⊗ Magnesium sulfate anhydrous **Fisher Scientific Catalog #M65500**

⊗ 12 N Hydrochloric acid

⊗ Potassium dihydrogen orthophosphate **ACP Chemicals Catalog #P-4550**

⊗ 18M sulfuric acid

⊗ Ammonium molybdate **Merck MilliporeSigma (Sigma-Aldrich) Catalog #09878-100G**

⊗ Ascorbic acid **Merck MilliporeSigma (Sigma-Aldrich) Catalog #A5960-100G**

## Troubleshooting

## Safety warnings

⚠ Polycarbonate filter can release toxic gas and smoke during combustion. An exhaust system is required for muffle furnace while using the X-HTDC method.

## Before start

We have found that crucibles may lose their temperature resistance after acid-washing or long soaks in alkaline detergent. Crucibles tended to shatter in the oven during the initial increase in temperature from room temperature to 500 °C, even when the ramp rate was carefully controlled at 150 °C/h. We recommend not soaking crucibles in acid but instead we suggest the crucibles be filled with 0.2 M HCl and then incubated at 90 °C for 30 minutes as the acid-washing step. It is necessary to inspect the temperature resistance of newly acquired crucibles by combusting them at 500 °C for 6 h (ramp rate: 150 °C/h) after acid-washing. We found that crucibles that pass this inspection do not usually shatter when heated to 800 °C.



## Sampling

### 1 **Sampling microalgae for total particulate phosphorus (i.e. intracellular phosphorus and adsorbed phosphorus)**

- 1.1 Filter microalgae in liquid media onto polycarbonate filters, using gentle vacuum pressure (130 mmHg).

#### Equipment

##### Filter forceps

NAME

blunt end, stainless steel

TYPE

Millipore

BRAND

XX6200006P

SKU

- 1.2 Rinse samples with filtered seawater

- 1.3 Place sample filters in 2 mL Cryogenic Vials.

#### Equipment

##### Cryogenic Vials with Closures

NAME

Polypropylene, 2 mL

TYPE


Corning®

BRAND

66021-974

SKU

- 1.4 Filter blank media (without cells) through polycarbonate filter as blank.

1.5 Flash freeze filters and store at  -20 °C .

## 2 Sampling microalgae for intracellular particulate phosphorus

2.1 Filter microalgae in liquid media onto polycarbonate filters , using gentle vacuum pressure (5 inches Hg).

### Equipment

Filter forceps

NAME

blunt end, stainless steel



TYPE

Millipore

BRAND

XX6200006P

SKU

2.2 Add  5 mL oxalate reagent onto the filter, and let oxalate reagent sit in the filter funnel for  00:05:00

5m

### Protocol



NAME

Preparation of oxalate reagent


CREATED BY

Ying-Yu Hu

Preview



2.2.1 Add  50 mL MilliQ water in a 250 mL beaker.

2.2.2 Weigh  40 g NaOH and slowly pour into the beaker.

2.2.3 Use squeeze bottle to rinse the weighing boat and transfer rinse water into the same beaker.

2.2.4 Use glass rod to gently stir and fully dissolve NaOH.

**Note**

The solution is very hot and corrosive. It can cause skin burns and eye damage.

2.2.5 Carefully transfer NaOH solution into 100 mL volumetric flask by using glass rod.












2.2.6 Rinse beaker with small amount of MilliQ water three times, transfer rinse water into the flask.

2.2.7 Mix the solution by gently shaking the capped volumetric flask and top to 100 mL with MilliQ water.

2.2.8 Transfer the prepared reagent into a 250 mL PP bottle.

2.2.9 Label the bottle with SDS pictogram.



- 2.2.10 In a 1000 mL beaker with stir bar, add  600 mL MilliQ water.
- 2.2.11 Add  18.6 g EDTA,  14.7 g sodium citrate,  0.74 g KCl and  5 g NaCl into the beaker, stir until all ingredients are dissolved.  5.7
- 2.2.12  10 Molarity (M) NaOH is added dropwise to bring pH in between 6 to 7 by using a transfer pipet
- 2.2.13 Add  12.6 g oxalic acid to the solution, stir the mixture while heating.
- 2.2.14 After oxalic acid is completely dissolved, stop heating and let it cool to room temperature. A water bath filled with tap water can be used to speed up cooling.  3.3
- 2.2.15 Add  10 Molarity (M) NaOH dropwise to bring pH to  8
- 2.2.16 Top to 1 L in volumetric flask with MilliQ water.
- 2.2.17 Filter oxalate reagent by rapid flow to a 1 L PP bottle.





### Equipment

Sterile Disposable Filter Units with PES Membrane		NAME
Thermo Scientific™ Nalgene™ Rapid-Flow™		BRAND
5964520		SKU

2.2.18 Label the bottle and keep it at  Room temperature .


2.3 Drain and then rinse the sample with filtered seawater once

2.4 Place sample filters in 2 mL Cryogenic Vials.

### Equipment

Cryogenic Vials with Closures		NAME
Polypropylene, 2 mL		TYPE
Corning®		BRAND
66021-974		SKU

2.5 Filter blank media (without cells) through polycarbonate filter as blank.

2.6 Flash freeze filters and store at  -20 °C .

## X-HTDC-ing



- 3 Mark number at the bottom of each crucible with pencil, log the following information:
- (1) The number of crucible
  - (2) The code of sample in the crucible

#### Equipment

##### Porcelain crucibles

40 mL

VWR

89037-996

NAME

TYPE

BRAND

SKU

#### Equipment

##### Crucible cover

VWR

71000-146

NAME

BRAND

SKU


- 4 Transfer sample to crucible with clean filter forceps and lay filter at the bottom.

- 5 [M] 0.17 M  $\text{MgSO}_4$  reagent:

Dissolve  1.023 g  $\text{MgSO}_4$  in 50 mL MilliQ water



Magnesium sulfate anhydrous **Fisher Scientific Catalog #M65500**

- 6 Add  200  $\mu\text{L}$  [M] 0.17 M  $\text{MgSO}_4$  directly onto each sample and blank filter.

#### Note

Sing-use pipet tip to avoid cross-contamination.

- 7 Cover the crucibles and place in the oven at **90 °C** until samples are completely dry.

#### Equipment

**Forced air oven**

NAME

VWR

BRAND

89511-410

SKU

#### Note

Remove samples out of the oven as soon as they are dried. If muffle furnace is not available, keep crucibles in the vacuum desiccator.

- 8 Combust dried samples at **800 °C** for **09:00:00**

9h

#### Equipment

**Muffle furnace**

NAME

F30428C

TYPE

Thermo

BRAND


10-505-13

SKU

**Note**

Map the location of crucibles in the oven, in case pencil mark disappears under 800°C.

**Note**

Ramp rate should be controlled at <  200 °C /hour or follow the instruction provided by manufacture, otherwise the crucibles might shatter.

SP.RAT: 150/PAMPU: hour

Or

SP.RAT: 2.5/PAMPU: minute

- 9 Allow samples to gradually cool down in the muffle furnace.
- 10 Pencil mark on crucibles should be still visible, however, it can be easily removed by water. Therefore, when removing samples out of the furnace, label the lid and crucible with sharpie immediately.

**Digesting**


- 11  0.2 M HCl reagent:

In a reagent bottle, dissolve one part of  12 N HCl in 59 parts of MilliQ water

 12 N Hydrochloric acid


**Note**

Volume of HCl\_0.2M\_mL = (5\_mL) X (#Sample + #Blank)


- 12 Preheat oven to  90 °C
- 13 Add 5 mL 0.2 M HCl to each crucible.
- 14 Gently swirl the crucible.



15 Cover the crucibles and place crucibles in the muffin tin pan for easier-handling.

16 Incubate in the oven for  00:30:00

30m

17 Cool samples down to  Room temperature

18 Gently swirl the crucible and then transfer 500 uL solution to 2 mL microtube. Duplicate each sample and blank.

#### Equipment

Maxymum Recovery® Snaplock Microcentrifuge Tube

NAME

2.0 mL, Polypropylene, Clear, Nonsterile,

TYPE

Axygen®

BRAND

MCT-200-L-C


SKU

## Preparing standard working solutions



2h

19 Standard working solutions and reagents can be prepared during sample digestion.

20  $\text{KH}_2\text{PO}_4$  primary standard stock solution ( $\approx 1 \text{ mM}$ )

 Potassium dihydrogen orthophosphate **ACP Chemicals Catalog #P-4550**

20.1 Transfer about 1 g  $\text{KH}_2\text{PO}_4$  into a beaker, cover the beaker with foil

20.2 Place the beaker into an oven, dry  $\text{KH}_2\text{PO}_4$  at  110 °C for at least  02:00:00

2h



20.3 Move  $\text{KH}_2\text{PO}_4$  into a vacuum desiccator, allow  $\text{KH}_2\text{PO}_4$  to cool to room temperature

20.4 Dissolve around 0.136 g dried  $\text{KH}_2\text{PO}_4$  in 1 L milliQ water.

- Use 1 L volumetric flask
- Take notes of the actual weight of  $\text{KH}_2\text{PO}_4$  for final concentration of standard stock solution

20.5 Transfer standard stock solution into a 1 L bottle and store in the fridge.

#### Note

This stock solution lasts quite a long time, unless there is evidence for growth of algae or other extraneous biotic material.

## 21 Standard working solution

	$\text{KH}_2\text{PO}_4$	Primary (ul)	Milli Q (ul)
	S1	0	1000
	S2	5	995
	S3	10	990
	S4	20	980
	S5	50	950
	S6	100	900
	S7	150	850
	S8	200	800

22 Transfer 500 uL of each standard working solution to 2 mL microtube.

## Preparing working reagents

2h


## 23

### Note



All reagents are freshly prepared before colorimetric measurement.

24 [M] 6 N (3 M) sulfuric acid reagent:


Carefully add 1 part [M] 18 M concentrated sulfuric acid into 5 part MilliQ water

 18M sulfuric acid



25 [M] 2.5 % ammonium molybdate reagent:


Weigh  0.25 g ammonium molybdate in a Falcon tube and top to  10 g with MilliQ water.

Cap and shake until totally dissolved.

 Ammonium molybdate **Merck MilliporeSigma (Sigma-Aldrich) Catalog #09878-100G**

26 [M] 10 % ascorbic acid reagent:

Weigh  1 g ascorbic acid in a Falcon tube and top to  10 g with MilliQ water; Cap and shake until all dissolved.

 Ascorbic acid **Merck MilliporeSigma (Sigma-Aldrich) Catalog #A5960-100G**

### Note

Wrap the tube with foil if the reagent is not used right after prepared.

27 Calculate the volume of molybdate-ascorbic reagent:

Total volume of reagent\_mL = (0.5 mL) X (#standard working solution + #samples + #blanks)

28 Mix the reagents into Falcon tube:

A	B
Reagent	Part(s) as in volume



A	B
MilliQ	2
6N sulphuric acid	1
2.5% ammonium molybdate	1
10% ascorbic acid	1

## Colorimetric measurement

2h

29 Preheat incubator/shaker to  37 °C

### Equipment

**SHAKING INCUBATOR**

NAME

71L


TYPE

Corning® LSE™

BRAND

6753

SKU

30 Add  500 µL reagent to each standard, sample and blank, starting from blanks, including blank for standards and blank for samples.

### Equipment

**Finntip Stepper Tips**

NAME

5 mL

TYPE

Thermo Scientific

BRAND

9404200

SKU



### Note

Before dispensing the reagent, wipe or dab the liquid drop on the outside of the tip, avoid wiping the open tip.

31 Vortex each tube.

32 Incubate at  37 °C for  03:00:00 while shaking at 200 rpm

3h

33 Load microplate with 250 uL reactant from each tube, duplicate.

	<u>1</u>	<u>2</u>	<u>3</u>	<u>4</u>	<u>5</u>	<u>6</u>	<u>7</u>	<u>8</u>	<u>9</u>	<u>10</u>	<u>11</u>	<u>12</u>
<u>A</u>	S1	S1	Samples and sample blanks: 40 with duplicate									
<u>B</u>	S2	S2										
<u>C</u>	S3	S3										
<u>D</u>	S4	S4										
<u>E</u>	S5	S5										
<u>F</u>	S6	S6										
<u>G</u>	S7	S7										
<u>H</u>	S8	S8										

Example of loading the microplate

### Equipment

96-Well Microplates, Polystyrene, Clear,

NAME

Greiner Bio-One

BRAND

655101

SKU

34 Read plate in microplate reader

A	B
Shake duration	00:00:05
Shaking type	Continuous
Shaking force	High
Shaking speed [rpm]	600
Wavelength [nm]	820
Use transmittance	No
Pathlength correction	No
Measurement Time [ms]	100

#### Equipment

Varioskan LUX Multimode Microplate Reader NAME

Thermo Fisher BRAND

VL0L00D0 SKU

## Calculating

- 35 Subtract the average absorbance at 820 nm of the blank standard replicates from the absorbance at 820 nm of all other standard working solutions.
- 36 Subtract the average absorbance at 820 nm of the blank sample (i.e. blank filter) replicates from the absorbance at 820 nm of all other individual samples.
- 37 Prepare a standard curve by plotting the average blank-corrected 820 nm absorbance for each standard working solution versus its concentration in uM.
- 38 Use the standard curve to determine the orthophosphate concentration of each unknown sample by using its blank-corrected 820 nm absorbance.



39 (P per sample)<sub>ug</sub> = (orthophosphate)<sub>uM</sub> X (V<sub>HCl</sub>)<sub>mL</sub> X (0.001) X (30.97)

## Citations

AntonioTovar-Sanchez, Sergio A Sañudo-Wilhelmy, Manuel Garcia-Vargas, Richard S Weaver, Linda C Popels, David A Hutchins. A trace metal clean reagent to remove surface-bound iron from marine phytoplankton  
[https://doi.org/10.1016/S0304-4203\(03\)00054-9](https://doi.org/10.1016/S0304-4203(03)00054-9)

P.S. Chen, T.Y. Toribara and Huber Warner. Microdetermination of Phosphorus  
<https://doi.org/10.1021/ac60119a033>