



Jun 23, 2023

Analysis of ER Flux in Cultured Induced Neurons using Keima ER reporters



Forked from [Analysis of Lysophagic Flux in Cultured Induced Neurons using RFP-GFP-galectin3](#)



[Nature Cell Biology](#)



In 3 collections

DOI

<https://dx.doi.org/10.17504/protocols.io.j8nlko51xv5r/v1>

Melissa Hoyer¹, Harper JW¹, Vinay Eapen¹

¹Department of Cell Biology, Harvard Medical School, Boston MA 02115

Harper

ASAP



Melissa Hoyer

Harvard Medical School

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.j8nlko51xv5r/v1>



Protocol Citation: Melissa Hoyer, Harper JW, Vinay Eapen 2023. Analysis of ER Flux in Cultured Induced Neurons using Keima ER reporters. **protocols.io** <https://dx.doi.org/10.17504/protocols.io.j8nlko51xv5r/v1>

Manuscript citation:

License: This is an open access protocol distributed under the terms of the **[Creative Commons Attribution License](#)**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: June 23, 2023

Last Modified: May 31, 2024

Protocol Integer ID: 83949

Keywords: ASAPCRN, axonal er accumulation within synaptic er bouton, autophagy for individual er protein cargo, quantitative landscape of er proteome remodeling, aberrant er accumulation in axon, axonal er accumulation, molecular inventory of er, er during cell state transition, individual er protein cargo, synaptic er bouton, keima er reporters the endoplasmic reticulum, er remodeling by autophagy, endoplasmic reticulum, extensive er remodeling during differentiation, specific autophagic capture of er, deficient neurons in vivo, phagy receptor fam134b, deficient neuron, er proteome remodeling, embedded receptor, er to degradative autophagy machinery, neuron, underlying er remodeling, aberrant er accumulation, refined tubular er network, phagy receptor mutant, cultured induced neuron, phagy receptor, selective autophagy, shaping protein, vast proteomic landscape, er proteome, protein, degradative autophagy machinery, receptor, monitoring extensive er remodeling, tractable induced neuron, lumenal protein, autophagy, dis

Funders Acknowledgements:

ASAP

Grant ID: ASAP-000282

Disclaimer

DISCLAIMER – FOR INFORMATIONAL PURPOSES ONLY; USE AT YOUR OWN RISK

The protocol content here is for informational purposes only and does not constitute legal, medical, clinical, or safety advice, or otherwise; content added to **[protocols.io](#)** is not peer reviewed and may not have undergone a formal approval of any kind. Information presented in this protocol should not substitute for independent professional judgment, advice, diagnosis, or treatment. Any action you take or refrain from taking using or relying upon the information presented here is strictly at your own risk. You agree that neither the Company nor any of the authors, contributors, administrators, or anyone else associated with **[protocols.io](#)**, can be held responsible for your use of the information contained in or linked to this protocol or any of our Sites/Apps and Services.



Abstract

The endoplasmic reticulum (ER) has a vast proteomic landscape to preform many diverse functions including protein and lipid synthesis, calcium ion flux, and inter-organelle communication. The ER proteome is remodeled in part through membrane-embedded receptors linking ER to degradative autophagy machinery (selective ER-phagy). A refined tubular ER network is formed in neurons within highly polarized dendrites and axons. Autophagy-deficient neurons *in vivo* display axonal ER accumulation within synaptic ER boutons, and the ER-phagy receptor FAM134B has been genetically linked with human sensory and autonomic neuropathy. However, mechanisms and receptor selectivity underlying ER remodeling by autophagy in neurons is limited. Here, we combine a genetically tractable induced neuron (iNeuron) system for monitoring extensive ER remodeling during differentiation with proteomic and computational tools to create a quantitative landscape of ER proteome remodeling via selective autophagy. Through analysis of single and combinatorial ER-phagy receptor mutants, we delineate the extent to which each receptor contributes to both magnitude and selectivity of ER clearance via autophagy for individual ER protein cargos. We define specific subsets of ER curvature-shaping proteins or luminal proteins as preferred clients for distinct receptors. Using spatial sensors and flux reporters, we demonstrate receptor-specific autophagic capture of ER in axons, which correlates with aberrant ER accumulation in axons of ER-phagy receptor or autophagy-deficient neurons. This molecular inventory of ER proteome remodeling and versatile genetic toolkit provides a quantitative framework for understanding contributions of individual ER-phagy receptors for reshaping ER during cell state transitions.



Materials

	A	B	C
	REAGENT or RESOURCE	SOURCE	IDENTIFIER
	Chemicals		
	Dulbecco's MEM (DMEM), F12	Thermo	11320033
	Phosphate Buffered Saline 1X	Corning	21-031-CV
	E8 components	See protocol dx.doi.org/10.17504/protocols.io.bsacnaaw	
	ND1 and ND2 components	See protocol dx.doi.org/10.17504/protocols.io.br9em93e	
	Hygromycin B	Life Technologies	10687-010
	Bafilomycin A1	Cayman Chemical	88899-55-2
	Sar405 Selective ATP-competitive inhibitor of Vps34	Apexbio	A8883
	pAC150 Keima-RAMP4		Addgene 201929
	pAC150 Keima-REEP5		Addgene 201929
	pCMV-hyPBase hyperactive piggyBac vector	Yusa et al 2011	Available upon request at the Sanger Institute Archives
	Software		



	A	B	C
	FlowJo Software	" Flowjo, v10.7 "	https://www.flowjo.com
	Hardware		
	ThermoFisher Neon Electroporator & Kit		
	AttuneNxt	Thermo Fisher Scientific,	Cat#A28993

Electroporation of PB vectors. Use ThermoFisher kit and ThermoFisher Neon Electroporator to electroporate ES cells with PB vector and PB helper vector.

- 1 Add 10ml buffer R to a sterile 1.5ml tube. Add 0.5µg of pAC150 ER Keima vector and 0.5µg of pCMV-hyPBBase hyperactive piggyBac vector. Pipet up and down to mix. Let it sit at RT for 10min. This is enough for 2 transfections (== one 6 well).   
- 2 Individualize cells with Accutase. Neutralize Accutase with 5 times volume E8 with Rock inhibitor. 
- 3 Count cells. You will need 2×10^5 for each transfection.
- 4 Spin down cells. Let it sit for a while so all the residue media can go down to the bottom of the tube. If the residue media is too much, take it out with a P200 pipet. 
- 5 Resuspend cells to a concentration of 2×10^5 per 5 ml (ie 4×10^7 per ml) using buffer R. You don't have to take all the residue media off but you will need to take into account the volume of residue media so you are not too much off. 
- 6 Prepare a 6 well matrigel coated plate. Add 2mL of E8+ rock inhibitor (1:1000) to the wells you will use. Two transfections go into one well. 
- 7 Wipe the Neon pipet station with EtOH and place it inside the hood.
- 8 Add 3ml of electrolytic buffer (buffer E) to the neon tube. Place the tube inside the station. You should feel a click before the tube is securely seated in the station.
- 9 Use program 13 from the optimization tab for electroporation parameter (Voltage: 1100. Pulse width: 20 Pulse number: 2). Program 9 should also work.
- 10 When everything is ready, mix 10-11ml of resuspended cells with the plasmid containing R buffer. The final volume should be in the range of 21-22ml. 
- 11 Take up a neon tip, pipet 10ml of the cell protein mix and electroporate with program 13. It is important to pipet slowly to avoid air bubble formation. It is also important to insert the pipet slowly into the station, especially during the end of the insertion when you will feel a click. Help the pipet down slowly during the clicking so there is no sudden movement of the tip, which might create tiny air bubbles. 



- 12 If you see air bubble in the tip, take it out, push everything out of the tip and repipet the mixture. 
- 13 If you see sparking during the electroporation, your efficiency will reduce significantly. 
- 14 Once electroporation is complete, push everything into one well of a 6 well plate. Do not pipet up and down with Neon tip. 
- 15 Repeat the same procedure with the same tip and the left over cell mixture. Place the second electroporated mixture into the same well. 
- 16 Disperse cells evenly in the well and place cells in a low O2 incubator. 17. Put electroporated cells into low oxygen incubator for 2 days 
- 17 Select cells with 50mg/mL hygromycin B 4 days post-electroporation. Grow cells in selection medium for 7-10 days until there is no longer any cell death and every cell has integration of the piggyBac Keima-ER construct. Also can sort for the same level of Keima in each ESC.

Differentiation of Stable Cell ES H9 line expressing Keima-ER reporter to induced neurons (iN)

- 18 Differentiation to induced neurons (iN) is done by following the protocol "Neural differentiation of AAVS1-TREG3-NGN2 pluripotent stem cells ([dx.doi.org/10.17504/protocols.io.br9em93e](https://doi.org/10.17504/protocols.io.br9em93e))

Analysis of Keima-ER via confocal microscopy

- 19 To help get spread out iNeurons for imaging, when cells get 90 percent confluent at any point in the day 5-7 range, cells are plated into onto 35 mm-glass bottom dishes (No. 1.5, 14 mm glass diameter, MatTek) at a low confluency (approximately 1×10^5 cells per 2mL dish). 
- 20 At day 12 or any day during differentiation, image the iNeurons 
- 21 5mL of fresh ND2/LLoMe is added to a final concentration of 500 mM of LLoMe in 2mL. Cells are treated for 1h.  



- 22 At day 12 or any day during differentiation, image the iNeurons. Cells are imaged at RT using a Yokogawa CSU-X1 spinning disk confocal on a Nikon Ti-E inverted microscope at the Nikon Imaging Center in Harvard Medical School. Nikon Perfect Focus System was used to maintain cell focus over time. The microscope is equipped with a Nikon Plan Apo 40x/1.30 N.A or 100x/1.40 N.A objective lens. 488nm (100mW) and 561nm (100mW) laser lines are controlled by AOTF. All images are collected with a Hamamatsu ORCA-ER cooled CCD camera (6.45 μm^2 photodiode) with MetaMorph image acquisition software.
- 23 Z series are displayed as maximum z-projections and saved using Fiji software.
- 24 Each channel z series are brightness and contrast adjusted equally and then converted to RGB for publication using FIJI software.

Analysis of Keima-ER via flow cytometry

- 25 iNeurons can be treated throughout the differentiation with a vps34i (1 microM in ND2 each day you add media) to block the ER flux. On the day of flow cytometry 4hr to collection for analysis treat one set of each condition with BAFA (100nM in ND2).
- 26 Detach cells from plate (1 well of 6-well dish of iN is sufficient) using accutase, add ND2 media and spin down, and resuspend carefully with 1XPBS so they do not clump. If required, filtered through a cell strainer cap tube. Place in vessel of choice for cytometry, for example tubes or 96 well plates.
- 27 Use dual-excitation (440nm for pH7 and 561nm for pH3) and collect in 620 nm range. Analyze at least 10,000 single, healthy cells. Calculate of acidic:neutral Keima ratio on a per-cell basis in FlowJo Software.
Use BafA treated as normalization for sample.