



Jul 21, 2025

# 🌐 Analysis of sphingosine, sphinganine and glucosylsphingosine from cells, animal tissues and plasma

DOI

[dx.doi.org/10.17504/protocols.io.bp2l6zkz1gqe/v1](https://dx.doi.org/10.17504/protocols.io.bp2l6zkz1gqe/v1)

Danielle te Vrugte<sup>1</sup>, Yuzhe Weng<sup>1</sup>, Mylene Huebecker<sup>1</sup>, Kerri-Lee Wallom<sup>1</sup>, Frances M Platt<sup>1,2,3,4</sup>

<sup>1</sup>Department of Pharmacology, University of Oxford, Mansfield Road, Oxford, OX1 3QT, United Kingdom;

<sup>2</sup>Aligning Science Across Parkinson's (ASAP) Collaborative Research Network, Chevy Chase, MD, 20815;

<sup>3</sup>Royal Society Wolfson Research Merit Award Holder; <sup>4</sup>Wellcome Trust Investigator in Science



Danielle te Vrugte

University of Oxford

## Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.bp2l6zkz1gqe/v1>

**Protocol Citation:** Danielle te Vrugte, Yuzhe Weng, Mylene Huebecker, Kerri-Lee Wallom, Frances M Platt 2025. Analysis of sphingosine, sphinganine and glucosylsphingosine from cells, animal tissues and plasma. **protocols.io**

<https://dx.doi.org/10.17504/protocols.io.bp2l6zkz1gqe/v1>

**License:** This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

**Protocol status:** Working

**We use this protocol and it's working**

**Created:** July 02, 2025

**Last Modified:** July 21, 2025

**Protocol Integer ID:** 221572

**Keywords:** sphingolipid metabolism, glucosylsphingosine from cell, analysis of sphingosine, glucosylsphingosine in cell line, extracted sphingosine, complex sphingolipid, glucosylsphingosine comprise, glucosylsphingosine, labelled c20 sphingosine, complex sphingolipids like ceramide, plasma sphingosine, sphingosine, c20 sphingosine, labelled c20 glucosylsphingosine standard, c20 glucosylsphingosine standard, glucosylceramide, sphinganine, glucose, sphingoid base, phthaldialdehyde, crucial molecule, cell line, cell, molar quantities of these molecule, molecule, fluorescent derivative

**Funders Acknowledgements:**

**The Michael J. Fox Foundation for Parkinson's Research (MJFF) and the Aligning Science Across Parkinson's (ASAP) Initiative**  
Grant ID: ASAP-000478

**Welcome Investigator Award**  
Grant ID: 202834/Z/16/Z

## Abstract

Sphingosine (SphO), sphinganine (SphA) and glucosylsphingosine (GlcSph) are crucial molecules involved in sphingolipid metabolism. Sphingosine and sphinganine are long-chain bases that serve as building blocks for more complex sphingolipids like ceramides and glucosylceramides (GlcCer). Glucosylsphingosine comprises of a sphingoid base with glucose attached.

Here we describe a method which can measure sphingosine, sphinganine and glucosylsphingosine in cell lines, animal tissue and plasma. The extracted sphingosine, sphinganine and glucosylsphingosine are reacted with Phthaldialdehyde (OPA) reagent, in the presence of the thiol compound mercaptoethanol, to form a fluorescent derivative which can be measured by reverse phase HPLC. With the inclusion of OPA-labelled C20 sphingosine or OPA-labelled C20 glucosylsphingosine standards, it is possible to obtain molar quantities of these molecules.

## Guidelines

This protocol requires the use of some hazardous materials. As such, users must be appropriately trained and hazardous materials stored, used and disposed of in accordance with your institution's health and safety policies, and approved laboratory policies, risk assessments and codes of practice.

## Materials

### **Acetone for HPLC, ≥99.8%**

34850-M Sigma-Aldrich

<https://www.sigmaaldrich.com/GB/en/product/sial/34850m>

### **Acetonitrile for HPLC, gradient grade, ≥99.9%**

34851 Sigma-Aldrich

<https://www.sigmaaldrich.com/GB/en/product/sigald/34851>

### **Bicinchonic Acid solution**

B9643-1L Sigma-Aldrich

[https://www.sigmaaldrich.com/GB/en/product/sigma/b9643?srltid=AfmBOog6n40M013YTRhV9rEDobVTfCL49\\_Z9ZPBpcNM3KWP8Pu7Zcnx7](https://www.sigmaaldrich.com/GB/en/product/sigma/b9643?srltid=AfmBOog6n40M013YTRhV9rEDobVTfCL49_Z9ZPBpcNM3KWP8Pu7Zcnx7)

### **Boric acid**

194810 MP Biomedicals, LLC

[https://www.mpbio.com/uk/boric-acid?queryID=27e466009c404622801dbc590258b64b&objectID=52001&indexName=MpBio\\_Production\\_Store\\_uk\\_products](https://www.mpbio.com/uk/boric-acid?queryID=27e466009c404622801dbc590258b64b&objectID=52001&indexName=MpBio_Production_Store_uk_products)

### **Chloroform Stabilized with Amylene, for HPLC, ≥99.8%**

C/4966/17 Fisher Scientific

<https://www.fishersci.co.uk/shop/products/chloroform-stabilized-amylene-hplc-fisher-chemical/10615492>

### **Chromolith Performance RP-18e 100-4.6 HPLC column (Merck, Darmstadt, Germany)**

1.02129.0001 VWR International Ltd.

<https://www.avantorsciences.com/uk/en/catalog-number/1.02129.0001>

### **Copper(II) sulfate solution**

C2284 Sigma-Aldrich

<https://www.sigmaaldrich.com/GB/en/product/sigald/c2284>

### **Ethanol, ≥99.8%**

32221-M Sigma-Aldrich

<https://www.sigmaaldrich.com/GB/en/product/sial/32221m>

### **C18-glucosylsphingosine standard (1 mM in EtOH, store at -80°C)**

860535 Glucosyl(β) Sphingosine (d18:1), Avanti Polar Lipids

<https://avantiresearch.com/product/860535>

### **C20-glucosylsphingosine standard (1 mM in EtOH, store at -80°C)**

860438 Glucosyl (β) Sphingosine (d20:1), Avanti Polar Lipids



<https://avantiresearch.com/product/860438>

**2-Mercaptoethanol, ≥99.0%**

63689 Sigma-Aldrich

[https://www.sigmaaldrich.com/GB/en/product/sigma/63689?  
srltid=AfmBOoqw0ieNmCN9Guq7aBo2lcVdf0TOeeGQL0V-AdRfGXVe2Gh1iCOF](https://www.sigmaaldrich.com/GB/en/product/sigma/63689?srltid=AfmBOoqw0ieNmCN9Guq7aBo2lcVdf0TOeeGQL0V-AdRfGXVe2Gh1iCOF)

**Methanol for HPLC grade, ≥99.9%**

34860 Sigma-Aldrich

<https://www.sigmaaldrich.com/GB/en/product/sigald/34860>

**OPA (Phthaldialdehyde), ≥99.0% (HPLC)**

79760 Sigma-Aldrich

<https://www.sigmaaldrich.com/GB/en/product/sigma/79760>

**Sodium chloride, AR**

S/3160/60 Fisher Scientific

<https://www.fishersci.co.uk/shop/products/sodium-chloride-certified-ar-analysis-meets-analytical-specification-ph-eur-fisher-chemical/10428420>

**Sodium hydroxide**

06203 Sigma-Aldrich

[https://www.sigmaaldrich.com/GB/en/product/sigald/06203?  
srltid=AfmBOorZjlns7uY\\_IQmW22fbzx3HWex8Bj9AgA7MtIEDMgAdZX2UH58V](https://www.sigmaaldrich.com/GB/en/product/sigald/06203?srltid=AfmBOorZjlns7uY_IQmW22fbzx3HWex8Bj9AgA7MtIEDMgAdZX2UH58V)

**C18-sphinganine standard (1 mM in EtOH, store at -80°C)**

860498 Sphinganine (d18:0), Avanti Polar Lipids

<https://avantiresearch.com/product/860498>

**C18-sphingosine standard (1 mM in EtOH, store at -80°C)**

860490 Sphingosine (d18:1), Avanti Polar Lipids

<https://avantiresearch.com/product/860490>

**C20-sphingosine standard (1 mM in EtOH, store at -80°C)**

860660 Sphingosine (d20:1), Avanti Polar Lipids

<https://avantiresearch.com/product/860660>

**Tris base**

sc-3715A Santa Cruz

[https://www.scbt.com/p/tris-base-77-86-1?  
srltid=AfmBOorCdFXWCa0RIGTvscM5u57kFWQMmJ\\_4FX1ABIIZcZ2giYjad8Zp](https://www.scbt.com/p/tris-base-77-86-1?srltid=AfmBOorCdFXWCa0RIGTvscM5u57kFWQMmJ_4FX1ABIIZcZ2giYjad8Zp)

**ISOLUTE® NH2 100 mg/1 mL**



#470-0010-A Biotage

<https://www.biotage.com/isolute-nh2>

### Equipment:

Pipettes

Homogeniser

Vortex mixer

Ultrasonic water bath

Centrifuge

Sample concentrator (Heatblock with a nitrogen supply)

RP-HPLC

### Troubleshooting

### Safety warnings

- ⚠ This protocol requires the use of some hazardous solvents, reagents and chemicals. Refer to the Safety Data Sheets (SDS) provided by supplier and applicable Control of Substances Harmful to Health (COSHH). The correct personal protective equipment must be worn, and incidents reported in line with your institution's policy and procedures.

### Before start

Check that you have the required reagents, solvents, chemicals, equipment and PPE.

## Cells homogenisation

- 1 Lyse and homogenise cell pellet ( $\geq 1 \times 10^6$  cells) in  200  $\mu\text{L}$  ddH<sub>2</sub>O through three cycles of freeze-thawing with vortexing after each cycle.
- 2 Measure protein concentration using the Bicinchonic acid assay (BCA) method in triplicate (repeat if triplicate values are not consistent).
- 3 Using determined protein concentrations, prepare samples as  200  $\mu\text{g}$  protein in  100  $\mu\text{L}$  ddH<sub>2</sub>O in a 2 ml screw-cap tube on ice.



## Animal tissue homogenisation

- 4 Homogenise tissue in ddH<sub>2</sub>O yielding tissue concentration of  25 mg tissue (wet weight) per  1 mL ddH<sub>2</sub>O.
- 5 Measure protein concentration using the Bicinchonic acid assay (BCA) method in triplicate (repeat if triplicate values are not in consistent).
- 6 Using determined protein concentrations, prepare samples as  200  $\mu\text{g}$  -  500  $\mu\text{g}$  protein in  100  $\mu\text{L}$  ddH<sub>2</sub>O in a 2 ml screw-cap tube on ice.



## Plasma

- 7 Use  50  $\mu\text{L}$  plasma and add  50  $\mu\text{L}$  ddH<sub>2</sub>O in a 2 ml screw-cap tube on ice.



## Extraction of Sph and GlcSph from cells, tissue and plasma

35m

### 8 Example of the experiment layout:

|  | Sample/Standard name | ddH <sub>2</sub> O | Sample           | Standard Master Mix      | Standard |
|--|----------------------|--------------------|------------------|--------------------------|----------|
|  | Sample 1             | 50 $\mu\text{l}$   | 50 $\mu\text{l}$ | 500 $\mu\text{l}$ C20 MM | -        |
|  | Sample 2             | 50 $\mu\text{l}$   | 50 $\mu\text{l}$ | 500 $\mu\text{l}$ C20 MM | -        |
|  | Sample 3             | 50 $\mu\text{l}$   | 50 $\mu\text{l}$ | 500 $\mu\text{l}$ C20 MM | -        |
|  | Sample 4             | 50 $\mu\text{l}$   | 50 $\mu\text{l}$ | 500 $\mu\text{l}$ C20 MM | -        |

|  | Sample/Standard name | ddH <sub>2</sub> O | Sample | Standard Master Mix                  | Standard               |
|--|----------------------|--------------------|--------|--------------------------------------|------------------------|
|  | Sample 5, etc.       | 50 µl              | 50 µl  | 500 µl C20 MM                        | -                      |
|  | Blank                | 100 µl             | -      | 500 µl CHCl <sub>3</sub> :MeOH (1:2) | -                      |
|  | C18 SphO             | 100 µl             | -      | 500 µl C20 SphO mix                  | 3 µl 0.1 mM C18 SphO   |
|  | C18 SphA             | 100 µl             | -      | 500 µl C20 SphO mix                  | 3 µl 0.1 mM C18 SphA   |
|  | C20 SphO             | 100 µl             | -      | 500 µl C20 SphO mix                  | -                      |
|  | C18 GlcSph           | 100 µl             | -      | 500 µl C20 GlcSph mix                | 3 µl 0.1 mM C18 GlcSph |
|  | C20 GlcSph           | 100 µl             | -      | 500 µl C20 GlcSph mix                | -                      |

See steps 9 and 10 below on how to make up the standards and standard master mixes.

- 9 5 standards and a blank are needed:
- 1 mM C18 sphingosine (C18 SphO),
  - 1 mM C18 sphinganine (C18 SphA),
  - 1 mM C20 sphingosine (C20 Sph),
  - 1 mM C18 glucosylsphingosine (C18 GlcSph),
  - 1 mM C20 glucosphingosine (C20 GlcSph),
  - Blank (with only CHCl<sub>3</sub>:MeOH (1:2)).

9.1 Make 0.1 mM C18 SphO, 0.1 mM C18 SphA and 0.1 mM C18 GlcSph standard solutions:

50 µL stock (1 mM) + 0.5 mL EtOH.

10 Make the different C20 (master) mixes as described below in Steps 10.1 to 10.3:

10.1 For the samples make a C20 master mix (**C20 MM**) in a 15 ml tube:

6 µL 1 mM C20 Sph and 6 µL 1 mM C20 GlcSph in 10 mL CH<sub>3</sub>Cl:MeOH (1:2) (need 500 µL MM C20 per sample, scale up if necessary).

10.2 For the Sph standards make a C20 SphO mix (**C20 SphO mix**) in a 15 ml tube:

3 µL 1 mM C20 Sph in 5 mL CHCl<sub>3</sub>:MeOH (1:2) (need 500 µL MM C20 SphO per Sph standard).

10.3 For the GlcSph standards make a C20 GlcSph mix (**C20 GlcSph mix**) in a 15 ml tube:



- 3  $\mu\text{L}$  1 mM C20 GlcSph in 5 mL  $\text{CHCl}_3\text{:MeOH}$  (1:2) (need 500  $\mu\text{L}$  MM C20 GlcSph per GlcSph standard).
- 11 For the blank and the standards put 100  $\mu\text{L}$  ddH<sub>2</sub>O in a 2 ml screw-cap tube.
- 12 Add 500  $\mu\text{L}$  C20 MM to the sample tubes,   
 500  $\mu\text{L}$  C20 SphO mix to the C18 SphO, C18 SphA and C20 SphO standard tubes,  
 500  $\mu\text{L}$  C20 GlcSph mix to the C18 GlcSph and C20 GlcSph standard tubes and just  
 500  $\mu\text{L}$  pure  $\text{CHCl}_3\text{: MeOH}$  (1:2) into the Blank tube.
- 13 Add 3  $\mu\text{L}$  of 0.1 mM C18 standards according to the standard tubes labels (see experiment layout).
- 14 Sonicate all tubes for 00:10:00 at Room temperature . 10m
- 15 Add 500  $\mu\text{L}$  NaCl (1 M).
- 16 Add 500  $\mu\text{L}$   $\text{CHCl}_3$ .
- 17 Add 100  $\mu\text{L}$  NaOH (3 M).
- 18 Vortex and leave for 00:15:00 at Room temperature (vortex every 5 min). 15m
- 19 Centrifuge at 16.000 x g for 00:10:00 at Room temperature . 10m
- 20 Transfer the lower, organic phase very carefully to a new 1.5 ml Eppendorf tube (**keep!**), while leaving the upper phase and protein interphase (which can be discarded).

## Purification of Sph and GlcSph



- 21 Condition the right number of ISOLUTE ®NH2 100 mg/1 mL columns with 2 x  1 mL 
- 22 Apply the extracted organic phase to the columns, let drip through gravity flow, discard flow-through. 
- 23 **Elute first sphingosine and sphinganine, then glucosylsphingosine from the same column:**
- 23.1 To elute sphingosine and sphinganine: Apply  900 µL acetone to the column and collect the eluent in a new 1.5 ml screw-cap tube. After elution, use syringe to push to empty the column. 
- 23.2 To elute glucosylsphingosine: Apply  900 µL acetone:EtOH (6:1) to the column and collect the eluent in a new 1.5 ml screw-cap tube. After elution, use syringe to push to empty the column. 
- 24 Dry down the solvents under a slow stream of N<sub>2</sub> gas.
- 25 The samples are now ready for OPA-labelling (or they can be stored at -20°C for OPA-labelling next day). 

## OPA-labelling of Sph and GlcSph

22m

- 26 Solutions: 
- 3% boric acid, pH 10.5:**  
dissolve  6 g boric acid in  180 mL H<sub>2</sub>O  
adjust pH with 10 M KOH to 10.5  
fill up to  200 mL (Filter sterilize).
- OPA-labelling solution:** weigh 1 mg – 5 mg OPA in 15 ml tube covered with aluminium foil.
- |                                     |                                              |
|-------------------------------------|----------------------------------------------|
| Dissolve OPA in EtOH (20 µl/mg OPA) | (e.g. 1 mg OPA: add 20 µl EtOH),             |
| add 2-mercaptoethanol (1 µl/mg OPA) | (e.g. 1 mg OPA: add 1 µl 2-mercaptoethanol), |
| dilute 1:2000 with 3% boric acid    | (e.g. 1 mg OPA: add 2 ml 3% boric acid).     |
- 27 Prepare the OPA-labelling solution just before labelling (light-sensitive, instable).



- 28 Resuspend the lipids in  50  $\mu\text{L}$  pre-warmed EtOH (  37  $^{\circ}\text{C}$  water bath/heat block), vortex gently. 
- 29 Add  50  $\mu\text{L}$  fresh OPA-labelling solution to each sample and vortex. 
- 30 Incubate at  Room temperature in the dark for  00:20:00 . Vortex again after 10 min.   
 
- 31 Add  100  $\mu\text{L}$  MeOH:5 mM Tris pH 7 (9:1) and vortex. 
- 32 Centrifuge at  3000 rpm for  00:02:00 at  Room temperature .   

- 33 Transfer  150  $\mu\text{L}$  of the supernatant (carefully pipette from the top) to a HPLC vial with plastic insert. 
- 34 For RP-HPLC inject  20  $\mu\text{L}$  -  50  $\mu\text{L}$  . 

## RP-HPLC

30m

- 35 The RP-HPLC system consists of a VWR Hitachi Elite LaChrom HPLC system with a L-2485 fluorescence detector set at Ex  $\lambda$ 340nm and Em  $\lambda$ 455nm.
- 36 The solid phase used is a Chromolith Performance RP-18e 100-4.6 HPLC column (Merck, Darmstadt, Germany).
- 37 The chromatographic flow rate is 2.0 mL/min, and total run time per sample is  00:30:00 . 

|  | Time | 100% MeOH | 100% ddH2O | 100% acetonitrile | 80% ddH2O/20% acetonitrile |
|--|------|-----------|------------|-------------------|----------------------------|
|  | 0.0  | 0.0       | 0.0        | 64.0              | 36.0                       |
|  | 2.0  | 0.0       | 0.0        | 64.0              | 36.0                       |
|  | 7.0  | 0.0       | 0.0        | 88.0              | 12.0                       |

|  | Time | 100% MeOH | 100% ddH <sub>2</sub> O | 100% acetonitrile | 80% ddH <sub>2</sub> O/20% acetonitrile |
|--|------|-----------|-------------------------|-------------------|-----------------------------------------|
|  | 8.0  | 5.0       | 0.0                     | 89.0              | 6.0                                     |
|  | 22.0 | 100.0     | 0.0                     | 0.0               | 0.0                                     |
|  | 24.0 | 100.0     | 0.0                     | 0.0               | 0.0                                     |
|  | 26.0 | 0.0       | 0.0                     | 64.0              | 36.0                                    |
|  | 30.0 | 0.0       | 0.0                     | 64.0              | 36.0                                    |

**Gradient conditions for reverse phase HPLC.** All chromatography was controlled and data were collected and processed using Waters Empower software.

- 38 The mobile phases are MeOH, acetonitrile and 80% ddH<sub>2</sub>O/20% acetonitrile mix.
- 39 The gradient conditions for the RP-HPLC are as described in the table above.
- 40 Identification of Sph-species: Individual sphingosine species are identified by their retention time in comparison to standards.
- 41 Approximate Retention Times:  
C18 GlcSph ~ 3.80 min  
C20 GlcSph ~ 6.06 min  
C18 SphO ~ 6.78 min  
C18 SphA ~ 8.00 min  
C20 SphO ~ 9.18 min
- 42 Quantification: Sph species are quantified by comparison of integrated peak areas with a known amount of spiked OPA-labelled C20 sphingosine standard (Avanti Polar Lipids) or OPA-labelled C20 glucosylsphingosine standard (Avanti Polar Lipids) in each sample.
- 43 Results should be normalized to protein content or per ml of plasma.

