

Feb 20, 2025

Ancient DNA Extraction from Bones and Teeth

DOI

<https://dx.doi.org/10.17504/protocols.io.x54v9ronpv3e/v1>

Michael V. Westbury¹, Alba Rey-Iglesia¹, Andrea A. Cabrera¹, Vanssy Li¹, Deon de Jager¹, Eline D. Lorenzen¹

¹Globe Institute, University of Copenhagen

Lorenzen Lab



Vanssy Li

Copenhagen university

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.x54v9ronpv3e/v1>

Protocol Citation: Michael V. Westbury, Alba Rey-Iglesia, Andrea A. Cabrera, Vanssy Li, Deon de Jager, Eline D. Lorenzen 2025. Ancient DNA Extraction from Bones and Teeth. [protocols.io https://dx.doi.org/10.17504/protocols.io.x54v9ronpv3e/v1](https://dx.doi.org/10.17504/protocols.io.x54v9ronpv3e/v1)

License: This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited



Protocol status: Working

We use this protocol and it's working

Created: February 06, 2025

Last Modified: February 20, 2025

Protocol Integer ID: 119694

Keywords: ancient dna, teeth ancient dna extraction protocol, ancient dna extraction from bone, extracting ancient dna, ancient dna from bone, ancient dna extraction, final dna extract volume, spin column protocol, using monarch spin column, extraction, monarch spin column, spin column, dna, tooth powder, bone, extraction protocol

Abstract

Ancient DNA Extraction Protocol (spin column-based)

- This protocol outlines the procedure for extracting ancient DNA from bone or tooth powder using Monarch spin columns. The final DNA extract volume will be 35 μ L.

Based on:

- The spin column protocol is based on Dabney et al. 2013 (DOI: [10.1073/pnas.1314445110](https://doi.org/10.1073/pnas.1314445110))
- The modified binding buffer is from Allentoft et al. 2015 extraction (DOI: [10.1038/nature14507](https://doi.org/10.1038/nature14507))

Protocol initially written up by Michael V. Westbury on 20-08-2018.

Annotations added by Andrea A. Cabrera and Alba Rey-Iglesia on 08-11-2019.

The protocols.io version was created by Wenxi Li and Deon de Jager.

Protocol materials

- ⊗ Buffer EB Qiagen Catalog #19086
- ⊗ Tween 20 Merck MilliporeSigma (Sigma-Aldrich) Catalog #P1379-500ml
- ⊗ 0.5M EDTA Merck MilliporeSigma (Sigma-Aldrich) Catalog #03690
- ⊗ Proteinase K, recombinant, PCR Grade Merck MilliporeSigma (Sigma-Aldrich) Catalog #RPROTKSOL-RO
- ⊗ Buffer PB Qiagen Catalog #19066
- ⊗ 5M Sodium acetate Merck MilliporeSigma (Sigma-Aldrich)
- ⊗ 5M Sodium chloride Fisher Scientific Catalog #10609823
- ⊗ Sodium acetate buffer solution 3M pH 5.2 Merck MilliporeSigma (Sigma-Aldrich) Catalog #S7899



Before start

- Ensure that ~50 mg of bone or tooth powder is prepared in advance (DO NOT treat with UV).
- UV-treat aliquots (PB buffer, PE buffer, EB buffer, ddH₂O) for 10 minutes before use.
- Do not UV Eppendorf Tubes/PCR TUBES
- PE buffer: Add Molecular Grade Ethanol as indicated on the bottle before use.
- Use Qiagen EB buffer directly if preferred.
- Clean all surfaces with 5% bleach solution followed by 70% ethanol before and after use.
- Clean all equipment with 70% ethanol before and after use.



Buffer Preparation

1 Digestion Buffer

Choose the appropriate option based on **[ProK]**.

Prepare fresh for each use.

⊗ 0.5M EDTA Merck MilliporeSigma (Sigma-Aldrich) Catalog #03690

⊗ Proteinase K, recombinant, PCR Grade Merck MilliporeSigma (Sigma-Aldrich) Catalog #RPROTKSOL-RO



Option 1: Use 10mg/mL ProK

	Reagent	Volume per sample (μL)	Final Concentration	Master Mix (14 samples)
	0.5M EDTA	900	0.45	12,600
	10 mg/mL ProK	25	0.25	350
	ddH ₂ O	75	N/A	1,050
	Total	1,000		14,000

Option 2: Use 19mg/mL ProK

	Reagent	Volume per sample (μL)	Final Concentration	Master Mix (14 samples)
	0.5M EDTA	900	0.45	12,600
	19 mg/mL ProK	13.15	0.25	184.1
	ddH ₂ O	86.85	N/A	1,215.9
	Total	1,000		14,000

2 Modified Binding Buffer

Prepare in bulk

Option 1 (with 5M NaOAc)

- 500 mL ⊗ Buffer PB Qiagen Catalog #19066



- **9 mL** 5M Sodium acetate **Merck MilliporeSigma (Sigma-Aldrich)**
- **1.25 mL** 5M Sodium chloride **Fisher Scientific Catalog #10609823**
- Adjust to **pH 4-5** (~2 mL Hydrochloric acid)

Option 2 (with 3M NaOAc, pH 5.2)

- **500 mL** Buffer PB **Qiagen Catalog #19066**
- **15 mL**
 Sodium acetate buffer solution 3M pH 5.2 **Merck MilliporeSigma (Sigma-Aldrich) Catalog #S7899**
- **2.5 mL** 5M Sodium chloride **Fisher Scientific Catalog #10609823**
- No pH adjustment necessary as the 3M sodium acetate is pH buffered at pH5.2.

Note

- The NaOAc and NaCl can be added directly to the 500 mL of Buffer PB. Make sure to label it as **Modified Buffer PB**.
- Make individual aliquots in 50 mL Falcon tubes.

3 PE Buffer

Remember to add Molecular Grade Ethanol to PE when opening a new bottle. See instructions on the PE buffer bottle.

4 EB Buffer

Appropriate for when using LoBind Eppendorf tubes.

- Buffer EB **Qiagen Catalog #19086** . Use as provided.

4.1 EBT Buffer *optional

Appropriate for when NOT using LoBind Eppendorf tubes.

- Buffer EB **Qiagen Catalog #19086**
- Tween 20 **Merck MilliporeSigma (Sigma-Aldrich) Catalog #P1379-500ml** at a final concentration of 0.05% (v/v)

*

Day 1 - Sample digestion

19h 2m

5 Optional: Pre-digestion

30m

*



- 5.1 In a 2 mL LoBind Eppendorf tube:
- Add  50 mg  Bone Powder to  1 mL **extraction buffer**.
 - Resuspend the bone powder by **vortexing**.
- 5.2 Incubate at  50 °C with **constant rotation** for  00:30:00 (use Parafilm to seal the lids).



30m



6 Digestion



Note

- The microcentrifuge tubes are the LoBind tubes.
- Use a rotation machine to keep the samples well mixed during digestion.

- 6.1 Remove Parafilm and **centrifuge** for  00:02:00 at **maximum speed** to pellet the bone powder.
- 6.2 Remove the supernatant and add  1 mL **fresh extraction buffer**.
- 6.3 Resuspend the bone powder by **vortexing**.
- 6.4 Incubate at  37 °C with **constant rotation** for ~  18:00:00 or overnight (use Parafilm to seal the lids).

2m

18h

Day 2 - DNA purification

41m

7 Concentration of DNA with Amicon vivaspin 30kD column



Note

- Amicon® vivaspin column are also called **centrifugal filter units**.
- When centrifuging the Amicon columns, do not close the tubes too tight.
- MinElute spin columns can be used instead of Monarch columns, though the latter are preferable for smaller volumes, like that used in the elution step.



7.1 Remove Parafilm from eppendorf and **centrifuge** for  00:02:00 min at **maximum speed** to pellet the bone powder.

2m

7.2 Transfer the **supernatant** (~1 mL) to an **Amicon® vivaspin 30kD column**.

7.3 **Centrifuge** at ~4,000g **until ~70 µL remains**. Start with a  00:10:00 min centrifugation, then adjust in  00:05:00 min intervals as needed.

15m

Note

- If the volume does not decrease to ~70 µL after several rounds of centrifuging, then repeat steps 7.4-7.5 until all the remaining sample is used.

7.4 **Mix** the **supernatant** with **10x** its volume of **modified binding buffer** in a **Monarch column**.

- Example: For 70 µL supernatant, add 700 µL of binding buffer.
- **Mix** by pipetting or inverting.

7.5 **Centrifuge** for  00:01:00 min at **8,000 RPM** and **discard the flowthrough**.

1m

- Tap the collection tube on a paper towel to get excess flowthrough out before reassembling.

7.6 Add  650 µL **PE buffer** to the spin column.

7.7 **Centrifuge** at **5,000 RPM** for  00:01:00 min and **discard the flowthrough**.

1m



7.8 Repeat steps 7.6-7.7.

7.9 Centrifuge at **maximum speed** for  00:01:00 min to dry the spin column.

1m



7.10 Place the column in a clean **1.5 mL LoBind** Eppendorf tube.



8 DNA Elution

10m 30s

8.1 Add  17.5 μ L **EB buffer** (or EBT buffer if not using LoBind tubes) to the silica membrane and **incubate** for  00:10:00 min at room temperature.

10m



8.2 **Centrifuge** at **maximum speed** for  00:00:30 seconds.

30s



8.3 **Repeat** steps 8.1-8.2 for a total of **35 μ L DNA extract**. Store DNA at -20°C.

9 Qubit Quantification

- If desired, measure the DNA concentration using a Qubit. See Qubit protocol for details.

Protocol references

Allentoft ME, Sikora M, Sjögren KG, et al. Population genomics of Bronze Age Eurasia. *Nature*. 2015;522(7555):167-172. doi:10.1038/nature14507

Dabney J, Knapp M, Glocke I, et al. Complete mitochondrial genome sequence of a Middle Pleistocene cave bear reconstructed from ultrashort DNA fragments. *Proc Natl Acad Sci U S A*. 2013;110(39):15758-15763. doi:10.1073/pnas.1314445110