

Jul 20, 2019

Anesthetizing Sea Cucumbers

DOI

dx.doi.org/10.17504/protocols.io.5nyg5fw



Jon Eilers¹

¹Walla Walla College



Jon Eilers

Walla Walla College

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN ACCESS



DOI: dx.doi.org/10.17504/protocols.io.5nyg5fw

Protocol Citation: Jon Eilers 2019. Anesthetizing Sea Cucumbers. **protocols.io** https://dx.doi.org/10.17504/protocols.io.5nyg5fw

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited



protocols.io Part of SPRINGER NATURE

Protocol status: Working

We use this protocol and it's working

Created: July 20, 2019

Last Modified: October 06, 2019

Protocol Integer ID: 26040

Keywords: anesthetizing sea cucumber, sea cucumbers this protocol, sea cucumber in the menthol, sea cucumber, minutes the sea cucumber, anesthetizing apostichopus japonicus, anesthestia, saltwaer solution, salt water

Abstract

This protocol is for anesthetizing sea cucumbers. After soaking the sea cucumber in the menthol/ethanol/saltwaer solution for 30-60 minutes the sea cucumber should be flacid and relatively unresponsive to poking. No sea cucumber has eviscerated during surgery after being anesthetized although they do contract and appear like they are trying. The anesthestia wears off quickly after being placed back in salt water.

This protocol is based on the method for anesthetizing Apostichopus japonicus as published in the citation below

Citation

YUSUKE YAMANA, TATSUO HAMANO AND KEN-ICHI YAMAMOTO (2005)

. Anesthetizer of the adult sea cucumber Apostichopus japonicus. NIPPON SUISAN GAKKAISHI 71(3):299-306.

10.2331/suisan.71.299

LINK

Materials

MATERIALS

Ethanol

Menthol

Salt Water



2

1 Add .5 grams of menthol crystals to 100 ml of ethanol and stir/shake until dissolved

Dilute to 10% menthol/ethanol solution by adding 900 ml of salt water

1m

In a separate container, add 400 ml diluted menthol/ethanol solution and 600 ml saltwater to create a 40% diluted solution

1m

Add sea cucumber to diluted solution and wait 30 minutes. After 30 minutes, gently poke and observe the reaction. If unresponsive, it is anesthetized. Otherwise, wait another 5-10 minutes and poke again.

30m

Expected result

Once anesthestized, the sea cucumber can be removed from the solution and will remain relatively unresponsive for 10-30 minutes.

Note

If feeding tentacles are relaxed and hanging out of the mouth, it is anesthetized.

Citations

YUSUKE YAMANA, TATSUO HAMANO AND KEN-ICHI YAMAMOTO. Anesthetizer of the adult sea cucumber Apostichopus japonicus

10.2331/suisan.71.299