

May 29, 2025

Animals: Nestin-Cre::DIO-SPOTlight transgenic mouse line generation

DOI

<https://dx.doi.org/10.17504/protocols.io.n2bvjb4pwgk5/v1>

Matt Oliver¹

¹Duke University



Matt Oliver

Duke University

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.n2bvjb4pwgk5/v1>

Protocol Citation: Matt Oliver 2025. Animals: Nestin-Cre::DIO-SPOTlight transgenic mouse line generation. **protocols.io**
<https://dx.doi.org/10.17504/protocols.io.n2bvjb4pwgk5/v1>

Manuscript citation:

Matthew L Oliver, Zachary F Caffall, Callie B Eatman, Timothy D Faw, Nicole Calakos
(2025) **DIO-SPOTlight Transgenic Mouse to Functionally Monitor Protein Synthesis Regulated by the Integrated Stress Response**
*eLife*14:RP104457<https://doi.org/10.7554/eLife.104457.1>

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited



Protocol status: Working

We use this protocol and it's working

Created: May 29, 2025

Last Modified: May 29, 2025

Protocol Integer ID: 219163

Keywords: ASAPCRN, DIO-SPOTlight , Transgenic mouse line, spotlight transgenic mouse line generation this protocol, spotlight transgenic mouse line generation, spotlight transgenic mouse line, nestin

Funders Acknowledgements:

ASAP

Grant ID: ASAP-020607

Abstract

This protocol describes the steps performed to generate Nestin-Cre::DIO-SPOTlight transgenic mouse line.

Nestin-Cre::DIO-SPOTlight transgenic mouse line generation

- 1 To produce the Rosa26-DIO-SPOTlight mouse clone the DIO-SPOTlight construct into the Ai9 vector backbone in place of the Lox-STOP-Lox-tdTomato cassette.
- 2 Transfer the linearized Ai9-DIO-SPOTlight plasmid to an available transgenic mouse core (rodent genetic engineering services).
- 3 Integrate the the linearized Ai9-DIO-SPOTlight plasmid into the Rosa26 locus using homologous recombination and electroporate into hybrid 129S6/C57BL/6N (G4) embryonic stem cells.
- 4 To achieve pan-neuronal expression cross Rosa26-DIO-SPOTlight hemizygous transgenic founders with Nestin-Cre mice.

Acknowledgements

- Ai9 vector backbone was a gift from Hongkui Zeng; Addgene plasmid # 22799; <http://n2t.net/addgene:22799>; RRID:Addgene_22799.
- Duke Transgenic Mouse Core performed generation of Rosa26-DIO-SPOTlight hemizygous transgenic founders
- Nestin-Cre mice were obtained from Jackson Laboratories; JAX: strain #003771, RRID: IMSR_JAX:003771 (Tronche, et al 1999)