

Dec 10, 2019

Anti-iNKT MicroBeads Isolation protocol

DOI

dx.doi.org/10.17504/protocols.io.bacuiaww

Ajit N Shah, Rachel Hatano¹

¹:



Ajit N Shah

OPEN  ACCESS



DOI: dx.doi.org/10.17504/protocols.io.bacuiaww

Protocol Citation: Ajit N Shah, Rachel Hatano 2019. Anti-iNKT MicroBeads Isolation protocol. **protocols.io**
<https://dx.doi.org/10.17504/protocols.io.bacuiaww>

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: December 10, 2019

Last Modified: December 10, 2019

Protocol Integer ID: 30836

Abstract

isolation of 6b11 positive cells

Materials

MATERIALS

 Anti-iNKT MicroBeads human **Miltenyi Biotec Catalog #130-094-842**

- 1 Determine cell number.
- 2 Centrifuge cell suspension at 300×g for 10 minutes. Aspirate supernatant completely.
- 3 Resuspend cell pellet in 400 µL of buffer per 10⁸ total cells. Add 100 µL of Anti-iNKT MicroBeads per 10⁸ total cells.
- 4 Mix well and incubate for 15 minutes in the refrigerator (2–8 °C).
- 5 Wash cells by adding 1–2 mL of buffer per 10⁸ cells and centrifuge at 300×g for 10 minutes. Aspirate supernatant completely.
- 6 Resuspend up to 10⁸ cells in 500 µL of buffer.
- 7 Place column in the magnetic field of a suitable MACS Separator.
- 8 Prepare column by rinsing with the appropriate amount of buffer: LS: 3 mL
- 9 Apply cell suspension onto the column. Collect flow-through containing unlabeled cells.
- 10 Wash column with the appropriate amount of buffer. Collect unlabeled cells that pass through and combine with the effluent. LS: 3×3 mL.
- 11 Remove column from the separator and place it on a suitable collection tube.
- 12 Pipette the appropriate amount of buffer onto the column.
- 13 Immediately flush out the magnetically labeled cells by firmly pushing the plunger into the column. LS: 5 mL



- 14 . To increase the purity of iNKT cells, the eluted fraction is enriched over a second MS or LS Column.