



Nov 21, 2025

Version 1

APOE Genotyping Using TaqMan + Data Processing V.1

DOI

<https://dx.doi.org/10.17504/protocols.io.n92ldzo68v5b/v1>

Mario Mosqueda¹

¹TGen, Precision Aging Network

Mario Mosqueda: On behalf of the Precision Aging Network



Mario A Mosqueda

TGen

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.n92ldzo68v5b/v1>

Protocol Citation: Mario Mosqueda 2025. APOE Genotyping Using TaqMan + Data Processing . **protocols.io**
<https://dx.doi.org/10.17504/protocols.io.n92ldzo68v5b/v1>

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working



Created: August 20, 2022

Last Modified: November 21, 2025

Protocol Integer ID: 68935

Keywords: protocol for apoe4, snps rs7412, apoe4, apo, using taqman, rs429358, protocol

Funders Acknowledgements:

Precision Aging Network

Grant ID: U19AG065169

Abstract

Protocol for APOE4 genotyping for SNPs rs7412 (C/T) and rs429358 (C/T)



Thaw DNA Samples

- 1 Thaw DNA samples and APOE Controls at room temp, gently vortex and briefly centrifuge to collect the DNA solution into the bottom of the tube or well-plate.. 20m
- 1.1 While the samples thaw, label a 96 well PCR plate "DNA Dilution"
- 1.2 Create a plate map, taking into account your APOE controls and Negative Control.

Controls:
 - 2/2 (DNA)
 - 2/3 (DNA)
 - 2/4 (DNA)
 - 3/3 (DNA)
 - 3/4 (DNA)
 - 4/4 (DNA)
 - Negative Control (NC) (UP H2O)
- 2 Dilute DNA in UP H2O (Sample and Control) to make 20 μ L at 0.2 - 2.0 ng/ μ L, seal and briefly centrifuge to collect the DNA solution into the bottom of each well. 1h

Thaw Reagents

- 3 Thaw TaqMan assay components at room temperature, flick and swirl Master Mix, and gently vortex SNPs and centrifuge as needed to collect assay components into the bottom of their respective tubes. 20m
- 3.1 While the assay components thaw, label two MicroAmp well plates with each SNP ID.

rs7412
rs429358

Make Reaction Mix

- 4 Get two clean microcentrifuge tubes and add the appropriate volume of Master Mix to both:

You can't use 'macro parameter character #' in math mode
- 4.1 Label each tube with one of the SNPs IDs.



- 5 Add the appropriate volume of each SNP assay to its corresponding tube:

15m

You can't use 'macro parameter character #' in math mode

- 6 Add 11 μ L of the Reaction Mix to each well that will be used in its corresponding MicroAmp plate. Seal and briefly centrifuge to collect solution into the bottom of each well.

20m

Continuation Check

- 7 Proceed To Step 8 If you've confirmed that the QuantStudio Is available and ready.

Add DNA to plate

- 8 Transfer 9 μ L of each sample from the DNA dilution plate to each well of the MicroAmp plate to have a total volume of 20 μ L per well (9 μ L of sample + 11 μ L of Reaction Mix added in the previous step).

30m

Seal with MicroAmp Adhesion Film. vortex gently to mix.

Centrifuge (300 g) briefly to collect the solution into the bottom of the well plate. Make sure no bubbles are present in the wells.

QuantStudio instrument run to perform qPCR SNP genotyping

- 9 Place MicroAmp plate in the QuantStudio instrument and set up experiment as directed by the [TaqMan SNP Genotyping Assays User Guide](#) and the [QuantStudio User Guide](#).
- 10 Once the Experiment has finished check to see that your Allelic Discrimination plot looks resembles:




Figure 1 Typical allelic discrimination plot

This figure is found in the TaqMan® SNP Genotyping Assays USER GUIDE linked in the step above. The guide also has observed issues and its possible cause if your allelic distribution graph doesn't resemble Figure 1.

- 11 Export Excel file of "Results" for genotype determinations. Make sure to Check "Well position" , "Call"

Download the Allelic Discrimination Plot

Data Processing

- 12 Using the calls made for rs7412 and rs429358 we can determine the APOE Genotype.

Download the Excel File below and fill the "Calls" Sheet.

Fill out your "qPCR Sample Name" and transfer the "Call" made from both the rs7412 and rs429358 "Results" excel files exported from the step above.

It should determine the APOE Genotype based on the calls provided along with the Allelic Frequencies in your sample population.

Add the Allelic Discrimination Plots



APOE Calculator template 8_4_202...

- 13 The calls from two SNPs, rs7412 (C/T) and rs429358 (C/T), determine each allele's APOE genotype. The combination of these SNPs, listed in the table below, defines the APOE status.

	A	B	C
	APOE Status	rs429358	rs7412
	e2e4	CT	CT
	e2e2	TT	TT
	e2e3	TT	CT
	e3e3	TT	CC
	e3e4	CT	CC
	e4e4	CC	CC