


Oct 20, 2019

## Aptamer conjugated beads affinity assay

DOI

[dx.doi.org/10.17504/protocols.io.8hxht7n](https://dx.doi.org/10.17504/protocols.io.8hxht7n)



Jorge Fernandez Mendez<sup>1</sup>

<sup>1</sup>Universidad Complutense de Madrid

AEGIS - Madrid iGEM 2019



Jorge Fernandez Mendez

Universidad Complutense de Madrid

### Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.8hxht7n>

**Protocol Citation:** Jorge Fernandez Mendez 2019. Aptamer conjugated beads affinity assay. **protocols.io**  
<https://dx.doi.org/10.17504/protocols.io.8hxht7n>



**License:** This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

**Protocol status:** Working

**We use this protocol and it's working**

**Created:** October 20, 2019

**Last Modified:** October 20, 2019

**Protocol Integer ID:** 28951

**Keywords:** affinity by an aptamer, aptamer, beads affinity, assay, treated well, target immobilization, bead, affinity, conjugated particle

## Abstract

The following protocol details how to test affinity by an aptamer-conjugated particles with their targets, via target immobilization in 96-treated wells.

## Guidelines

The steps for coating the wells with the target can differ depending on the molecule to immobilize in the plate.

## Materials

### MATERIALS



Nunc®; FluoroNunc®;LumiNunc®; 96-Well Plates, clear, C-shaped Bottom, MaxiSorp, Certified, Clear, For Time-Resolved Fluorescence **Thermo Fisher Catalog #437958**

- Streptavidin - biotylated DNA x18 thymine sequence at 8  $\mu$ M concentration.
- E. Coli suspension in Carbonate-Bicarbonate buffer pH = 9.6 (OD600 = 0.3 )
- Conjugated Latex beads Stock at 0.1 % wt in PBS-T 0.05% 1.4 mM  $MgCl_2$
- 5% BSA solution in PBS buffer.
- PBS-T (Tween-20 0.05 %) buffer for washings.

## Troubleshooting



## Well Coatings

- 1 Add to each well of the 96 maxisorp plate 100  $\mu$ L of E. Coli or Streptavidin-DNA suspensions. And incubate overnight at 4 °C.
- 2 After incubation, remove the liquid in the plate flipping the 96-plate.  
  
! Don't aspire the liquid in the wells with a pipette, to avoid damaging the coating.
- 3 Add 150  $\mu$ L of 3% BSA solution to each well. Make two additional BSA additions for negative controls. Incubate the 96-wells for 4h with mild agitation.

## Affinity Assay

- 4 Add 80  $\mu$ L of each conjugated latex beads stock to the coated wells. Remember to add at least in two wells, coated with the different targets respectively; E. Coli and Streptavidin-DNA conjugate. Incubate under mild agitation at room temperature for 1 hour.
- 5 Remove the liquid by plate flipping, as mentioned in step 2.
- 6 Wash the wells adding 200 $\mu$ L PBS-T (Tween-20 0.05 %), incubating 2 minutes under mild agitation, and flipping the 96-plate. Repeat this step twice.
- 7 Cover the wells with 50  $\mu$ L of distilled water and measure absorbance at the absorbtion peak for latex beads.

For the the fluorescent orange latex beads we have employed, the absorbsnce measurement have been performed at 570nm.