


Jun 07, 2023

Version 1

# Aqueous (SBIp) Delipidation of a Whole Mouse Brain V.1

DOI

[dx.doi.org/10.17504/protocols.io.n2bvj81mwgk5/v1](https://dx.doi.org/10.17504/protocols.io.n2bvj81mwgk5/v1)



Andrew Recknagel<sup>1</sup>, Kevin Cao<sup>2</sup>, Naveen Ouellette<sup>2</sup>, Molly Logsdon<sup>2</sup>, Judith Baka<sup>2</sup>, Jayaram Chandrashekar<sup>2</sup>

<sup>1</sup>Janelia Farm; <sup>2</sup>Allen Institute for Neural Dynamics

Allen Institute for Neural...



Naveen Ouellette

Allen Institute

## Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.n2bvj81mwgk5/v1>

**Protocol Citation:** Andrew Recknagel, Kevin Cao, Naveen Ouellette, Molly Logsdon, Judith Baka, Jayaram Chandrashekar 2023.

Aqueous (SBIp) Delipidation of a Whole Mouse Brain. **protocols.io**

<https://dx.doi.org/10.17504/protocols.io.n2bvj81mwgk5/v1>



**License:** This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

**Protocol status:** Working

**We use this protocol and it's working**

**Created:** March 24, 2023

**Last Modified:** June 19, 2024

**Protocol Integer ID:** 79363

**Keywords:** delipidation, aqueous delipidation, clearing, brain delipidation, lipid removal via phase separation, whole mouse brain aqueous strategy, organic delipidation, delipidation, lipid removal, aqueous biphasic buffer, sbip, phosphate buffer, solvent, lipid, aqueous strategy, forming micelle, phase separation, antibody, whole mouse brain, residual membrane, b1n buffer, antibody labeling, brain, detergent, buffer

**Funders Acknowledgements:**

Allen Institute

## Abstract

Aqueous strategies for whole-brain delipidation involve lipid removal via phase separation and detergent washes. SBiP is an aqueous biphasic buffer that extracts lipids at polar-nonpolar solvent interfaces. The brain is then washed with the detergent-based B1n buffer, which further disrupts residual membranes, forming micelles that can be washed out. When paired with organic delipidation, these steps can return a solvent-shrunken brain to normal size in phosphate buffer, suitable for post-delipidation antibody labeling.

## Guidelines

It is recommended to etch identification information on the glass vial. Ink or label adhesive may dissolve if exposed to organic solvents.



## Materials

⊗ Sodium Dodecyl Sulfate **Merck MilliporeSigma (Sigma-Aldrich) Catalog #74225**

⊗ Sodium phosphate dibasic **Merck MilliporeSigma (Sigma-Aldrich) Catalog #7558-79-4**

⊗ Sodium phosphate monobasic monohydrate **Merck MilliporeSigma (Sigma-Aldrich) Catalog #S9638**

⊗ 2-methyl-2-butanol **Merck MilliporeSigma (Sigma-Aldrich) Catalog #152463**

⊗ 2-propanol **Merck MilliporeSigma (Sigma-Aldrich) Catalog #278475**

⊗ Glycine **Fisher Scientific Catalog #BP381-500**

⊗ PBS - Phosphate-Buffered Saline (10X) pH 7.4 **Invitrogen - Thermo Fisher Catalog #AM9625**

⊗ Triton X-100 **Merck MilliporeSigma (Sigma-Aldrich) Catalog #T8787-50ML**

⊗ 5% Sodium Azide **Fisher Scientific Catalog #71448-16**

⊗ 10N NaOH **Merck MilliporeSigma (Sigma-Aldrich) Catalog #SX0607N-6**

## Equipment

**WHEATON® Liquid Scintillation Vials, Caps Attached to Vials, Glass, Polyethylene Cone, 22-400, 20 mL** NAME

Vial TYPE

Wheaton BRAND

DWK986546 SKU

<https://www.dwk.com/na/wheaton-liquid-scintillation-vials-caps-attached-to-vials-glass-polyethylene-cone-22-400-20-ml-986546> LINK

20 mL Glass Vial with Polyethylene cone Caps SPECIFICATIONS



## Equipment

**Fisherbrand™ Multi-Purpose Tube Rotator** NAME

Carousel TYPE

FisherBrand BRAND

88-861-049 SKU

<https://www.fishersci.com/shop/products/multi-purpose-tube-rotators/88861049> LINK





## RECIPES

**10 mM Phosphate Buffer pH 8.3**

Combine the following reagents, adjust pH to 8.3.



	Reagent	Amount	Final Concentration
	1M Phosphate Buffer	5 mL	10 mM
	Milli-Q water	495 mL	

### **SBiP Solution: 0.08% SDS, 16% 2-Methyl-2-Butanol, 8% 2-Propanol, in H<sub>2</sub>O**

Combine the following reagents on ice. Use a fume hood when adding 2-methyl-2-butanol and 2-propanol. Mix  On ice until solution is uniform and clear. Store immediately at  4 °C until ready for use. Use each batch within a month for best effect.

	Reagent	Volume
	Milli-Q water (ice cold)	350 mL
	50mM Na <sub>2</sub> HPO <sub>4</sub>	2 mL
	4% SDS (in H <sub>2</sub> O, pH 7.4)	10 mL
	2-Methyl-2-butanol	80 mL
	2-Propanol	40 mL
	<b>Total</b>	<b>482 mL</b>

### **B1n Buffer: 0.1% Triton X-100, 2% Glycine, 0.02% NaN<sub>3</sub> in H<sub>2</sub>O**

Combine the following reagents and stir at  Room temperature until fully dissolved. Store for a few months at  Room temperature .


	A	B
	Milli-Q Water	up to 500 mL
	Triton X-100	500 uL



	A	B
	Glycine	10 g
	10N NaOH	50 $\mu$ L
	5% Sodium azide	2 mL
	<b>Total</b>	<b>500 mL</b>

## Troubleshooting

## Safety warnings

 2-Methyl-2-butanol and 2-propanol are corrosive and flammable. Perform the steps that involve these reagents under the fume hood. Wear a lab coat, safety goggles or glasses, and gloves.

## Before start

If performing this delipidation as part of the **Whole Mouse Brain Delipidation, Immunolabeling, and Expansion Microscopy** protocol, first delipidate the whole mouse brain with organic solvents using the **THF and DCM Delipidation of a Whole Mouse Brain** protocol.



## SBiP Delipidation


5d 12h








### 1 Delipidate with aqueous SBiP solution

Start with a whole mouse brain perfused with 4% PFA, post-fixed, stored in PBS in a 20 mL vial.


### 2 Use a 20 mL vial for processing an adult mouse brain. All steps in this section are carried out on a carousel rotator 10 rpm, Room temperature .

4d 9h

Replace solution in vial with  20 mL SBiP for each of the following steps:

- SBiP for  03:00:00
- SBiP for  06:00:00
- SBiP  Overnight
- SBiP for  24:00:00
- SBiP for  24:00:00
- SBiP for  24:00:00
- SBiP for  24:00:00








#### Note

It is important that the SBiP solution is mixed thoroughly during the aqueous steps in order to effectively delipidate the brain. When it is thoroughly mixed, the solution should look uniformly turbid at  Room temperature .

#### Safety information

2-Methyl-2-butanol and 2-propanol are corrosive and flammable. Perform the steps that involve these reagents under the fume hood. Wear a lab coat, safety goggles or glasses, and gloves.



- 3 Wash the brain with  20 mL B1n buffer at  Room temperature  Overnight .
- 4 Replace B1n buffer to complete ~  24:00:00 wash. 1d
- 5 Wash the brain with 1X PBS, rotating at  Room temperature for the following steps: 3h
  - PBS for  01:00:00
  - PBS for  02:00:00 +
- 6 Aqueous (SBIp) delipidation complete. Store in 1X PBS 0.05% Azide for up to 6 months.