


Jul 23, 2019

# 🌍 Arabidopsis growth medium with varying sulfur and nitrogen content

DOI

[dx.doi.org/10.17504/protocols.io.5q6g5ze](https://dx.doi.org/10.17504/protocols.io.5q6g5ze)



Meike Burow<sup>1</sup>

<sup>1</sup>University of Copenhagen



Verena Jeschke

OPEN  ACCESS



DOI: [dx.doi.org/10.17504/protocols.io.5q6g5ze](https://dx.doi.org/10.17504/protocols.io.5q6g5ze)

**Protocol Citation:** Meike Burow 2019. Arabidopsis growth medium with varying sulfur and nitrogen content. **protocols.io** <https://dx.doi.org/10.17504/protocols.io.5q6g5ze>

**License:** This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

**Protocol status:** Working

**We use this protocol and it's working**

**Created:** July 23, 2019

**Last Modified:** July 23, 2019

**Protocol Integer ID:** 26110

**Keywords:** medium, low sulfur, low nitrogen, seedling

## Abstract

Media recipe to grow *A. thaliana* seedlings on plates.  
Variation in S and N concentrations.





## stock solutions

1

Media solution	stock solutions
<b>+S/+N</b>	A+C+G
<b>+S/-N</b>	A+D+G
<b>-S/+N</b>	B+E+H
<b>-S/-N</b>	B+F+H

 media recipes.xlsx

## for 1 L medium

2

General media recipe	Volume
Minor salts (100x)	10 mL
Major salts (10x)	100 mL
Ferrous solution (100x)	10 mL
Vitamin mix (1000x, Murashige and Skoog, Prod. No. M0409.0100, Duchefa Biochemie)	1 mL
water to 1L	
adjust pH to 5.8 with 1M KOH	
divide 1L to 2× 500 mL (into 1L bottles)	
add 4g agar per 500mL	
autoclave	

## Concentrations of S and N in final medium

3

Concentration of Sulphate
in +S: 1.68 mM
in -S: 0.015 mM
Concentration of Nitrate

in +N: 3 mM
in -N: 0.3 mM