

Mar 14, 2020

Version 2

Artificial Cerebrospinal Fluid IV (ACSF.IV) V.2

DOI

dx.doi.org/10.17504/protocols.io.bdpbi5in

Allen Institute for Brain Science¹

¹Allen Institute

BICCN / BICAN

Allen Institute for Brain S...



Dillan Brown

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.bdpbi5in>

Protocol Citation: Allen Institute for Brain Science 2020. Artificial Cerebrospinal Fluid IV (ACSF.IV). **protocols.io**
<https://dx.doi.org/10.17504/protocols.io.bdpbi5in>

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: March 14, 2020

Last Modified: March 14, 2020

Protocol Integer ID: 34243

Keywords: Artificial Cerebrospinal Fluid, ACSF, RP0222, incubation, electrophysiological, recording, artificial cerebrospinal fluid iv, making artificial cerebrospinal fluid iv, incubation of fresh mouse brain slice, fresh mouse brain slice, mouse brain, electrophysiological recording, mouse, national institutes of health, brain slice, including incubation, national institute of mental health,

Abstract

This protocol describes making Artificial Cerebrospinal Fluid IV (ACSF.IV). ACSF.IV is used for multiple applications including incubation of fresh mouse brain slices prior to electrophysiological recording.

Note: Research reported in this publication was supported by the National Institute Of Mental Health of the National Institutes of Health under Award Number U19MH114830. The content is solely the responsibility of the authors and does not necessarily represent the official views of the National Institutes of Health.

Attachments



RP0222_Artificial_Ce...

54KB

Troubleshooting

