

Mar 02, 2020

Version 1

Artificial Cerebrospinal Fluid VIII (ACSF.VIII) V.1

DOI

<https://dx.doi.org/10.17504/protocols.io.bctqiwmw>

Allen Institute for Brain Science¹

¹Allen Institute

BICCN / BICAN

Allen Institute for Brain S...



Dillan Brown

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.bctqiwmw>

Protocol Citation: Allen Institute for Brain Science 2020. Artificial Cerebrospinal Fluid VIII (ACSF.VIII). **protocols.io**
<https://dx.doi.org/10.17504/protocols.io.bctqiwmw>

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited



Protocol status: Working

We use this protocol and it's working

Created: February 20, 2020

Last Modified: March 02, 2020

Protocol Integer ID: 33360

Keywords: Artificial Cerebrospinal Fluid, ACSF, incubation, recording, electrophysiological, RP0228, artificial cerebrospinal fluid viii,

Abstract

This protocol describes how to make Artificial Cerebrospinal Fluid VIII (ACSF.VIII). ACSF.VIII is used for multiple applications including incubation of fresh human brain slices prior to electrophysiological recording.

Attachments



RP0228_Artificial_Ce...

60KB

