

Oct 12, 2016

Version 2

ASP12A Recipe for culturing *Aureococcus anophagefferens* V.2

DOI

<https://dx.doi.org/10.17504/protocols.io.f3ybgpw>

Eric Gann

Protist Research to Opti...

The Aquatic Microbial E...



Ashley A Humphrey

University of Tennessee, Knoxville

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.f3ybgpw>

External link: <http://www.nefsc.noaa.gov/nefsc/publications/crd/crd0517/crd0517.pdf>

Protocol Citation: Eric Gann 2016. ASP12A Recipe for culturing *Aureococcus anophagefferens*. protocols.io <https://dx.doi.org/10.17504/protocols.io.f3ybgpw>

License: This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working



Created: October 12, 2016

Last Modified: June 11, 2019

Protocol Integer ID: 3928

Keywords: culturing aureococcus anophagefferen, aureococcus anophagefferen, asp12a recipe, recipe

Abstract

Defined Medium for Aureococcus anophagefferens

Guidelines

	Anhydrous Salts						
	Constituent	Final Conc.	Molar Mass	g/1.0L	g/0.75L	g/0.5L	g/0.25L
	NaCl	4.28E-01	58.44	25.01232	18.75924	12.50616	6.25308
	KCl	1.10E-02	74.5513	0.8200643	0.615048225	0.410032	0.205016
	Na2SO4	1.00E-02	142.04	1.4204	1.0653	0.7102	0.3551
	NaHCO3	2.00E-03	84.007	0.168014	0.1260105	0.084007	0.042004
	Hydrous Salts						
	Constituent	Final Conc.	Molar Mass	g/1.0L	g/0.75L	g/0.5L	g/0.25L
	CaCl2·2H2O	1.25E-02	147.02	1.83775	1.3783125	0.918875	0.459438
	MgCl2·6H2O	4.00E-02	203.3027	8.132108	6.099081	4.066054	2.033027
	Macronutrients						
	Constituent	Final Conc.	Molar Mass	g/1.0L	g/0.75L	g/0.5L	g/0.25L
	NaNO3	1.47E-03	84.9947	0.124942209	0.093706657	0.062471	0.031236
	Ca2·glycerophosphate	2.00E-05	210.14	0.0042028	0.0031521	0.002101	0.001051
	K2HPO4	4.00E-05	174.2	0.006968	0.005226	0.003484	0.001742
	Other						
	Constituent	Final Conc.	Molar Mass	g/1.0L	g/0.75L	g/0.5L	g/0.25L
	Tris Base	4.13E-03	121.14	0.5003082	0.37523115	0.250154	0.125077
	Stock Solutions						
	Solution	Final Conc.	g/L	mL/1.0L	mL/0.75L	mL/0.5L	mL/0.25L
	Fe-EDTA sodium Salt	1.79E-06	0.65702	1	0.75	0.5	0.25

	TSM - I			1	0.75	0.5	0.25
	TSM - II			1	0.75	0.5	0.25
	TSM - III			1	0.75	0.5	0.25
	S3 Vitamin Solution			1	0.75	0.5	0.25
	Trace Metal Solution - TMS - I						
	Constituent	Final Conc.	Molar Mass	g/1.0L 1000X Stock			
	Na2EDTA·2H2O	8.53E-05	372.24	31.752072			
	FeCl3·6H2O	1.79E-06	270.3	0.483837			
	H3BO3	2.77E-04	61.83	17.12691			
	MnCl2·4H2O	9.10E-06	197.91	1.800981			
	ZnSO4·7H2O	1.07E-06	287.56	0.3076892			
	NaMoO4·2H2O	5.21E-06	218.9779	1.140874859			
	CoCl2·6H2O	2.21E-07	237.93	0.05258253			
	Trace Metal Solution - TMS - II						
	Constituent	Final Conc.	Molar Mass	g/1.0L 1000X Stock			
	KBr	1.25E-04	119.002	14.87525			
	SrCl2 · 6H2O	3.42E-05	266.62	9.118404			
	RbCl	2.34E-06	120.92	0.2829528			
	Trace Metal Solution - TMS - III						
	Constituent	Final Conc.	Molar Mass	g/1.0L 1000X Stock			
	Na2SeO3	2.00E-08	172.94	0.0034588			
	NaF	5.26E-05	41.98871	2.208606146			
	NaAsO2	4.00E-08	129.91	0.0051964			
	Al2(SO4)3·18H2O	1.11E-06	666.43	0.7397373			
	NH4VO3	2.00E-08	116.98	0.0023396			
	S3 Vitamin Solution						



	Constituent	Final Conc.	Molar Mass	g/1.0L 1000X Stock			
	Thiamine HCL	2.96E-07	337.27	0.09983192			
	biotin	4.10E-09	244.31	0.001001671			
	vitamin B12	7.40E-10	1.36E+03	0.0010064			



- 1 Prepare Trace Metal Solutions I - III, Fe-EDTA Sodium Salt Solution, and Vitamin Solution. Add constituents to 750 mL MilliQ water, shake to dissolve. Add MilliQ water to 1.0 L. Filter sterilize. Trace Metal Solutions and Fe-EDTA Sodium Salt Solution are stored at 4°C in 15 - 50 mL aliquots, while the Vitamin solution is stored at -20°C in 5 - 10 mL aliquots.
- 2 Add anhydrous salts, hydrous salts, macronutrients and Tris Base to 75% desired MilliQ water while stirring to dissolve. Add MilliQ water to final volume, autoclave.
- 3 Once cooled, add in 1mL of the Trace Metal Solutions, Fe-EDTA Sodium Salt Solution, and Vitamin Solution per 1L of media made.