

Jun 19, 2025

Assay for lysosomal activity using Magic Red

 [Nature Communications](#)

DOI

<https://dx.doi.org/10.17504/protocols.io.5jyl8q437l2w/v1>

Bishal Basak^{1,2}, Erika Holzbaur^{1,2}

¹Department of Physiology, Perelman School of Medicine, University of Pennsylvania, Philadelphia, PA, 19104;

²Aligning Science Across Parkinson's (ASAP) Collaborative Research Network, Chevy Chase, MD, 20815

University of Pennsylv...



Bishal Basak

University of Pennsylvania

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.5jyl8q437l2w/v1>

Protocol Citation: Bishal Basak, Erika Holzbaur 2025. Assay for lysosomal activity using Magic Red . **protocols.io**
<https://dx.doi.org/10.17504/protocols.io.5jyl8q437l2w/v1>

Manuscript citation:

<https://www.biorxiv.org/content/10.1101/2024.06.12.598718v1.full>

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: June 19, 2025

Last Modified: June 19, 2025

Protocol Integer ID: 220604

Keywords: lysosomal activity, using magic red, assay, magic red

Funders Acknowledgements:

ASAP

Grant ID: ASAP-000350

Abstract

Assay for lysosomal Cathepsin B activity in neurons using Magic Red

Protocol materials

 Magic Red substrate **Bio-Rad Laboratories Catalog #ICT937**



Plate neurons

- 1 Plate 120,000–150,000 neurons in an imaging dish at days in vitro (DIV) 0

Note

Protocol for primary neuron isolation and in vitro culture can be obtained from:
dx.doi.org/10.17504/protocols.io.81wgby723vpk/v1

(DIV) 0

Preparing reagents

- 2 Reconstitute Magic Red substrate upon arrival in  50 μL of DMSO

 Magic Red substrate **Bio-Rad Laboratories Catalog #ICT937**

- 3 On the day of imaging prepare 1:10 fresh dilution of reconstituted Magic Red solution in water
- 4 On the day of imaging equilibrate 4 ml of fresh maintenance media (MM) and 4 ml of imaging media (IM) per imaging dish in an incubator

Note

Maintenance media composition can be obtained from:
dx.doi.org/10.17504/protocols.io.81wgby723vpk/v1

Imaging media composition can be obtained from:
dx.doi.org/10.17504/protocols.io.36wgqgkz3gk5/v1

Imaging steps

15m

- 5 At DIV 6 or 7, aspirate out conditioned media from the imaging dish, wash the neurons once with fresh MM
- 6 Add  2 mL of fresh MM and then add  5 μL of the diluted Magic Red solution (prepared in step 3)
- 7 Incubate the neurons in an incubator at  37 °C for  00:15:00

15m



- 8 Post incubation, aspirate out the media, and wash once with fresh IM
- 9 Add  2 mL of fresh IM and then add  5 μ L of the diluted Magic Red solution (prepared in step 3)
- 10 Image neurons live under a confocal microscope with a heating chamber kept at  37 °C