

Dec 11, 2025

## Assessing Growth-Limiting Factors of *C. elegans* in 96-Well Liquid Culture

DOI

<https://dx.doi.org/10.17504/protocols.io.36wgqped5vk5/v1>

John Bergqvist<sup>1</sup>

<sup>1</sup>MRC Laboratory of Medical Sciences, Imperial College London

Behavioural Genomics



John Bergqvist

MRC Laboratory of Medical Sciences, Imperial College London

### Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



**DOI:** <https://dx.doi.org/10.17504/protocols.io.36wgqped5vk5/v1>

**Protocol Citation:** John Bergqvist 2025. Assessing Growth-Limiting Factors of *C. elegans* in 96-Well Liquid Culture. protocols.io <https://dx.doi.org/10.17504/protocols.io.36wgqped5vk5/v1>

**License:** This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited



**Protocol status:** Working

**We use this protocol and it's working**

**Created:** October 02, 2025

**Last Modified:** December 11, 2025

**Protocol Integer ID:** 228866

**Keywords:** C. elegans, liquid assay, paraformaldehyde, growth conditions, full growth of the strain, worm length measurement, growth of the strain, shaking incubation, well liquid culture the purpose, static incubation, incubating round, well liquid culture, worm, liquid medium, strain

## Abstract

The purpose of this protocol is to determine what factor or factors limit the growth of a C. elegans strain from L1 to adulthood in liquid medium in 96-well plates. The hypotheses to be tested are 1) PFA-killed OP50 limits the growth of the strain compared to alive OP50, 2) static incubation limits the growth of the strain compared to shaking incubation, and 3) wells contain too little S-medium to support full growth of the strain. These hypotheses will be tested by incubating round, F-bottom 96-well plates that are either shaking or static with PFA-killed OP50 or alive OP50, and with worms being dispensed in M9-medium or S-medium. The plates are recorded at 25 fps for 30s, where the first 10s are without bluelight stimulation and the final 10s are with bluelight stimulation. The readout of this protocol is a pixel variation analysis and worm length measurement.

## Materials

- 6 mL PFA-killed OP50 (OD = 2.0)
- 6 mL alive OP50 (OD = 2.0)
- 10,000x L1 N2
- 10,000x L1  $\Delta$ nhf-8
- 8x 96-WPs
- 6x seeded 100 mm plates
- 15 mL falcon tubes
- Pasteur pipettes
- M9 medium
- S.H2O
- 20°C shaking incubator over 76h



## Make sure the following are prepared

### 1 0.5% PFA-treated OP50

#### Protocol

NAME

**PFA treatment of OP50**

CREATED BY

Bonnie Evans

[Preview](#)

### 2 Make S-medium

#### Note

#### 1 L Pre-autoclave S-medium (S-basal)

|                                 |        |           |
|---------------------------------|--------|-----------|
| NaCl                            | 5.85 g | 4th floor |
| K <sub>2</sub> HPO <sub>4</sub> | 1 g    | 4th floor |
| KH <sub>2</sub> PO <sub>4</sub> | 6 g    | 4th floor |
| MilliQ water                    | 1 L    | 4th floor |

#### Note

#### 1 L Post-autoclave S-medium

|                           |             |            |      |           |
|---------------------------|-------------|------------|------|-----------|
| Cholesterol               | [M] 5 mg/mL | in ethanol | 1 mL | 4th floor |
| 1 M potassium citrate Ph6 | 10 mL       |            |      | 4th floor |
| trace metals solution     | 10 mL       |            |      | 4th floor |
| 1 M CaCl <sub>2</sub>     | 3 mL        |            |      | 4th floor |
| 1 M MgSO <sub>4</sub>     | 3 mL        |            |      | 4th floor |



## Day -8

- 3 Pick 10-20 adult worms onto 3 seeded 100 mm plates per strain in afternoon.

## Day -4

- 4 Wash off the adult worms with M9.
- 5 Add a couple of mL of M9 to the plates and streak a plastic sterile loop over the agar to loosen eggs.
- 6 Pour the eggs into a 15 mL falcon tube.
- 7 Repeat step 3-4 until few eggs are left on the plate.
- 8 Wash the eggs in M9 three times.
- 9 Leave eggs in ~9 mL M9 in empty 60 mm plates over-night at 20°C.

## Day -3

- 10 Wash L1s once in M9 in a 15 mL falcon tube.
- 11 Half the population into another 15 ml falcon tube, spin the worms down and replace the M9 medium with S-medium. You should have one tube with L1 worms in M9 medium and one tube with worms in S-medium now.
- 12 Adjust the worm number to obtain 20 worms per 20  $\mu$ L.
- 13 Add 10  $\mu$ L of Sterile H<sub>2</sub>O to all eight 96-well plates with the ViaFill. After dispensing S.H<sub>2</sub>O, prime once with MilliQ water and then prime the ViaFill without MilliQ water so



that no more liquid is left in the tubes.

- 14 Divide up the tubes so that the tubes dispensing in rows A to D are for dispensing PFA-killed OP50 and tubes dispensing in rows E to H are for dispensing live OP50.
- 15 Dispense 30  $\mu\text{L}$  of PFA-killed OP50 (pellet resuspended in S-medium; OD=2.0; rows A-D) and alive OP50 (pellet resuspended in S-medium; OD=2.0, rows E-H). After dispensing bacteria, prime 2x with MilliQ water and then prime the ViaFill without MilliQ water so that no more liquid is left in the tubes.
- 16 Divide the ViaFill tubes into two groups so that the tubes covering rows A, B, E and F are sourcing worms from Strain 1, and tubes covering rows C, D, G, H are sourcing worms from Strain 2.
- 17 Shake the two M9-containing falcon tubes for Strain 1 and Strain 2, and stick the dedicated groups of tubes into the respective falcon tube. Dispense 20  $\mu\text{L}$  of worms into Plate 1 and Plate 2. After dispensing the worms, prime 2x with MilliQ water and then prime the ViaFill without MilliQ water so that no more liquid is left in the tubes.
- 18 Shake the two S medium-containing falcon tubes for Strain 1 and Strain 2, and stick the dedicated groups of tubes into the respective falcon tube. Dispense 20  $\mu\text{L}$  of worms into Plate 3 and Plate 4. After dispensing the worms, prime 2x with MilliQ water and then prime the ViaFill without MilliQ water so that no more liquid is left in the tubes.
- 19 Place Plate 1 (M9 medium-filled) and Plate 3 (S medium-filled) in the 20°C incubator (static) inside a box with a lid on.
- 20 Wrap Plate 2 (M9 medium-filled) and Plate 4 (S medium-filled) in tissue paper that has been wet with MilliQ water and squeezed. Place the plates in a plastic plate rack and keep them in a 20°C shaking incubator at 140 rpm.
- 21 Leave plates incubating for 76h.

## Day 0.

- 22 Record the plates.