

Dec 18, 2019

© ASSESSING RTTA ACTIVITY (Support Protocol 6)



In 1 collection

DOI

dx.doi.org/10.17504/protocols.io.5w8g7hw



Michael S. Fernandopulle¹, Ryan Prestil¹, Christopher Grunseich¹, Chao Wang², Li Gan², Michael E. Ward¹

¹National Institute of Neurological Disorders and Stroke, National Institutes of Health, Bethesda, Maryland;

²Gladstone Institute of Neurological Disease, Gladstone Institutes, San Francisco, California

Neurodegeneration Method Development Community Tech. support email: ndcn-help@chanzuckerberg.com



Julia Rossmanith

protocols.io

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN ACCESS



DOI: https://dx.doi.org/10.17504/protocols.io.5w8g7hw

External link: https://doi.org/10.1002/cpcb.51



Protocol Citation: Michael S. Fernandopulle, Ryan Prestil, Christopher Grunseich, Chao Wang, Li Gan, Michael E. Ward 2019. ASSESSING RTTA ACTIVITY (Support Protocol 6) . **protocols.io** https://dx.doi.org/10.17504/protocols.io.5w8g7hw

Manuscript citation:

Fernandopulle, M. S., Prestil, R., Grunseich, C., Wang, C., Gan, L., & Ward, M. E. (2018). Transcription-factor mediated differentiation of human iPSCs into neurons. Current Protocols in Cell Biology, e51. doi:https://doi.org/10.1002/cpcb.51

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: July 26, 2019

Last Modified: December 18, 2019

Protocol Integer ID: 26304

Keywords: i3LMN, i3Neurons, iPSC, iPSC-derived neurons, transcription factor-mediated differentiation, assessing rtta activity, rtta activity, support protocol, protocol

Attachments



fernandopulle2018.pd...

1.7MB

Materials

- Cell line of interest and appropriate culture medium
- pBI-MCS-EGFP plasmid (Addgene, cat. no.)

pBI-MCS-EGFP addgene Catalog #16542

- Doxycycline (2 mg/ml in PBS; 1000×; Sigma, cat. no. D9891)
- Doxycycline hyclate Merck MilliporeSigma (Sigma-Aldrich) Catalog #D9891
- Fluorescent microscope
- Additional reagents and equipment for general iPSC culture (<u>Basic Protocol 1</u>) and transfection (<u>Basic Protocol</u> 2)

Troubleshooting



Safety warnings

• Please see SDS (Safety Data Sheet) for hazards and safety warnings.



1 Prepare four wells of the cell line of interest by EDTA split (Basic Protocol 1).

Note

For at least the first time performing this protocol, include positive and negative control cells in order to compare with cells of interest.

- 2 1 to 2 days after plating, follow <u>Basic Protocol 2</u> (steps 6 to 9) for two of the wells in order to transfect with the pBI-MCS-EGFP plasmid.

Note

For AAVS1: Δ 0.75 μ g of pTALdNC-AAVS1_T2 (Addgene, <u>cat. no. 80496</u>) and Δ 0.75 μ g of pTALdNC-AAVS1_T1 (Addgene, <u>cat. no. 80495</u>) per transfection.

For CLYBL: $\underline{\underline{A}}$ 0.75 μg of pZT-C13-R1 (Addgene, <u>cat. no. 62197</u>) and $\underline{\underline{A}}$ 0.75 μg of pZT-C13-L1 (Addgene, <u>cat. no. 62196</u>) per transfection.

- 2.2 Combine the contents of the two tubes and vortex again for 00:00:02 to 00:00:03. Incubate this mixture for 00:10:00 at 8 Room temperature.
- 2.3 Retrieve the cells plated in steps 3 to 4, and, using a P200 pipet tip and repeat pipettor, add $\stackrel{\square}{\underline{}}$ 200 $\stackrel{\square}{\underline{}}$ of the complete transfection solution from step 6 **dropwise**, evenly across the surface of the well. Return the cells to the incubator overnight.
- 2.4 (24:00:00) after transfection, aspirate transfection medium and replace with fresh E8.



If applicable, evaluate transfection efficiency by fluorescence microscopy.

Note

All cells transfected with the hNGN2 (Addgene, <u>cat. no. 105840</u> or <u>110492</u>) and hNIL constructs (Addgene, <u>cat. no. 105841</u> or <u>105842</u>) will transiently express mCherry for **3 to 4 days**, while only those cells with transgene insertion will maintain **stable expression** of mCherry for longer periods of time.

See <u>Basic Protocols 3</u> and <u>4</u> for options for enrichment and clonal isolation.

3 One day later, aspirate medium and replace with fresh E8. In one trans- fected well and one un-transfected well, supplement medium with 2 μg/ml doxycycline.

Note

If cells are more than 30% confluent when transfected, doxycycline may be added at the time of transfection in order to permit imaging the next day.

4 One day later, aspirate medium, rinse twice with PBS, and image with a fluorescent microscope.

Note

eGFP typically has a half-life of 24 hr, so fluorescence should persist for at least 2 to 3 days after doxycycline treatment.

Since transfection is heterogeneous, only a subset of cells should express GFP. In addition, levels of expression may vary by the number of plasmid copies delivered to each cell. However, in cells with active tetracycline transactivator, the average intensity of GFP will be notably higher in the well exposed to doxycycline than in the well which was not. Some leaky expression is expected due to the sequence of the tetracycline operator in this plasmid, so un-transfected cells thus serve as a further negative control.