

Jan 11, 2026

Assessment of Mitochondrial OCR For SH-SY5Y Cell using Agilent Seahorse XFe Cell Mito Stress Test

DOI

<https://dx.doi.org/10.17504/protocols.io.4r3l2qdypl1y/v1>

Ali Ghoochani^{1,2,3,4}, Monther Abu-Remaileh^{1,2,3,4}

¹Department of Chemical Engineering, Stanford University, Stanford, CA 94305, USA;

²Department of Genetics, Stanford University, Stanford, CA 94305, USA;

³The Institute for Chemistry, Engineering and Medicine for Human Health (Sarafan ChEM-H), Stanford University, Stanford, CA 94305, USA;

⁴Aligning Science Across Parkinson's (ASAP) Collaborative Research Network, Chevy Chase, MD, 20815, USA

asap



Ali Ghoochani

Stanford University

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.4r3l2qdypl1y/v1>



Protocol Citation: Ali Ghoochani, Monther Abu-Remaileh 2026. Assessment of Mitochondrial OCR For SH-SY5Y Cell using Agilent Seahorse XFe Cell Mito Stress Test . **protocols.io** <https://dx.doi.org/10.17504/protocols.io.4r3l2qdypl1y/v1>

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: October 13, 2024

Last Modified: January 11, 2026

Protocol Integer ID: 109767

Keywords: Mitochondrial function, Mitochondrial OCR, Seahorse assay, using agilent seahorse xfe cell mito stress test, agilent seahorse xfe cell mito stress test, assessment of mitochondrial function, assessment of mitochondrial ocr, indicator of mitochondrial function, mitochondrial ocr, mitochondrial function, sy5y cell, oxygen species, effects of iron supplementation, oxygen consumption rate, iron supplementation

Abstract

This protocol describes the assessment of mitochondrial function in SH-SY5Y cells using the Agilent Seahorse XFe Cell Mito Stress test, which measures the oxygen consumption rate (OCR) as an indicator of mitochondrial function. Additionally, it evaluates the effects of iron supplementation and reactive oxygen species (ROS) inhibition.



Cell Seeding and Treatment

2d

- 1 Seed 80,000 cells (SH-SY5Y) in 180 μ L media in each well of a Seahorse XFe96 microplate.
Note: Leave well A blank empty as background control.
- 1.1 For iron supplementation, treat cells with 0.2 mg/mL ferric ammonium citrate (FAC)
For ROS inhibition, treat cells with 250 nM Liproxstatin-1
- 1.2 Incubate the plate at 37°C with 5% CO_2 for 48:00:00



2d

Preparation of Sensor Cartridge

- 2 Hydrate the sensor cartridge
- 2.1 The day prior to measuring OCR, hydrate the sensor cartridge with Seahorse XF Calibrant $200\ \mu\text{L}$ /well and $400\ \mu\text{L}$ PBS/ edge well and keep the sensor cartridge at 37°C in a non- CO_2 incubator overnight

Preparing the Assay Medium

1h

- 3 On the day of the assay, replace the growth medium with pre-warmed assay medium ($180\ \mu\text{L}$) containing Seahorse XF base medium, 1 mM pyruvate, and 20 mM glucose.
- 3.1 Incubate the microplate in a non- CO_2 incubator at 37°C for 01:00:00 to equilibrate the temperature and pH.

1h

Preparing Modulating Agents

- 4 In the meantime, prepare following modulating agents (10X) in assay medium for injection:



	A	B	C	D
	Port Loading Design	10X concentration	Add to Port volume	Final well concentration (μM)
	Port A: Oligomycin	10 μM	20 μL	1 μM
	Port B : 2,4-Dinitrophenol (DNP)	500 μM	25 μL	50 μM
	Port C : Rotenone	5 μM	27 μL	0.5 μM

4.1 Load 10X concentration of agents into port of the cartridge (From step 2.1)

Calibrating the Analyzer

15m

- 5 Place the loaded sensor cartridge into analyzer and run the program. The calibration usually takes  00:15:00 .

15m

Measuring the OCR

- 6 After calibration is done, transfer the Seahorse XFe96 cell microplate (from step 3.1) to the analyzer and press "continue"

Post-Assay Procedures

- 7 Protein Normalization:
- 7.1 After the assay, aspirate the assay medium and lyse the cells with a suitable lysis buffer and measure protein concentration using a BCA assay.
- 7.2 Use this data for normalizing the OCR results.