

Dec 04, 2024

ATAC-seq protocol

DOI

dx.doi.org/10.17504/protocols.io.14egn924ml5d/v1

Renhe Luo¹, Michael Beer², Danwei Huangfu¹

¹Sloan Kettering Institute; ²Johns Hopkins University

IGVF



Michael Beer

Johns Hopkins University

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.14egn924ml5d/v1>

Protocol Citation: Renhe Luo, Michael Beer, Danwei Huangfu 2024. ATAC-seq protocol. **protocols.io**
<https://dx.doi.org/10.17504/protocols.io.14egn924ml5d/v1>

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited



Protocol status: Working

We use this protocol and it's working

Created: December 04, 2024

Last Modified: December 04, 2024

Protocol Integer ID: 113699

Keywords: seq protocol atac, seq protocol for esc, atac, seq protocol, esc, seq, protocol, de differentiation

Funders Acknowledgements:

NHGRI

Grant ID: U01HG012051

Abstract

ATAC-seq protocol for ESC-DE differentiation

Troubleshooting



ATAC-seq protocol

- 1 50K cryopreserved cells were washed in cold PBS and lysed.
- 2 The transposition reaction containing TDE1 Tagment DNA Enzyme (Illumina; 20034198) was incubated at 37°C for 30 minutes.
- 3 The DNA was purified with the MinElute PCR Purification Kit (QIAGEN; 28004) and amplified for 5 cycles using NEBNext High-Fidelity 2X PCR Master Mix (New England Biolabs; M0541L).
- 4 After evaluation by real-time PCR, 3-14 additional PCR cycles were done.
- 5 The final product was cleaned by AMPure XP beads (Beckman Coulter; A63882) at a 1X ratio, and size selection was performed at a 0.5X ratio.
- 6 Libraries were sequenced on a HiSeq 4000 or NovaSeq 6000 platform in a PE50 run, using the HiSeq 3000/4000 SBS Kit or NovaSeq 6000 S1 Reagent Kit (100 Cycles) (Illumina).