



Aug 29, 2024

ATG9A-vesicle isolation protocol

 [Nature Cell Biology](#)

DOI

<https://dx.doi.org/10.17504/protocols.io.dm6gpzok8lzp/v1>

Daniel Bernklau¹, Justyna Sawa-Makarska¹, Julia Romanov¹

¹Sascha Martens lab, University of Vienna, Max Perutz Labs - Vienna (AT)



Elias Adriaenssens

Sascha Martens lab, University of Vienna, Max Perutz Labs - ...

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.dm6gpzok8lzp/v1>

External link: <https://doi.org/10.1038/s41556-025-01712-y>

Protocol Citation: Daniel Bernklau, Justyna Sawa-Makarska, Julia Romanov 2024. ATG9A-vesicle isolation protocol. [protocols.io](#) <https://dx.doi.org/10.17504/protocols.io.dm6gpzok8lzp/v1>

Manuscript citation:

Adriaenssens, E., Schaar, S., Cook, A.S.I. *et al.* Reconstitution of BNIP3/NIX-mitophagy initiation reveals hierarchical flexibility of the autophagy machinery. *Nat Cell Biol* **27**, 1272–1287 (2025). <https://doi.org/10.1038/s41556-025-01712-y>



License: This is an open access protocol distributed under the terms of the **[Creative Commons Attribution License](#)**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: July 04, 2024

Last Modified: August 29, 2024

Protocol Integer ID: 103561

Keywords: ASAPCRN, vesicle isolation protocol this protocol detail, vesicle isolation protocol, vesicle isolation, atg9, atg9a, protocol detail, isolation, protocol

Funders Acknowledgements:

Aligning Science Across Parkinson's (ASAP)

Grant ID: ASAP-000350

Marie Skłodowska-Curie MSCA Postdoctoral fellowship

Grant ID: 101062916

Abstract

This protocol details about the ATG9-vesicle isolation.

Materials

Sodium-Orthovanadate (Applichem, prepared as 100 millimolar (mM) in MQ, frozen)

b-Glycerophosphate (Sigma, prepared as 1 Molarity (M) in MQ, frozen)

NaF (Sigma, prepared as 500 millimolar (mM) in MQ, frozen)

Complete EDTA free protease inhibitor tablets, Sigma 5056489001 (prepared as 50x stock 1 tbl in 1 ml H2O)

Anti-FLAG M2 Affinity Gel from Sigma A2220-1ml

D(+)-Saccharose (sucrose) from Roth 4661.1

FLAG Peptide from Sigma F32920-4MG

Dissolved at 4 mg /ml in SEC buffer (25 millimolar (mM) NaCl, 1 millimolar (mM) DTT), aliquoted, kept at -20 °C

3 mL Luer Lock HENKE-JECT syringe

Braun Sterican 26 G x 1'' Gr. 18 needle

MilliQ water (filtered additionally through a 0.2µm membrane)

All buffers used were filtered prior to use

0.2 µm celluloseacetate syringe filter from Chromafil CA-20/25 (S)

	A	B	C	D	E	F
			Vesicle Isolation Buffer		Elution Base Buffer	
	Reagent	Stock conc	Final conc	Vol for 20ml	Final Conc	Vol for 20ml
	HEPES pH7.5	1M	20mM	400µl	20mM	400µl
	NaCl	5M	150mM	600µl	150mM	600µl
	Sucrose	1M	250mM	5ml	-	-
	Complete EDTA free protease inhibitors (PI)	50x	1x	400µl	1x	400µl
	Beta-Glycerophosphate	1M	20mM	400µl	20mM	400µl
	Sodium-Orthovanadate	100mM	1mM	200µl	1mM	200µl
	NaF	500mM	1mM	40µl	1mM	40µl
	EDTA pH8.0	0.5M	1mM	40µl	-	-
	In MilliQ H2O					

	A	B	C	D	E	F
	Sterile filtered, precooled					

Before start

Hap1 cells were previously CRISPR engineered in our lab to contain ATG9A endogenously tagged on the C-terminus with mEGFP-3C-FLAG.



Preparation of starting material (Hap1 ATG9A-mEGFP-3C-FLAG cell pellets)

- 1 Grow and expand cells in IMDM medium supplied with 10% FCS and Pen/Strep until the required number of cells is reached.
Prep size: For one preparation of vesicles the required amount is 100×10^6 to 150×10^6 cells which is equivalent to 4 - 5 $\times$ 15 cm dishes of confluent cells.
Detach cells with trypsin, count and aliquot to prep size (see above), wash with cold PBS, spin down  200 rcf, 4°C , remove supernatant, flash freeze the cell pellets and keep at  -80 °C until use.

2w

Day1: Cell lysis and binding to Flag beads

1h

- 2 **Equilibrate Flag beads**
Always pellet the beads with  4000 rcf, 4°C cooled 5415R centrifuge for 3'. 
- 2.1 Take  70 μ L of FLAG bead slurry per prep.
- 2.2 Wash 3x with MilliQ water 
- 2.3 Wash 3x with Vesicle Isolation Buffer. 
- 2.4 Prepare a 1 : 1 slurry in **Vesicle isolation Buffer**
- 3 **Cell lysis**
Always work  On ice 
">> sample name" means: take a sample for Western Blot e.g. 10 μ l sample + 2 μ l 6xProtein loading dye
- 3.1 Thaw cells  On ice
- 3.2 **Resuspend** the pellet in  1665 μ L of **Vesicle Isolation Buffer** and transfer to a fresh 2ml Eppi. 



- 3.3 Incubate On ice for 00:20:00 20m
- 3.4 Lyse the cells by pushing them **30 times (up and down) through a 26 x g needle** 10m
- 3.5 Let them chill for 00:10:00 On ice 10m
- 3.6 Repeat the lysis step (**push 30 times up and down through a 26 x g needle**) 10m
- 3.7 6m
Pellet the cell debris and nuclei (200 rcf, 4°C, 00:06:00)
- 3.8 Transfer the supernatant to a fresh 2 ml Eppi (>> Sample Supernatant)
- 4 Add 70 µL of equilibrated **FLAG beads** (slurry). 10m
Incubate Overnight at 4 °C (on the roller in the cold room in a 50ml Falcon tube)

Day2 Elution of ATG9A vesicles

3h 15m

- 5 Pellet the beads (in a 5415R centrifuge) at 200 rcf, 00:03:00 and remove the supernatant (>> Sample Unbound). 3m
- 6 Resuspend the beads in 1 mL of **Vesicle Isolation Buffer** and transfer to a fresh 2ml Eppi (+ 665 µL **Vesicle Isolation Buffer** already pre-pipetted in the Eppi).
- 7 Incubate 00:01:00 rolling at 4 °C or on ice with occasionally inverting the Eppi carefully 1m
- 8 Pellet the beads at 2000 rcf, 00:03:00 and remove the supernatant 3m



- 9 Resuspend the beads in  1 mL of Vesicle Isolation Buffer and transfer into a fresh 2 ml LowBind tube 
- 9.1 Add  665 μL of **Elution Base Buffer** very slowly to the tube (decreases the sucrose concentration; now be even more careful than before!) 
- 10 Incubate  00:01:00 rolling at  4 °C or on ice with occasionally inverting the tube carefully 

- 11 Pellet the beads at  2000 rcf, 00:03:00 and remove the supernatant 

- 12 Resuspend the beads in  1 mL of **Elution Base Buffer** and transfer into a fresh 2 ml LowBind tube  665 μL **Elution Base Buffer** already in pre-pipetted in the tube). 
- 13 Incubate  00:01:00 rolling at  4 °C or on ice with occasionally inverting the tube carefully 

- 14 Pellet the beads at  2000 rcf, 00:03:00 and remove the supernatant. 

- 15 Resuspend the beads in  666 μL of **Elution Base Buffer** and transfer to a fresh 1.5 ml LowBind tube (>> Sample Elution Input). 
- 16 Add  16.65 μL of  4 mg /ml **FLAG peptide**
- 17 Incubate for  02:45:00 at 4°C rolling 

- 18 Pellet the beads at  2000 rcf, 00:03:00 and transfer the supernatant to a fresh LowBind tube (>> Sample Elution) 
This is your vesicle prep! Handle with care. 



- 19 Resuspend the beads again in  666 μL of **Elution Base Buffer** (>>Sample Beads after Elution)