

Mar 11, 2020


Version 2

🌐 ATR FTIR Spectroscopy of Aqueous Cell Cultures V.2

📁 In 1 collection

DOI

<https://dx.doi.org/10.17504/protocols.io.bdidi4a6>



Emrulla Spahiu¹, Senol Dogan², Jörg chnauß², Mayda Gursel¹, Feride Severcan³

¹Middle East Technical University;

²Leipzig University, Faculty of Physics and Earth Sciences, Peter Debye Institute for Soft Matter Physics;

³Department of Biophysics, Faculty of Medicine, Altinbas University



Emrulla Spahiu

Middle East Technical University

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.bdidi4a6>

Protocol Citation: Emrulla Spahiu, Senol Dogan, Jörg chnauß, Mayda Gursel, Feride Severcan 2020. ATR FTIR Spectroscopy of Aqueous Cell Cultures. **protocols.io** <https://dx.doi.org/10.17504/protocols.io.bdidi4a6>



License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: March 11, 2020

Last Modified: March 11, 2020

Protocol Integer ID: 34085

Keywords: atr, ftir, infrared, spectroscopy, aqueous, cell culture, data processing, PCA, spectra, ftir spectroscopy of aqueous cell culture sample, aqueous cell cultures atr, aqueous cell culture sample, ftir spectroscopy, infrared spectra acquisition,

Abstract

ATR-FTIR Spectroscopy of aqueous cell culture samples is explained step-by-step. Infrared spectra acquisition, processing and analysis is included briefly.



Spectra acquisition

- 1 Arrange the temperature of the working environment to 🌡️ 25 °C
- 2 Clean the ATR crystal using a cellulose wipe and water-ethanol-water sequence. Scratching the crystal surface must be avoided.

Equipment

ATR FTIR

NAME

Spectroscope

TYPE

Perkin Elmer Inc., Norwalk, CT, USA

BRAND

X

SKU

- 3 Leave the crystal to air-dry
- 4 Select "Scan and Instrument Setup" and input sample name, scan range (650-4000 cm^{-1}), scan number (128), and spectral resolution (2 cm^{-1}). This step needs optimization for specific sample being used.

Software

Spectrum One

NAME

Perkin Elmer

DEVELOPER

- 5 Scan the background spectrum to overcome the probable atmospheric interference that can be triggered by H_2O and CO_2 molecules present in the air, subtracted automatically by the software.



- 6 Additionally, record PBS as sample (subtract this manually from each spectrum)
This eliminates interference of PBS and flattens the water band located around 2125 cm⁻¹
- 7 Put  1 µL cell samples and let to dry with a mild N₂ flux for  00:05:00 . This allows cells to settle on the crystal. 5m
- 8 Scan the sample. Process takes about  00:07:00 to finish. 7m
As the scanning process goes on, one must be sure that the quality check sign is everytime green.
- 9 Scan each sample under the same conditions in three independent scans (the average constitutes a replicate n). Replicates (n=5) are used in data analysis.
This minimizes intra-sample variability and eliminates possible variation that might arise from experimental conditions.

Spectral pre-processing

- 10 The average of triplicates is taken manually.

Software

Spectrum One

NAME

Perkin Elmer

DEVELOPER

- 11 13-point Savitzky-Golay smoothing
- 12 Baseline correction at 3800-2750-1800-900 cm⁻¹points.

Principal Component Analysis (PCA)

- 13 Import spectra as data file into



Software

Unscrambler

NAME

Camo, NO

DEVELOPER

14 Define the data sets to 3000-650 cm^{-1} region.

15 Make the PCA Model.
Define the regions with high variation using the Loadings plots.
Check how the replicates are grouping and if there are outliers in Scores plots.

Spectral Calculations

16 Peak labelling can be done automaticaly. For increased precision we measured manually:
Measure the center of weight of each band. Read the frequency.
Compare the peaks you obtained, to the literature:

Citation

Zanyar Movasaghi, Shazza Rehman, and Ihtesham ur Rehman (2026)
. Fourier Transform Infrared (FTIR) Spectroscopy of Biological Tissues.
Applied Spectroscopy Reviews.

[10.1080/05704920701829043](https://doi.org/10.1080/05704920701829043)

LINK

Make an excel file of your peaks and their assignment to the functional groups.
For side-peak band assignment, the second derivative spectrum is used.

17 Concentration studies are carried out by integrating the curves under the peaks.

18 For conformational alteration studies, bandwidth measurements are done manually using the 75% height of the peaks.



Citations

Step 16

Zanyar Movasaghi, Shazza Rehman, and Ihtesham ur Rehman. Fourier Transform Infrared (FTIR) Spectroscopy of Biological Tissues

10.1080/05704920701829043