

Dec 19, 2019

Attachment 1: Preparation of Fluorescein Diacetate and Propidium Iodide Stock Solutions (FDA/PI)

DOI

dx.doi.org/10.17504/protocols.io.5bgg2jw



Integrated Islet Distribution Program¹

¹[Integrated Islet Distribution Program, City of Hope, Duarte, CA]

Integrated Islet Distribution Program
Tech. support email: iidp-email@coh.org



Integrated Islet Distribution Program

Integrated Islet Distribution Program, City of Hope

OPEN  ACCESS



DOI: dx.doi.org/10.17504/protocols.io.5bgg2jw

Protocol Citation: Integrated Islet Distribution Program 2019. Attachment 1: Preparation of Fluorescein Diacetate and Propidium Iodide Stock Solutions (FDA/PI). **protocols.io** <https://dx.doi.org/10.17504/protocols.io.5bgg2jw>

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: July 09, 2019

Last Modified: June 11, 2020

Protocol Integer ID: 25672

Keywords: Inclusion and Exclusion Fluorescent Dyes, Fluorescein Diacetate, Propidium Iodide, Viability Assay



Abstract

This Standard Operating Procedure is adapted from the work of the '*National Institutes of Health-Sponsored Clinical Islet Transplantation Consortium Phase 3 Trial: Manufacture of a Complex Cellular Product at Eight Processing Facilities*' following the SOP cited in the document '*Purified Human Pancreatic Islet - Viability Estimation of Islet Using Fluorescent Dyes (FDA/PI): Standard Operating Procedure of the NIH Clinical Islet Transplantation Consortium*'

This SOP defines the procedure for assessment of viability of human isolated islet preparations, which include endocrine and exocrine tissue, for use in the National Institute of Diabetes and Digestive and Kidney Diseases (NIDDK) sponsored research in the Integrated Islet Distribution Program (IIDP). This protocol is written to assist the participating islet isolation centers and investigators who are part of this program.

Fluorescein Diacetate/ Propidium Iodide (FDA)/(PI) Viability Assay is a rapid fluorometric method to test the integrity of the plasma membrane simultaneously using inclusion and exclusion dyes; the assay differentiates between viable and nonviable cells and is, consequently, used for determination of viability of islet preparations.

Guidelines

Fluorescein Diacetate/ Propidium Iodide (FDA)/(PI) Viability Assay is a rapid fluorometric method to test the integrity of the plasma membrane simultaneously using inclusion and exclusion dyes; the assay differentiates between viable and nonviable cells and is, consequently, used for determination of viability of islet preparations.

- The inclusion dye is Fluorescein Diacetate (FDA) and the exclusion dye is Propidium Iodide (PI). The final concentrations are as follows:
 - ◆ FDA: 0.46 μ M
 - ◆ PI: 14.34 μ M
- Fluorescein Diacetate is a nonpolar ester, which passes through plasma membranes and is hydrolyzed by intracellular esterases to produce free fluorescein. The polar fluorescein is confined within cells with an intact plasma membrane and can be observed under appropriate excitation conditions. FDA functions as an inclusion dye, i.e., viable cells will appear bright green fluorescent using FDA.
- Propidium iodide functions as an exclusion dye that cannot penetrate living cells but readily enters dead or dying cells. Once PI penetrates through the cell membrane, it binds to nucleic acids and causes them to fluoresce bright orange/red. PI absorbs in green light and fluoresces orange/red.
- Both of the fluorescent dyes used in this assay are light sensitive and must be kept in the dark, covered with aluminum foil.
- The fluorescent dyes are temperature sensitive and must be stored as follows:
 - ◆ FDA: $\leq -20^{\circ}\text{C}$
 - ◆ PI: $2-8^{\circ}\text{C}$

Materials


MATERIALS


Fluorescein Diacetate
Merck MilliporeSigma (Sigma-Aldrich) Catalog #F7378


Propidium Iodide
Merck MilliporeSigma (Sigma-Aldrich) Catalog #P4170


Acetone
Merck MilliporeSigma (Sigma-Aldrich) Catalog #179124


Thermo Scientific™ Nunc™ Cell Culture/Petri Dishes
Fisher Scientific Catalog #174926



Corning® Dulbecco's Phosphate-Buffered Saline 1X without calcium and magnesium
Corning Catalog #21-031-CM

Equipment

OHAUS™ Explorer™ Semi-Micro Balance or equivalent	NAME
Semi-Micro Balance	TYPE
Ohaus	BRAND
01-919-370	SKU
https://www.fishersci.com/shop/products/ohaus-explorer-semi-micro-balances/p-6541003#	LINK
0.01 mg to 52 gm capacity	SPECIFICATIONS



Equipment	
Snap Cap Microcentrifuge Tube or equivalent	NAME
Polypropylene Microcentrifuge Tube	TYPE
Corning Costar Snap Cap Microcentrifuge Tube	BRAND
07200210	SKU
https://www.fishersci.com/shop/products/costar-microcentrifuge-tubes-6/07200210	LINK
2 mL snap cap polypropylene micro tube	SPECIFICATIONS

Equipment	
5 mL Transport Tubes	NAME
5 mL Transport Tubes	TYPE
Globe Scientific Self-Standing Transport Tubes or	BRAND
22-010-1227	SKU
https://www.fishersci.com/shop/products/self-standing-transport-tubes-separate-screw-cap-7/p-7112488#?keyword=5+ML+CENTRIFUGE+TUBE+WITH+CAP	LINK
5 mL Polypropylene tube with Polyethylene Screw cap	SPECIFICATIONS
	

Equipment

Fisherbrand™ Pipette-Specific Tips	NAME
Pipette Tips, 100 to 1000uL or equivalent	TYPE
Pipet Tips	BRAND
02-681-182	SKU
https://www.fishersci.com/shop/products/fisherbrand-pipetter-specific-tips-natural-100-1000-l-101-6mm-long-bulk-pack/02681182#?keyword=true	LIN K
100 to 1000uL	SPECIFICATIONS



Safety warnings

Always wear gloves and observe standard chemical procedures:

Fluorescein Diacetate: FDA MSDSAction.pdf

- Protect from light. Avoid contact and inhalation. Nitrile gloves are recommended in the MSDS when handling FDA.

Propidium Iodide: PI MSDSAction.pdf

- Use personal protective equipment. Product may be toxic if inhaled, swallowed, or splashed on skin. Avoid dust formation. Avoid breathing vapors, mist, or gas. Ensure adequate ventilation. Wear gloves and observe Safety Data Sheet. Suspected of causing genetic defects.

Acetone: AcetoneMSDSAction.pdf

- Solvent/Flammable. Keep away from heat, spark, and open flame. Keep container tightly closed. Use with adequate ventilation. Avoid contact with eyes. Avoid prolonged or repeated breathing of vapor. Handle with gloves. Gloves must be inspected prior to use. Use proper glove removal technique (without touching glove's outer surface) to avoid skin contact with this product. Dispose of contaminated gloves after use in accordance with applicable laws and good laboratory practices. Wash and dry hands. Store aliquots in tightly sealed in glass or polypropylene tubes with polypropylene or polyethylene closure.

Preparation of Stock Fluorescein Diacetate (FDA) Solution

- 1 Prepare **1M** 24 micromolar (μM) Fluorescein Diacetate (FDA) according to the formulation in sub-steps below.

Note

The expiration date, for both PI and FDA stains, is six months from the date of preparation.

Both of the fluorescent dyes used in this assay are light sensitive and must be kept in the dark, covered with aluminum foil.

1.1 Calculation:

FDA FW = 416.4

Stock Concentration = 24 μM

Volume required = 200 mL

Formula: FW X Concentration X Volume = $(416.4) \times (24 \times 10^{-6}) \times (200 \times 10^{-3}) = 0.00199$ g FDA

Sigma Aldrich has a Mass Molarity Calculator that can be used to determine smaller quantities: <https://www.sigmaaldrich.com/chemistry/stockroom-reagents/learning-center/technical-library/mass-molarity-calculator.html>

- 1.2 Dissolve 0.00199 g of FDA in 200 mL of acetone in a glass bottle and cover with aluminum foil.

0.00199 g FDA

- FDA Supplier _____
- FDA Lot # _____
- Expiration Date _____

200 mL Acetone

- Acetone Supplier _____
- Acetone Lot # _____
- Expiration Date _____

- 1.3
 - Store Stock FDA solution tightly sealed in 4 mL or smaller aliquots in polypropylene tubes or glass tubes at \leq -20 °C , **protected from light**, for up to six months.



- Record the following on each aliquot tube:

FDA (24 µM) Stock Solution Assigned Lot number: _____

Expiration Date: _____

Prepared By/Date: _____

Reviewed By/Date: _____

Preparation of Stock Propidium Iodide (PI)

- Prepare **1M** 750 micromolar (µM) Propidium Iodide (PI) according to the formulation in sub-steps below.

2.1 Calculation:

PI FW = 668.4

Stock Concentration = 750 µM

Volume required = 25 mL

$$\text{FW} \times \text{Concentration} \times \text{Volume} = (668.4) \times (750 \times 10^{-6}) \times (25 \times 10^{-3}) = 0.0125 \text{ g PI}$$



- Dissolve 0.0125 g of PI in 25 mL of DPBS and cover with aluminum foil.

 0.0125 g Propidium Iodide

- PI Supplier _____
- PI Lot # _____
- Expiration Date _____

 25 mL DPBS

- PBS Supplier _____
- PBS Lot # _____
- Expiration Date _____

- Store Stock PI solution tightly sealed in 0.5 mL or smaller aliquots in polypropylene Microcentrifuge Snap Cap tubes or glass tubes at  2 °C to  8 °C , **protected from light**, for up to six months.

- Record the following on each aliquot tube:

PI (750 µM) Stock Solution Assigned Lot number: _____

Expiration Date: _____

Prepared By/Date: _____

Reviewed By/Date: _____