

Dec 18, 2019

## Attachment 1: Preparation of Fluorescein Diacetate and Propidium Iodide Stock Solutions (FDA/PI)

DOI

[dx.doi.org/10.17504/protocols.io.5bgg2jw](https://dx.doi.org/10.17504/protocols.io.5bgg2jw)



Integrated Islet Distribution Program<sup>1</sup>

<sup>1</sup>[Integrated Islet Distribution Program, City of Hope, Duarte, CA]

Integrated Islet Distribution Program  
Tech. support email: [iidp-email@coh.org](mailto:iidp-email@coh.org)



Integrated Islet Distribution Program

Integrated Islet Distribution Program, City of Hope

### Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.5bgg2jw>

**Protocol Citation:** Integrated Islet Distribution Program 2019. Attachment 1: Preparation of Fluorescein Diacetate and Propidium Iodide Stock Solutions (FDA/PI). **protocols.io** <https://dx.doi.org/10.17504/protocols.io.5bgg2jw>

**License:** This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

**Protocol status:** Working

**We use this protocol and it's working**

**Created:** July 09, 2019

**Last Modified:** June 11, 2020

**Protocol Integer ID:** 25672

**Keywords:** Inclusion and Exclusion Fluorescent Dyes, Fluorescein Diacetate, Propidium Iodide, Viability Assay, nih clinical islet transplantation consortium, sponsored clinical islet transplantation consortium phase, integrated islet distribution program, research in the integrated islet distribution program, purified human pancreatic islet, viability estimation of islet, isolated islet preparation, islet preparation, viability assay, fluorescein diacetate, preparation of fluorescein diacetate, participating islet isolation center, islet isolation center, assay differentiate, fda, propidium iodide stock solution, assay, islet, diabetes, using fluorescent dye, complex cellular product, nonviable cell, integrity of the plasma membrane, national institutes of health

## Abstract

This Standard Operating Procedure is adapted from the work of the '*National Institutes of Health-Sponsored Clinical Islet Transplantation Consortium Phase 3 Trial: Manufacture of a Complex Cellular Product at Eight Processing Facilities*' following the SOP cited in the document '*Purified Human Pancreatic Islet - Viability Estimation of Islet Using Fluorescent Dyes (FDA/PI): Standard Operating Procedure of the NIH Clinical Islet Transplantation Consortium*'

This SOP defines the procedure for assessment of viability of human isolated islet preparations, which include endocrine and exocrine tissue, for use in the National Institute of Diabetes and Digestive and Kidney Diseases (NIDDK) sponsored research in the Integrated Islet Distribution Program (IIDP). This protocol is written to assist the participating islet isolation centers and investigators who are part of this program.

Fluorescein Diacetate/ Propidium Iodide (FDA)/(PI) Viability Assay is a rapid fluorometric method to test the integrity of the plasma membrane simultaneously using inclusion and exclusion dyes; the assay differentiates between viable and nonviable cells and is, consequently, used for determination of viability of islet preparations.

## Guidelines

**Fluorescein Diacetate/ Propidium Iodide (FDA)/(PI) Viability Assay** is a rapid fluorometric method to test the integrity of the plasma membrane simultaneously using inclusion and exclusion dyes; the assay differentiates between viable and nonviable cells and is, consequently, used for determination of viability of islet preparations.

- The inclusion dye is Fluorescein Diacetate (FDA) and the exclusion dye is Propidium Iodide (PI). The final concentrations are as follows:
  - ◆ FDA: 0.46  $\mu$ M
  - ◆ PI: 14.34  $\mu$ M
- Fluorescein Diacetate is a nonpolar ester, which passes through plasma membranes and is hydrolyzed by intracellular esterases to produce free fluorescein. The polar fluorescein is confined within cells with an intact plasma membrane and can be observed under appropriate excitation conditions. FDA functions as an inclusion dye, i.e., viable cells will appear bright green fluorescent using FDA.
- Propidium iodide functions as an exclusion dye that cannot penetrate living cells but readily enters dead or dying cells. Once PI penetrates through the cell membrane, it binds to nucleic acids and causes them to fluoresce bright orange/red. PI absorbs in green light and fluoresces orange/red.
- Both of the fluorescent dyes used in this assay are light sensitive and must be kept in the dark, covered with aluminum foil.
- The fluorescent dyes are temperature sensitive and must be stored as follows:
  - ◆ FDA:  $\leq -20^{\circ}\text{C}$
  - ◆ PI:  $2-8^{\circ}\text{C}$

# Materials

## MATERIALS

- Fluorescein Diacetate **Merck MilliporeSigma (Sigma-Aldrich) Catalog #F7378**
- Propidium Iodide **Merck MilliporeSigma (Sigma-Aldrich) Catalog #P4170**
- Acetone **Merck MilliporeSigma (Sigma-Aldrich) Catalog #179124**
- Thermo Scientific™ Nunc™ Cell Culture/Petri Dishes **Fisher Scientific Catalog #174926**
- Corning® Dulbecco's Phosphate-Buffered Saline 1X without calcium and magnesium **Corning Catalog #21-031-CM**

## Equipment

|   |                |
|---|----------------|
| OHAUS™ Explorer™ Semi-Micro Balance or equivalent   | NAME           |
| Semi-Micro Balance  | TYPE           |
| Ohaus   | BRAND          |
| 01-919-370  | SKU            |
| <a href="https://www.fishersci.com/shop/products/ohaus-explorer-semi-micro-balances/p-6541003#">https://www.fishersci.com/shop/products/ohaus-explorer-semi-micro-balances/p-6541003#</a> | LINK           |
| 0.01 mg to 52 gm capacity   | SPECIFICATIONS |



## Equipment

### Snap Cap Microcentrifuge Tube or equivalent

NAME

Polypropylene Microcentrifuge Tube

TYPE

Corning Costar Snap Cap Microcentrifuge Tube

BRAND

07200210

SKU

<https://www.fishersci.com/shop/products/costar-microcentrifuge-tubes-6/07200210>

LINK

2 mL snap cap polypropylene micro tube

SPECIFICATIONS

## Equipment

### 5 mL Transport Tubes

NAME

5 mL Transport Tubes

TYPE

Globe Scientific Self-Standing Transport Tubes or

BRAND

22-010-1227

SKU

<https://www.fishersci.com/shop/products/self-standing-transport-tubes-separate-screw-cap-7/p-7112488#?keyword=5+ML+CENTRIFUGE+TUBE+WITH+CAP>

LI  
NK

5 mL Polypropylene tube with Polyethylene Screw cap

SPECIFICATIONS



Equipment

|   |                |
|---|----------------|
| Fisherbrand™ Pipette-Specific Tips  | NAME           |
| Pipette Tips, 100 to 1000uL or equivalent   | TYPE           |
| Pipet Tips  | BRAND          |
| 02-681-182  | SKU            |
| <a href="https://www.fishersci.com/shop/products/fisherbrand-pipetter-specific-tips-natural-100-1000-l-101-6mm-long-bulk-pack/02681182#?keyword=true">https://www.fishersci.com/shop/products/fisherbrand-pipetter-specific-tips-natural-100-1000-l-101-6mm-long-bulk-pack/02681182#?keyword=true</a> | LIN<br>K       |
| 100 to 1000uL   | SPECIFICATIONS |



Troubleshooting

## Safety warnings

### **Always wear gloves and observe standard chemical procedures:**

#### **Fluorescein Diacetate:** FDA MSDSAction.pdf

- Protect from light. Avoid contact and inhalation. Nitrile gloves are recommended in the MSDS when handling FDA.

#### **Propidium Iodide:** PI MSDSAction.pdf

- Use personal protective equipment. Product may be toxic if inhaled, swallowed, or splashed on skin. Avoid dust formation. Avoid breathing vapors, mist, or gas. Ensure adequate ventilation. Wear gloves and observe Safety Data Sheet. Suspected of causing genetic defects.

#### **Acetone:** AcetoneMSDSAction.pdf

- Solvent/Flammable. Keep away from heat, spark, and open flame. Keep container tightly closed. Use with adequate ventilation. Avoid contact with eyes. Avoid prolonged or repeated breathing of vapor. Handle with gloves. Gloves must be inspected prior to use. Use proper glove removal technique (without touching glove's outer surface) to avoid skin contact with this product. Dispose of contaminated gloves after use in accordance with applicable laws and good laboratory practices. Wash and dry hands. Store aliquots in tightly sealed in glass or polypropylene tubes with polypropylene or polyethylene closure.

## Preparation of Stock Fluorescein Diacetate (FDA) Solution

- 1 Prepare **1M** 24 micromolar ( $\mu\text{M}$ ) Fluorescein Diacetate (FDA) according to the formulation in sub-steps below.

### Note

*The expiration date, for both PI and FDA stains, is six months from the date of preparation.*

*Both of the fluorescent dyes used in this assay are light sensitive and must be kept in the dark, covered with aluminum foil.*

### 1.1 Calculation:

FDA FW = 416.4

Stock Concentration = 24  $\mu\text{M}$

Volume required = 200 mL

Formula: FW X Concentration X Volume =  $(416.4) \times (24 \times 10^{-6}) \times (200 \times 10^{-3}) = 0.00199$  g FDA

Sigma Aldrich has a Mass Molarity Calculator that can be used to determine smaller quantities: <https://www.sigmaaldrich.com/chemistry/stockroom-reagents/learning-center/technical-library/mass-molarity-calculator.html>

- 1.2 Dissolve 0.00199 g of FDA in 200 mL of acetone in a glass bottle and cover with aluminum foil.

0.00199 g FDA

- FDA Supplier \_\_\_\_\_
- FDA Lot # \_\_\_\_\_
- Expiration Date \_\_\_\_\_

200 mL Acetone

- Acetone Supplier \_\_\_\_\_
- Acetone Lot # \_\_\_\_\_
- Expiration Date \_\_\_\_\_

- 1.3
  - Store Stock FDA solution tightly sealed in 4 mL or smaller aliquots in polypropylene tubes or glass tubes at  $\leq$  -20  $^{\circ}\text{C}$  , **protected from light**, for up to six months.





- Record the following on each aliquot tube:

FDA (24 µM) Stock Solution Assigned Lot number: \_\_\_\_\_

Expiration Date: \_\_\_\_\_

Prepared By/Date: \_\_\_\_\_

Reviewed By/Date: \_\_\_\_\_

## Preparation of Stock Propidium Iodide (PI)

- 2 Prepare **1M** 750 micromolar (µM) Propidium Iodide (PI) according to the formulation in sub-steps below.

### 2.1 Calculation:

PI FW = 668.4

Stock Concentration = 750 µM

Volume required = 25 mL

$$\text{FW} \times \text{Concentration} \times \text{Volume} = (668.4) \times (750 \times 10^{-6}) \times (25 \times 10^{-3}) = 0.0125 \text{ g PI}$$

- 2.2 Dissolve 0.0125 g of PI in 25 mL of DPBS and cover with aluminum foil.

 0.0125 g Propidium Iodide

- PI Supplier \_\_\_\_\_
- PI Lot # \_\_\_\_\_
- Expiration Date \_\_\_\_\_

 25 mL DPBS

- PBS Supplier \_\_\_\_\_
- PBS Lot # \_\_\_\_\_
- Expiration Date \_\_\_\_\_

- 2.3
- Store Stock PI solution tightly sealed in 0.5 mL or smaller aliquots in polypropylene Microcentrifuge Snap Cap tubes or glass tubes at **2 °C** to **8 °C** , **protected from light**, for up to six months.

- Record the following on each aliquot tube:

PI (750 µM) Stock Solution Assigned Lot number: \_\_\_\_\_

Expiration Date: \_\_\_\_\_

Prepared By/Date: \_\_\_\_\_

Reviewed By/Date: \_\_\_\_\_