


Jun 19, 2023

Autofluorescence Microscopy QC for Multimodal Molecular Imaging

DOI

dx.doi.org/10.17504/protocols.io.e6nvwm5dlmk/v1



Katerina V Djambazova¹, Nathan Heath Patterson¹, Lukasz Migas², Angela R.S. Kruse¹, Melissa Farrow¹, Raf Van De Plas², Jeff Spraggins¹

¹Vanderbilt University; ²Delft University of Technology

VU Biomolecular Multim...

KPMP



Katerina V Djambazova

Vanderbilt University

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.e6nvwm5dlmk/v1>

Protocol Citation: Katerina V Djambazova, Nathan Heath Patterson, Lukasz Migas, Angela R.S. Kruse, Melissa Farrow, Raf Van De Plas, Jeff Spraggins 2023. Autofluorescence Microscopy QC for Multimodal Molecular Imaging . **protocols.io**

<https://dx.doi.org/10.17504/protocols.io.e6nvwm5dlmk/v1>



License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: February 08, 2023

Last Modified: October 18, 2023

Protocol Integer ID: 76686

Keywords: autofluorescence microscopy qc for multimodal molecular imaging, autofluorescence microscopy qc, autofluorescence qc, fluorescence, fluorescence intensity, fluorescence channel, multimodal molecular imaging, microscope performance, imaging, tissue af, egfp

Abstract

Autofluorescence QC is performed prior to collecting on-tissue AF. Fluorescence intensity can be tracked for each fluorescence channel (DAPI, EGFP, and DsRed), allowing the microscope performance and the LED light source to be benchmarked and tracked over time.

Materials

StarLight Calibration Slides (Bang Laboratories) - Glacial Blue (SL1GB), Dragon Green (SL1DG), Envy Green (SL1EG)

Zeiss AxioScan Z1 Microscope

Troubleshooting

Autofluorescence Microscopy QC

- 1 Fluorescence QC is collected prior to any on-tissue AF collection.

Use fluorescence calibration slides - StarLight Calibration Slides - Bang Laboratories DAPI (excitation-360nm, emission-450nm), Dragon green (excitation-480nm, emission-520nm), and Envy green (excitation-525nm, emission-565nm)

- 2 Load calibration slides into Zeiss AxioScan.Z1 slide scanner

- 3 Select the same ~ 25 × 25 mm square region in each calibration slide

Note

To reproducibly sample the exact region each time, the square region can be saved as a "region of interest", which can be imported for each slide.

- 4 Acquire fluorescence microscopy data for each calibration slide.

Note

Example set-up and relative ranges:

DAPI QC: DAPI (ex. 360 nm, em. 450 nm) Calibration Slide

LED Light: 385 nm
Light Source Intensity: 90%
Exposure time: 20 ms

eGFP QC: Dragon Green (ex. 480 nm, em. 520 nm) Calibration Slide

LED Light: 475 nm
Light Source Intensity: 90%
Exposure time: 60 ms

DsRed QC: Envy Green (ex. 525 nm, em. 565 nm) Calibration Slide

LED Light: 567 nm
Light Source Intensity: 90%
Exposure time: 250 ms

- 5 Data Analysis:

Once data are collected, background pixels are subtracted, and the average intensity of each channel is calculated and monitored over time.