




Aug 28, 2024

Version 1

🌐 Automated immunodetection of proteinase-K-resistant Lewy Pathology in FFPE human colon tissue V.1

DOI

dx.doi.org/10.17504/protocols.io.kxygxy8ekl8j/v1



Elizabeth Videlock¹

¹UCLA

Neurodegeneration Method Development Community
Tech. support email: ndcn-help@chanzuckerberg.com



Elizabeth Videlock

UCLA

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.kxygxy8ekl8j/v1>

Protocol Citation: Elizabeth Videlock 2024. Automated immunodetection of proteinase-K-resistant Lewy Pathology in FFPE human colon tissue. **protocols.io** <https://dx.doi.org/10.17504/protocols.io.kxygxy8ekl8j/v1>

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: August 28, 2024

Last Modified: August 28, 2024

Protocol Integer ID: 106608

Keywords: gut lewy pathology, tissue resource for parkinsons disease, resistant lewy pathology, arizona parkinsons disease consortium, p30 ag19610 arizona alzheimers disease core center, translational pathology core laboratory, resistant lewy pathology in ffpe, parkinsons research, human colon tissue automated detection, parkinsons disease, automated immunodetection of proteinase, arizona alzheimers research center, human colon tissue, provision of human colon tissue, banner sun health research institute brain, automated immunodetection, arizona biomedical research commission, leica biosystem, body donation program, proteinase, tissue resource, tissue processing

Funders Acknowledgements:

Chan Zuckerberg Initiative

Grant ID: 2021-235154

Abstract

Automated detection of gut Lewy Pathology using Leica Bond RX processor based on Protocol F using Bond Polymer Refine Detection kit (Leica Biosystems, Cat#: DS9800). Adapted from "Method 5" as described in Beach TG et al., J Parkinsons Dis. 2016.


Citation

Beach TG, Corbillé AG, Letournel F, Kordower JH, Kremer T, Munoz DG, Intorcchia A, Hentz J, Adler CH, Sue LI, Walker J, Serrano G, Derkinderen P (2016)

. Multicenter Assessment of Immunohistochemical Methods for Pathological Alpha-Synuclein in Sigmoid Colon of Autopsied Parkinson's Disease and Control Subjects..

<https://doi.org/>

LINK



a. Adjacent sections from a healthy control stained for phosphorylated α -Syn (MJFR-13, abcam ab168381) are shown with (left) and without (right) proteinase K treatment demonstrating robust detection of insoluble LP. b. Higher power view of Lewy Pathology (solid arrow) in a PD myenteric plexus ganglion. Samples obtained from Banner Health Brain and Body Donation Program. Scale bars are 50 μ m.


Acknowledgments

We are grateful to the Banner Sun Health Research Institute Brain and Body Donation Program of Sun City, Arizona, for the provision of human colon tissue and data. The Brain and Body Donation Program is supported by the National Institute of Neurological Disorders and Stroke (U24 NS072026 National Brain and Tissue Resource for Parkinsons Disease and Related Disorders), the National Institute on Aging (P30 AG19610 Arizona Alzheimers Disease Core Center), the Arizona Department of Health Services (contract 211002, Arizona Alzheimers Research Center), the Arizona Biomedical Research Commission (contracts 4001, 0011, 05-901 and 1001 to the Arizona Parkinsons Disease Consortium) and the Michael J. Fox Foundation for Parkinsons Research.

We wish to acknowledge Yunfeng Li and the Translational Pathology Core Laboratory (TPCL) at UCLA for assistance with tissue processing, IHC staining, and scanning.

Materials

Note

- 
BOND Polymer Refine Detection Kit
Leica Biosystems Catalog #DS9800
contains
- Peroxide Block (30 mL) 3–4% (v/v) Hydrogen peroxide.
 - Polymer (30 mL) Anti-rabbit Poly-HRP-IgG (<25 µg/mL) containing 10% (v/v) animal serum in tris-buffered saline/0.1% ProClin™ 950.
 - DAB Part 1 (2.4 mL) 66 mM 3,3'-Diaminobenzidine tetrahydrochloride hydrate, in a stabilizer solution.
 - DAB Part B (30 mL) ≤0.1% (v/v) Hydrogen Peroxide in a stabilizer solution.
 - Hematoxylin (30 mL) <0.1% Hematoxylin

Equipment

new equipment	NAME
Leica BOND RX	BRAND
3342171	SKU
https://www.leicabiosystems.com/ihc-ish-fish/ihc-ish-instruments/products/leica-bond-rx/	SPECIFICATIONS



Protocol materials


- ✕ Bond™ Wash Solution 10X Concentrate **Leica Biosystems Catalog #AR9590**
- ✕ Bond™ Wash Solution 10X Concentrate **Leica Biosystems Catalog #AR9590**
- ✕ Anti-Alpha-synuclein (phospho S129) antibody [MJF-R13 (8-8)] **Abcam Catalog #AB168381**
- ✕ EnVision System-HRP, Labelled Polymer (Rabbit) **Agilent Technologies Catalog #K4003**
- ✕ Histo-Clear **National Diagnostics Catalog #HS2001GLL**
- ✕ Chemical PermMount Mounting Medium **Fisher Scientific Catalog #SP15-100**
- ✕ BOND Polymer Refine Detection Kit **Leica Biosystems Catalog #DS9800**
- ✕ Dewax Solution **Leica Biosystems Catalog #AR9222**
- ✕ BOND Enzyme Pretreatment Kit **Leica Biosystems Catalog #AR9551**
- ✕ Epitope Retrieval Solution 1 BOND **Leica Biosystems Catalog #AR9961**

Troubleshooting




1 Leica Bond RX routine factory based "Bake and Dewax" protocol.

 Dewax Solution **Leica Biosystems Catalog #AR9222**


2 **Enzyme Antigen Retrieval:** Incubate in Enzyme 1 solution for  00:20:00 at

20m


 Room temperature

 BOND Enzyme Pretreatment Kit **Leica Biosystems Catalog #AR9551**

Note

For comparison, this step can be replaced with Heat Induced Antigen Retrieval using ER1 buffer for  00:20:00

 Epitope Retrieval Solution 1 BOND **Leica Biosystems Catalog #AR9961**

3 Incubate in Peroxide Block for  00:05:00


5m

4 Incubate in Bond Wash buffer three times  00:00:00 each.

 Bond™ Wash Solution 10X Concentrate **Leica Biosystems Catalog #AR9590**

5 **Primary antibody:** Incubate in 1:10,000

1h

 Anti-Alpha-synuclein (phospho S129) antibody [MJF-R13 (8-8)] **Abcam Catalog #AB168381**


for  01:00:00 at  Room temperature

6 Incubate in Bond Wash buffer three times  00:00:00 each.

 Bond™ Wash Solution 10X Concentrate **Leica Biosystems Catalog #AR9590**

7 Incubate with labeled polymer  00:10:00

10m

 EnVision System-HRP, Labelled Polymer (Rabbit) **Agilent Technologies Catalog #K4003**











8 Incubate in Bond Wash buffer three times  00:02:00 each.

2m

9 Incubate with polymer for  00:08:00 .

8m



- 10 Incubate in Bond Wash buffer five times  00:00:00 each.
- 11 Deionized water wash  00:00:00 .
- 12 Incubate with Mixed DAB Refine for  00:00:00 followed by Mixed DAB Refine for  00:10:00 . 10m
- 13 Deionized water wash for three times  00:00:00 each.
- 14 Incubate with Hematoxylin for  00:10:00 . 10m
- 15 Incubate in Bond Wash buffer three times  00:00:00 each.
- 16 Deionized water wash  00:00:00 .
- 17 Slides are dehydrated in series of alcohols, cleared with
 Histo-Clear **National Diagnostics Catalog #HS2001GLL** and mounted with
 Chemical Permunt Mounting Medium **Fisher Scientific Catalog #SP15-100**

Protocol references

Beach TG, Corbillé AG, Letournel F, Kordower JH, Kremer T, Munoz DG, Intorcía A, Hentz J, Adler CH, Sue LI, Walker J, Serrano G, Derkinderen P. Multicenter Assessment of Immunohistochemical Methods for Pathological Alpha-Synuclein in Sigmoid Colon of Autopsied Parkinson's Disease and Control Subjects. *J Parkinsons Dis*. 2016 Oct 19;6(4):761-770. doi: 10.3233/JPD-160888. PMID: 27589538; PMCID: PMC5501392.

Citations

Beach TG, Corbillé AG, Letournel F, Kordower JH, Kremer T, Munoz DG, Intorcía A, Hentz J, Adler CH, Sue LI, Walker J, Serrano G, Derkinderen P. Multicenter Assessment of Immunohistochemical Methods for Pathological Alpha-Synuclein in Sigmoid Colon of Autopsied Parkinson's Disease and Control Subjects.
<https://doi.org/>