



Mar 29, 2024

Autophagosome live cell imaging

DOI

<https://dx.doi.org/10.17504/protocols.io.5qpvnob14o/v1>

Dan Tudorica¹, bishal.penn²

¹Hurley Lab, QB3, UC Berkeley; ²University of Pennsylvania



Dan Tudorica

Hurley Lab, QB3, UC Berkeley

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.5qpvnob14o/v1>

Protocol Citation: Dan Tudorica, bishal.penn 2024. Autophagosome live cell imaging. **protocols.io**
<https://dx.doi.org/10.17504/protocols.io.5qpvnob14o/v1>

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: March 29, 2024



Last Modified: March 29, 2024

Protocol Integer ID: 97564

Keywords: live cell imaging for quantification, live cell imaging, cell

Abstract

For quantification of PINK-PARKIN mitophagy

Cell Culture

- 1 Seed Parkin-expressing HeLa cells for ~80% confluence on day of transfection on tissue-culture coated glass slides, suitable for use on a confocal microscope.
- 2 On day of transfection, transfect pCIG2-LAMP1-mNeon Green, pHaloTag-LC3B, and an autophagy regulating protein tagged with mcherry via lipofectamine 3,000 (ThermoFisher). Briefly, dilute plasmids in Opti-mem reduced serum medium and add P3000 reagent in accordance with manufacturer instructions. Also dilute lipofectamine 3000 reagent in optimem, then mix the dilute plasmid and lipofectamine solutions together. Incubate for ~10 minutes, then add dropwise to adhered cells.
- 3 Allow cells to recover overnight
- 4 On day of experiment swap cells into growth medium containing 10 μ M Oligomycin A and 5 μ M Antimycin A. Return to incubator for 4 H

Imaging and Analysis

- 5 When ready to image, swap medium with Leibovitz's L-15 medium (Gibco 11415064) supplemented with 10% fetal bovine serum and 1% Glutamax along with Antimycin A and Oligomycin A.
- 6 Calibrate laser power and exposure for best image quality. Keep these parameters consistent across treatments or conditions. Collect Z stacks of each target cells.
- 7 To analyze, threshold LC3 and LAMP1 channels, then binarize using Yen thresholding in ImageJ. To quantify area of colocalization, use the 'AND' function from the Image Calculator function. Then calculate area of colocalization using Analyze Particles.