

Jun 19, 2025

Autophagy detection assay

 [Nature Communications](#)

DOI

<https://dx.doi.org/10.17504/protocols.io.36wggqkz3gk5/v1>

Bishal Basak^{1,2}, Erika Holzbaur^{1,2}

¹Department of Physiology, Perelman School of Medicine, University of Pennsylvania, Philadelphia, PA, 19104;

²Aligning Science Across Parkinson's (ASAP) Collaborative Research Network, Chevy Chase, MD, 20815

University of Pennsylv...



Bishal Basak

University of Pennsylvania

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.36wggqkz3gk5/v1>

Protocol Citation: Bishal Basak, Erika Holzbaur 2025. Autophagy detection assay. **protocols.io**
<https://dx.doi.org/10.17504/protocols.io.36wggqkz3gk5/v1>



Manuscript citation:

<https://www.biorxiv.org/content/10.1101/2024.06.12.598718v1.full>

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: May 13, 2025

Last Modified: June 19, 2025

Protocol Integer ID: 218190

Keywords: autophagy detection, autophagosome detection, lysotracker green in neuron, using dojindo dapred dye, dojindo dapred dye, lysotracker green, autophagosome

Funders Acknowledgements:

ASAP

Grant ID: ASAP-000350

Abstract

Autophagosome detection using Dojindo DAPRed dye and Lysotracker Green in neurons

Protocol materials

 DAP Red **Dojindo Catalog #D677**

 Lysotracker Green **Thermo Fisher Scientific Catalog #L7526**



Plate neurons

- 1 Plate 120,000–150,000 neurons in an imaging dish at *days in vitro* (DIV) 0

Note

Protocol for primary neuron isolation and in vitro culture can be obtained from:
[dx.doi.org/10.17504/protocols.io.81wgby723vpk/v1](https://doi.org/10.17504/protocols.io.81wgby723vpk/v1)

Preparing reagents

- 2 On the day of imaging, equilibrate 8 ml of fresh maintenance media (MM) per imaging dish in an incubator

at  37 °C

Note

Maintenance media composition can be obtained from:
[dx.doi.org/10.17504/protocols.io.81wgby723vpk/v1](https://doi.org/10.17504/protocols.io.81wgby723vpk/v1)

- 3 Prepare intermediate dilution of DAP Red in MM: for each imaging dish, make

 500 µL of fresh MM containing  1 µL of DAP Red

 DAP Red **Dojindo Catalog #D677**

- 4 Prepare intermediate dilution of LysoTracker Green in MM: for each imaging dish make

 50 µL of fresh MM containing  0.1 µL of LysoTracker Green

 LysoTracker Green **Thermo Fisher Scientific Catalog #L7526**

- 5 Prepare imaging media (IM) with LysoTracker Green- for each imaging dish make

 50 µL of fresh MM containing  0.1 µL of LysoTracker Green



Note

Imaging media composition for 50 ml:

45% Glucose 660 μ L

B-27 1 mL

Hibernate E Up to 50 mL

Imaging steps

1h

- 6 At DIV 6 or 7, aspirate out conditioned media from the imaging dish, wash the neurons once with fresh MM
- 7 Add  1.5 mL of fresh MM and then add  0.5 mL of MM with the intermediate dilution of DAP Red (prepared in step 3)
- 8 Incubate the neurons in an incubator at  37 °C for  00:30:00 30m
- 9 Post incubation, aspirate out the media, and wash once with fresh MM
- 10 Add  2 mL of fresh MM followed by  50 μ L of MM with the intermediate dilution of LysoTracker Green (prepared in step 4)
- 11 Incubate the neurons in an incubator at  37 °C for  00:30:00 30m
- 12 Post incubation, aspirate out the media, and wash once with fresh IM.
- 13 Add  2 mL of fresh IM followed by  50 μ L of IM with the intermediate dilution of LysoTracker Green (prepared in step 5)
- 14 Image neurons live under a confocal microscope with a heating chamber kept at  37 °C

