


Oct 20, 2019

Bacteria Staining

DOI

dx.doi.org/10.17504/protocols.io.8gshtwe



Guillermo Fernández Rodríguez¹

¹AEGIS - Madrid iGEM 2019

AEGIS - Madrid iGEM 2019



Guillermo Fernández Rodríguez

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.8gshtwe>

Protocol Citation: Guillermo Fernández Rodríguez 2019. Bacteria Staining. **protocols.io**
<https://dx.doi.org/10.17504/protocols.io.8gshtwe>

License: This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited



Protocol status: Working

We use this protocol and it's working

Created: October 19, 2019

Last Modified: October 20, 2019

Protocol Integer ID: 28914

Keywords: Staining, Crystal violet, Coating, Maxisorp-96well, staining protocol, bacteria, replicates per dilution, well plate, dilution

Abstract

A bacteria staining protocol has been automated by OT-2. It allows to check the amount of target we had coated on the 96 well plate.

(We used 5 replicates per dilution)

Materials

MATERIALS

⊗ Crystal violet **Gold Biotechnology Catalog #C-328**

⊗ Nuclease-free water or water filtered using a Milli-Q filtering system **Ambion Catalog #AM9932**

⊗ Sodium bicarbonate **Merck MilliporeSigma (Sigma-Aldrich) Catalog #S6014**

⊗ PBST (PBS 1:1000 Tween-20)

⊗ Centrifuge **Eppendorf Catalog #5415D**

⊗ LB **Research Products International Corp (RPI) Catalog #L24400-2000.0**

⊗ White 96-Well Immuno Plates, Maxisorp, Flat-Bottom, MaxiSorp, 350 μ L **Thermo Fisher Catalog #436110**
















Troubleshooting

Before start

Clean all the working surface with ethanol.



Staining bacteria

- 1 Inoculate a single colony of E.Coli DH5 α from LB agar plate in  10 mL of LB. Use a sterile pipette tip, selecting a single colony from LB agar plate. The liquid culture is incubated overnight  37 °C .
- 2 Spin at  4000 rpm for  00:05:00 . Discard the supernatant, collect pellet and re-suspend in  10 mL of NaHCO₃-Na₂CO₃, 50mM, pH 9,6. Mix by inverting the tube.
- 3 Spin at 4000 rpm for  00:05:00 Discard the supernatant, collect pellet and re-suspend in  8 mL of NaHCO₃-Na₂CO₃, 50mM, pH 9,6. Mix by inverting the tube.
- 4 Read the absorbance (600nm). Dilute the sample with NaHCO₃-Na₂CO₃, 50mM, pH 9,6. and adjust the absorbance to 1.
- 5 Make the following dilutions with -Na₂CO₃, 50mM, pH 9,6:1:5, 1:10, 1:30, 1:50 and 1:100.
- 6 Add  200 μ L the sample into 96 well-plate Nunc MaxiSorp. Incubate overnight at  4 °C .
- 7 Wash 3x  200 μ L PBS Tween 0,1%, pH 7,4. Remove the drops after the last wash.
- 8 Add  150 μ L of crystal violet/well. Incubate for  00:15:00 at  Room temperature
- 9 Wash 4x  250 μ L distilled water. Wait for  01:30:00 for air-drying before counting colonies in the microscope.