Oct 18, 2020

Bacterial Transformation

DOI

dx.doi.org/10.17504/protocols.io.bnifmcbn

Jiaxin Li¹

¹South China University of Technology







DOI: dx.doi.org/10.17504/protocols.io.bnifmcbn

Protocol Citation: Jiaxin Li 2020. Bacterial Transformation. protocols.io https://dx.doi.org/10.17504/protocols.io.bnifmcbn

License: This is an open access protocol distributed under the terms of the <u>Creative Commons Attribution License</u>, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working We use this protocol and it's working

Created: October 16, 2020

Last Modified: October 18, 2020

Protocol Integer ID: 43303

- 1 Thaw a tube of DH5a Competent E. coli cells on ice for 5-10 minutes
- 2 Add 1µl of plasmid DNA to per 100µl cell mixture. Flick the tube 4-5 times to mix cells and DNA but do not vortex.
- 3 Place the mixture on ice for 30 minutes. Do not mix.
- 4 Heat shock at exactly 42°C for exactly 90 seconds. Immediately place on ice for 2-3 minutes. Do not mix.
- 5 Pipette 800 µl of 42°C LB liquid medium into the mixture.
- 6 Place at 37°C, 180-220 rpm for 45-60 minutes.
- 7 Centrifuge the bacterium at 3000r for 3 minutes.
- 8 Discard 300 µl supernatant and resuspend the sediment.
- 9 Spread 50-100 μ l of each dilution onto a selection plate and incubate for 12 hours at 37°C.