

Oct 18, 2020

## Bacterial Transformation

DOI

[dx.doi.org/10.17504/protocols.io.bnifmcbn](https://dx.doi.org/10.17504/protocols.io.bnifmcbn)

Jiaxin Li<sup>1</sup>

<sup>1</sup>South China University of Technology



Jiaxin Li

### Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.bnifmcbn>

**Protocol Citation:** Jiaxin Li 2020. Bacterial Transformation. **protocols.io** <https://dx.doi.org/10.17504/protocols.io.bnifmcbn>

**License:** This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

**Protocol status:** Working

**We use this protocol and it's working**

**Created:** October 16, 2020

**Last Modified:** October 18, 2020

**Protocol Integer ID:** 43303

**Keywords:** transformation

## Troubleshooting



- 1 Thaw a tube of DH5a Competent E. coli cells on ice for 5-10 minutes
- 2 Add 1µl of plasmid DNA to per 100µl cell mixture. Flick the tube 4-5 times to mix cells and DNA but do not vortex.
- 3 Place the mixture on ice for 30 minutes. Do not mix.
- 4 Heat shock at exactly 42°C for exactly 90 seconds. Immediately place on ice for 2-3 minutes. Do not mix.
- 5 Pipette 800 µl of 42°C LB liquid medium into the mixture.
- 6 Place at 37°C, 180-220 rpm for 45-60 minutes.
- 7 Centrifuge the bacterium at 3000r for 3 minutes.
- 8 Discard 300 µl supernatant and resuspend the sediment.
- 9 Spread 50-100 µl of each dilution onto a selection plate and incubate for 12 hours at 37°C.