


Jun 24, 2017

Bacteriophage Cryopreservation

DOI

dx.doi.org/10.17504/protocols.io.iyyccfw



Dr. Steven Wilhelm

The Aquatic Microbial E...



Steven W Wilhelm

The University of Tennessee, Knoxville

OPEN  ACCESS



DOI: dx.doi.org/10.17504/protocols.io.iyyccfw

Protocol Citation: Dr. Steven Wilhelm 2017. Bacteriophage Cryopreservation. **protocols.io**

<https://dx.doi.org/10.17504/protocols.io.iyyccfw>

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

Created: June 19, 2017

Last Modified: November 29, 2017

Protocol Integer ID: 6456

Abstract

Please contact Dr. Steven Wilhelm (wilhelm@utk.edu) for additional information regarding this protocol.

Materials

STEP MATERIALS

 Glycerol growcells.com Catalog #MRGE-4002

 Glycerol growcells.com Catalog #MRGE-4002

Protocol materials





 Glycerol growcells.com Catalog #MRGE-4002

 Glycerol growcells.com Catalog #MRGE-4002

 Glycerol growcells.com Catalog #MRGE-4002



Preparation

- 1 Vortex a freshly lysed culture for 10 sec
 00:00:10
- 2 Add 1 mL of phage into each cryogenic vial
 1 mL
- 3 Add 1 mL of 20% glycerol to each vial with phage
 1 mL
 Glycerol growcells.com Catalog #MRGE-4002
- 4 Invert tubes several times
- 5 Place immediately in a -80°C freezer

Transferring frozen virus particles

- 6 Remove cryogenic vials one at a time from the freezer
- 7 Using a sterile pipette tip, scrape a small amount of frozen virus particles into a 20 mL mid-log culture
- 8 Place cryogenic vials immediately back into the freezer
- 9 Allow to incubate and lyse culture under constant light conditions appropriate for your specific culture