

Jan 19, 2024

Version 1

BAF_Protocol_001 In-gel Digestion V.1

DOI

dx.doi.org/10.17504/protocols.io.14egn7r5mv5d/v1

Nicholas Sherman¹

¹University of Virginia

Nicholas Sherman: Biomolecular Analysis Facility Core



Nicholas Sherman

University of Virginia Biomolecular Analysis Facility Core

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.14egn7r5mv5d/v1>

Protocol Citation: Nicholas Sherman 2024. BAF_Protocol_001 In-gel Digestion. **protocols.io**

<https://dx.doi.org/10.17504/protocols.io.14egn7r5mv5d/v1>

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working



Created: June 21, 2022

Last Modified: May 05, 2025

Protocol Integer ID: 65051

Keywords: proteomics, mass spectrometry, digestion, gel digestion this protocol, gel digestion, gel digestion of protein, gel method, clean peptide digest, baf-protocol-001 in, gel, gel band, peptides for mass spectrometry analysis, peptide, protein, digest, mass spectrometry analysis, protocol

Abstract

This protocol is for in-gel digestion of proteins, including large mixtures, to produce peptides for mass spectrometry analysis. The gel method can be useful for getting rid of detergents or small molecules that might interfere with minimal loss. The input is a gel band up to 1cm x 1cm. The output is a relatively clean peptide digest that is ready for quick cleanup by C18 tip.

Guidelines

1. Wash hands before starting. Refrain from wearing heavy lotions or perfumes.
2. Do not wear clothing that contains wool or a lot of loose fibers.
3. Use nitrile gloves (not latex).
4. Perform all procedures under PCR clean hood cabinet.
5. All sample preparation must use polypropylene tubes cleaned by addition/removal of EtOH (pure ethanol), then H₂O (water), then EtOH. Dry in clean hood for several hours upside down. This strips off polymer coating.
7. All reagents must be prepared in small glass bottles (polypropylene lined caps - no glue foil).
8. Solutions should be prepared fresh.



Materials

Pre-cleaned microtubes 1.5 mL - USA scientific, SEAL-RITE® 1.5 mL MICROCENTRIFUGE TUBES color: natural.
Pre-cleaned microtubes 0.5 mL - USA Scientific, SEAL-RITE® 0.5 mL MICROCENTRIFUGE TUBES color: natural.
Screw cap microtubes 1.5 mL - Fisher Brand, Conical Screw Cap Tube, color: natural, part number: 02-681-339.
20mL liquid scintillation vials PE cone lined cap - Wheaton DWK986756
Pipette tips - Fisher Brand, color: yellow, part number: 02-681-151.
Micropipettes
ABC - Fluka analytical, Ammonium Bicarbonate, Sigma Aldrich, part number: 09830.
DTT - Fisher BioReagents, Dithiothreitol, part number: BP 172-5.
Iodoacetamide - Sigma Aldrich, Iodoacetamide, part number: 1149.
Promega Trypsin: Promega, Sequencing grade modified, frozen, part number: V511C.
FA - Fisher Chemical, Formic Acid, Optima™ LC/MS Grade, part number: A117-50.
MeOH - Fisher Chemical, Methanol, Optima™ LC/MS Grade, part number: A456-4.
ACN - Fisher Chemical, Acetonitrile, Optima™ LC/MS Grade, part number: A955-4.
TFA - Pierce™, Trifluoroacetic Acid, part number: 28903.
Water - Fisher Chemical, Optima™ LC/MS Grade, part number: W6-4.
Syringe - Unimetrics PKS250, 250µL Peek Laboratory Syringe
2 to 20 µL Micropipette - Gilson™ F144056MT
10 to 100 µL Micropipette - Gilson™ F144057MT
20 to 200 µL Micropipette - Gilson™ F144058MT
100 to 1000 µL Micropipette - Gilson™ F144059MT

Troubleshooting

Before start

REAGENTS: (All reagents to be prepared fresh for each digestion)

1. Destain solution: 10 mL MeOH (methanol), 10 mL H₂O (water), 500 µL Acetic Acid (glacial)
2. 100 mM ABC (ammonium bicarbonate): 0.158 g in 20 mL distilled H₂O
3. 50 mM ABC: 0.079g in 20 mL distilled water
4. ACN (Acetonitrile)
5. 10mM DTT: 0.0015 g in 1 mL of 100 mM ABC (Use screw caps microtubes. DO NOT mix until directly before you are ready to use)
6. 50 mM Iodoacetamide: 0.01 g in 1 mL of 100 mM ABC (Use screw caps microtubes. DO NOT mix until directly before you are ready to use)
7. Trypsin solution: Keep on ice. Promega sequencing grade trypsin (cat. # V5113 (porcine)) dilute when ready to use in 50 mM ABC. It comes frozen in 40µL of 50 mM Acetic Acid
8. Extraction solution 1: 5% formic acid in H₂O (950/50)
9. Extraction solution 2: 5% formic acid in 50% MeCN/H₂O (500/450/50)

**DAY 1:**

19h 30m

- 1 Cut gel band into small pieces (~2 mm). 10m
- 2 Add the pieces in a pre-cleaned 1.5 mL polypropylene tube. 5m
- 3 Destain the pieces in 0.5 mL destain solution for ~2hrs. 2h
- 4 Remove destain (using 1 gel load tip for all tubes in batch) and replace with 0.5 mL destain solution for 30min. It's okay if some stain remains. It will be removed during further wash steps. 30m
- 5 Remove the destain solution and dehydrate gel slices in 200 μ L ACN (Acetonitrile) for 5 min. 5m
- 6 Remove ACN and repeat. 5m
- 7 Reduce the gel pieces in 30 μ L of 10 mM DTT (dithiothreitol) for 0.5 h at room temperature. 30m
- 8 Remove the DTT solution. 1m
- 9 Alkylate in 30 μ L 50 mM IA (iodoacetamide) at room temperature, in dark for 0.5 h. (get ice for 50 mM ABC). 30m
- 10 Remove the IA solution. 1m
- 11 Wash gel pieces with 100 μ L 100 mM ABC (ammonium bicarbonate) for 10 min and remove. 10m
- 12 Dehydrate gel pieces in 200 μ L ACN for approximately 5 min. Remove ACN. 5m
- 13 Rehydrate gel pieces in 200 μ L 100 mM ABC for 10 min.



- | | | |
|----|--|-----|
| | | 10m |
| 14 | Remove ABC supernatant. | 1m |
| 15 | Dehydrate gel slices in 200 μ L ACN approximately 5 min. Remove ACN and repeat. | 5m |
| 16 | Prepare trypsin.
For Gel Pieces: Add 960 μ L of cold 50 mM ABC to 20 μ g of Promega trypsin (Part number V5113) (it comes diluted in 40 μ L of 50 mM acetic acid) on ice (20 μ g/mL). | 2m |
| 17 | Add 30-50 μ L of the trypsin solution to cover the gel pieces and rehydrate on ice for 30 min. | 30m |
| 18 | Microfuge and remove any excess trypsin solution and add 5-20 μ L of 50 mM ABC. React overnight at 37°C. | 16h |

Day 2:

3h 40m

- | | | |
|----|---|-----|
| 19 | Prepare Extraction Solutions:
a. Extraction Solution 1: H ₂ O/Formic Acid (950/50)
b. Extraction Solution 2: ACN/H ₂ O/Formic Acid (500/450/50) | 10m |
| 20 | DO NOT remove excess ABC supernatant from yesterday!!! | |
| 21 | Add 10 μ L Extraction Solution 1 to each tube (or more if needed to cover the gel pieces) Let rest for 10 minutes. Take off supernatant to a clean 0.5 mL tube. | 10m |
| 22 | Add 10 μ L Extraction Solution 2. Incubate for 10 min. Take off supernatant to the same 0.5 mL tube. | 10m |
| 23 | Repeat the 10 μ L Extraction Solution procedure. Take off supernatant to the same 0.5 mL tube. | 10m |
| 24 | Evaporate the sample via speed vac and reconstitute to 10-100 μ L total volume (depending on the amount of total protein) with 0.1% Formic acid. | 2h |
| 25 | Perform desalting using C18 tips (BAF_Protocol_003):
– 10 μ L tips for a small amount of proteins | 1h |



– 100 μ L for a higher amount of proteins.