

Apr 05, 2024

## BaseScope In Situ Hybridization

DOI

[dx.doi.org/10.17504/protocols.io.5qpvo364zv4o/v1](https://dx.doi.org/10.17504/protocols.io.5qpvo364zv4o/v1)

madalynn.erb Erb<sup>1</sup>

<sup>1</sup>Van Andel Research Institute

Madalynn's Workspace



madalynn.erb Erb

Van Andel Research Institute

### Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.5qpvo364zv4o/v1>

**Protocol Citation:** madalynn.erb Erb 2024. BaseScope In Situ Hybridization. **protocols.io**  
<https://dx.doi.org/10.17504/protocols.io.5qpvo364zv4o/v1>

**License:** This is an open access protocol distributed under the terms of the **[Creative Commons Attribution License](#)**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited



**Protocol status:** Working

**We use this protocol and it's working**

**Created:** January 29, 2024

**Last Modified:** May 31, 2024

**Protocol Integer ID:** 94340

**Keywords:** ASAPCRN, basescope in situ hybridization basescope, situ hybridization basescope, basescope, situ hybridization, mouse brain section, mouse brain

## Disclaimer

The [protocols.io](https://www.protocols.io) team notes that research involving animals and humans must be conducted according to internationally-accepted standards and should always have prior approval from an Institutional Ethics Committee or Board.

## Abstract

BaseScope *in situ* hybridization on mouse brain sections

## Troubleshooting



## Day 1

- 1 Use 35 $\mu$ M floating brain sections
- 2 Wash sections 3 times (5 minutes) in PBS to remove cryoprotectant solution
- 3 Mount tissue onto Superfrost plus slides
- 4 Dry slides at room temp and then dip in mQ H<sub>2</sub>O to remove salt
- 5 Dry slides at 60°C for 2 hours
- 6 Let tissue dry at room temp overnight – make sure tissue is completely stuck to slides

## Day 2

- 7 Dry slides at 60°C for 30 minutes
- 8 Dehydrate slides in **50%, 70%, 100% EtOH**, for 5 min each at room temp

### Note

Use sterile distilled H<sub>2</sub>O for the rest of the protocol

- 8.1 Let slides air dry for 5 min at room temp.
- 9 Make 200mL Target Retrieval Reagent (1:10 180mL H<sub>2</sub>O + 20mL 10X Target Retrieval Reagent)
- 10 On the bench – add 5-8 drops of Hydrogen Peroxide to tissue sections



- 10.1 Incubate at room temp for 10 min
- 11 Rinse slides 2 times in H<sub>2</sub>O
  - 11.1 Submerge in dish of H<sub>2</sub>O and move slides up and down 3-5 times for each wash
- 12 In the steamer heat one dish of H<sub>2</sub>O and one dish of Target Retrieval Reagent
  - 12.1 Check temp of Target retrieval reagent (must be at least 95°C)
  - 12.2 After preheating – place slides into H<sub>2</sub>O for 10 sec to acclimate.
  - 12.3 Move slides to Target Retrieval Reagent and steam for 5 min
- 13 Remove the water container and the container with slides from the steamer and place on the bench.
  - 13.1 Let slides cool on the bench for 20 -30 minutes (wait until temperature of water container has dropped to 30°C)
- 14 Wash slides in H<sub>2</sub>O for 15 sec at room temp
- 15 Wash slides in 100% EtOH for 3 min at room temp
- 16 Dry slides at 60°C for ≥5 min
- 17 Use an Immedge hydrophobic barrier pen to draw a barrier around your tissue section(s)



- 17.1 Let the barrier dry  $\geq 1$  min or overnight at room temp (you can take a break at this point)
- 18 Turn on RNA Scope oven – set to 40°C
- 18.1 Wet humidifying paper (filter paper) completely with H<sub>2</sub>O and place in the Humidity Control Tray
- 18.2 Warm the covered tray in the oven for 30 minutes before use
- 19 Place slides in slide rack and add ~5 drops of Protease IV
- 19.1 Incubate in pre-warmed RNAscope oven in the humidity tray at 40°C for 30 min
- 20 Wash slides in H<sub>2</sub>O at room temp
- 21 Prepare 1L of wash buffer (20mL of 50X Wash buffer + 980mL H<sub>2</sub>O)
- 21.1 *Heat 50X wash buffer to 40°C for 10-20 min before making 1X wash buffer*
- 22 Equilibrate each BaseScope reagent and probe to room temp for 30 min before use
- Important:** *Do not let sections dry out between incubation steps*  
  
*All incubation steps are in the slide rack – remove slides from slide rack to wash between incubations*  
  
*For the incubations in the HybEZ oven – cover the slide rack with the lid and make sure to lock the oven so that the water in the humidity chamber does not evaporate.*
- 23 Add ~4 drops of BaseScope probe to completely cover tissue



- 23.1 Incubate in the oven at 40°C for **2 hours**
- 23.2 Wash 2 times in Wash Buffer at room temp - 2 min each wash
- 24 Add ~4 drops of **AMP1** to completely cover tissue
- 24.1 Incubate in the oven at 40°C for **30 min**
- 24.2 Wash 2 times in Wash Buffer at room temp - 2 min each wash
- 25 Add ~4 drops of **AMP2** to completely cover tissue
- 25.1 Incubate in the oven at 40°C for **30 min**
- 25.2 Wash 2 times in Wash Buffer at room temp - 2 min each wash
- 26 Add ~4 drops of **AMP3** to completely cover tissue
- 26.1 Incubate in the oven at 40°C for **15 min**
- 26.2 Wash 2 times in Wash Buffer at room temp - 2 min each wash
- 27 Add ~4 drops of **AMP4** to completely cover tissue
- 27.1 Incubate in the oven at 40°C for **30 min**
- 27.2 Wash 2 times in Wash Buffer at room temp - 2 min each wash



28 Add ~4 drops of **AMP5** to completely cover tissue

28.1 Incubate in the oven at 40°C for **30 min**

28.2 Wash 2 times in Wash Buffer at room temp - 2 min each wash

29 Add ~4 drops of **AMP6** to completely cover tissue

29.1 Incubate in the oven at 40°C for **15 min**

29.2 Wash 2 times in Wash Buffer at room temp - 2 min each wash

29.3 *We are done with the oven – but keep using the slide rack for incubations at room temp*

30 Add ~4 drops of **AMP7** to completely cover tissue

30.1 Incubate at **room temp** for **15 min**

30.2 Wash 2 times in Wash Buffer at room temp - 2 min each wash

30.3 *Staining intensity can be modified by adjusting the AMP7 incubation time*

31 Add ~4 drops of **AMP8** to completely cover tissue



- 31.1 Incubate on the bench at **room temp** for **15 min**
- 31.2 Wash 2 times in Wash Buffer at room temp - 2 min each wash
- 32 Prepare BaseScope Fast RED working solution (1:60 ratio)
- 32.1 Spin down BaseScope Fast RED-B before using
- 32.2 2uL RED-B + 120uL RED-A : mix well
- 32.3 Use the Fast Red solution **within 5 minutes** – do not expose to direct sunlight or UV light
- 33 Add ~120μL Fast Red solution to each slide / section
- 33.1 Cover sections on tray to protect from light
- 33.2 Incubate at room temp for **10 min**
- 33.3 Rinse 2 times in H<sub>2</sub>O at room temp
- 34 Dry slides at 60°C for ≥ 15 min (until slides are completely dry)
- 34.1 *The red substrate is alcohol sensitive. DO NOT dehydrate the slides in alcohol, and make sure your reagents are not contaminated with alcohol.*
- 35 Place 1-2 drops of VectaMount or Ecomount on the slide and coverslip the tissue sections.
- 36 Dry slides at room temp for ≥ 5 min



