

Feb 17, 2025

Beam Crossing

 Forked from [Beam test](#)

DOI

<https://dx.doi.org/10.17504/protocols.io.n92ldr2jxg5b/v1>

Marina Lorente Picón¹, Núria Peñuelas¹, Ariadna Laguna¹, Miquel Vila¹, Chelsie Steele², Yujie Fan²

¹Vall d'Hebron Research Institute; ²Caltech



Chelsie Steele

Caltech

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.n92ldr2jxg5b/v1>

Protocol Citation: Marina Lorente Picón, Núria Peñuelas, Ariadna Laguna, Miquel Vila, Chelsie Steele, Yujie Fan 2025. Beam Crossing. [protocols.io](#) <https://dx.doi.org/10.17504/protocols.io.n92ldr2jxg5b/v1>

License: This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working



Created: January 08, 2025

Last Modified: February 17, 2025

Protocol Integer ID: 117859

Keywords: beam crossing beam test, mice

Funders Acknowledgements:

Aligning Science Across Parkinson's

Grant ID: ASAP-020495

Abstract

Beam test for mice

Beam Crossing

- 1 Mice are placed at the beginning of an elevated horizontal plexiglass beam bar covered on scotch tape to avoid limb slip (1m length – 2.5cm width- 0.5m elevation). At the beginning of the beam we place a little lamp and at the end of the beam we place a dark box to motivate the animal to go towards the dark side.
- 2 Mice are allowed to explore the beam bar and receive four assisted trials to cross it a week before the test. This is called the "training period"
- 3 The week after training, the beam crossing test starts and animals can traverse the entire length of the beam unassisted for three times. Each of the three trials are filmed for later analysis in BORIS software.
- 4 Clean the beam bar between animals.

Video Analysis

- 5 For each video, the time it takes the animal to cross the beam is recorded, then the average time of the three trials is taken in BORIS Software (DOI: [10.1111/2041-210X.12584](https://doi.org/10.1111/2041-210X.12584)) for analysis.
- 6 The limb errors per video are also counted. An error is counted when, during a forward movement, a limb (forelimb or hindlimb) slips through the beam. The average number of foot slips is then taken.