


Jul 11, 2019

BG11 medium (working group Wilde)

 [Journal of Bacteriology](#)

DOI

dx.doi.org/10.17504/protocols.io.wj5fcq6

Christin Köbler¹, Annegret Wilde¹

¹Institute of Biology III, Faculty of Biology, University of Freiburg, 79104 Freiburg, Germany

CyanoWorld



Christin Köbler

OPEN  ACCESS



DOI: dx.doi.org/10.17504/protocols.io.wj5fcq6

External link: <https://doi.org/10.1128/JB.00478-19>

Protocol Citation: Christin Köbler, Annegret Wilde 2019. BG11 medium (working group Wilde). **protocols.io**
<https://dx.doi.org/10.17504/protocols.io.wj5fcq6>

Manuscript citation:

Wiegand A, Köbler C, Oyama K, Dörrich AK, Azai C, Terauchi K, Wilde A, Axmann IM, Array. Journal of Bacteriology 202(4). doi: [10.1128/JB.00478-19](https://doi.org/10.1128/JB.00478-19)

License: This is an open access protocol distributed under the terms of the **[Creative Commons Attribution License](#)**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: December 21, 2018

Last Modified: July 11, 2019

Protocol Integer ID: 18781

Abstract

100x BG11 (500 ml)

74.79 g NaNO_3

3.75 g $\text{MgSO}_4 \times 7\text{H}_2\text{O}$

1.8 g $\text{CaCl}_2 \times 2\text{H}_2\text{O}$

0.30 g Citric Acid

0.28 ml 0,5 M Na_2EDTA (pH 8)

Add *Aqua bidest* (**ddH₂O**) to a final volume of 500 ml. Aliquot (5×100 ml) and autoclave. Store at 4°C.

Trace-Metal-Mix (500 ml)

1.43 g H_3BO_3

900 mg $\text{MnCl}_2 \times 4\text{H}_2\text{O}$

110 mg $\text{ZnSO}_4 \times 7\text{H}_2\text{O}$

195 mg $\text{Na}_2\text{MoO}_4 \times 2\text{H}_2\text{O}$

39.5 mg $\text{CuSO}_4 \times 5\text{H}_2\text{O}$

24.7 mg $\text{Co}(\text{NO}_3)_2 \times 6\text{H}_2\text{O}$

Put some **ddH₂O** into the bottle; dissolve the salts separately and fill up to 500 ml. Sterile filter and aliquot the solution (10 × 50 ml). Store at 4°C.

2x BG11 (1 L)

20 ml 100x BG11

2 ml $\text{K}_2\text{HPO}_4 \times 3\text{H}_2\text{O}$ (30 mg/ml)

2 ml Na_2CO_3 (20 mg/ml)

2 ml Trace-Metal-Mix

Add 900 ml deionized water, afterwards add the solutions, fill up to 1 L with deionized water and autoclave. Store at RT.

Before use add 2 ml of Fe-Ammonium-Citrate (6 mg/ml).

1xBG11 (1 L)

10 ml 100xBG11

10 ml TES pH 8



1 ml K_2HPO_4 (30 mg/ml)

1 ml Na_2CO_3 (20 mg/ml)

1 ml Trace-Metal-Mix

Add 900 ml deionized water, afterwards add the solutions, fill up to 1 L with deionized water and autoclave. Store at RT.

Before use add 1 ml of Fe-Ammonium-Citrate (6 mg/ml).

