

Jan 24, 2025

Bi-parental conjugation protocol for non-model bacteria using E.coli WM3064 DAP auxotroph donor

DOI

<https://dx.doi.org/10.17504/protocols.io.kqdg3qbzqv25/v1>

Atharva Karde^{1,2}, Hans Christopher Bernstein^{1,2,3}

¹The Norwegian College of Fishery Science, Faculty of Biosciences, Fisheries and Economics, UiT - The Arctic University of Norway, Tromsø, Norway;

²Microalgae & Microbiomes Research Group (M2RG), UiT - The Arctic University of Norway, Tromsø, Norway;

³The Arctic Centre for Sustainable Energy (ARC), UiT - The Arctic University of Norway, Tromsø, Norway



Atharva Karde

UiT The Arctic University of Norway

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.kqdg3qbzqv25/v1>

Protocol Citation: Atharva Karde, Hans Christopher Bernstein 2025. Bi-parental conjugation protocol for non-model bacteria using E.coli WM3064 DAP auxotroph donor. **protocols.io** <https://dx.doi.org/10.17504/protocols.io.kqdg3qbzqv25/v1>



License: This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We have used this protocol successfully with Arctic marine diatom-associated bacterial recipients.

Created: January 16, 2025

Last Modified: January 24, 2025

Protocol Integer ID: 118568

Keywords: conjugation, bi-parental mating, bi-parental, plasmid, transfer, DAP auxotroph, non-model bacteria, dap auxotroph donor strain, bacteria, plasmid

Funders Acknowledgements:

Kometen – Nord-Norges forskningsfond


Abstract

This is a protocol for using bi-parental conjugation to transfer a plasmid into non-model recipient bacteria from a DAP auxotroph donor strain.

Guidelines


1. **Recipients : donor ratio.** For fast-growing commercial strains, the recipient : donor ratio can be 1 : 1, but non-model/slow growing strains will be rapidly out-competed by the donor. In such cases, a larger recipient : donor ratio will give significantly better results. 10 : 1 or even 100 : 1 are some variations of mixtures worth trying.
2. **Conjugation time.** This may need to be optimized based on the recipient strain used. Some may require 24+ hours of incubation. Tweak the mixture ratio and conjugation time such that the donor strain does not out-compete the recipient cells entirely. The incubation temperature can also be regulated to control this.


Materials

1. **Donor strain** carrying the desired plasmid to be transferred. Your plasmid must contain an oriT sequence in order to be transferred. The donor strain (if using a different one than WM3064) must carry rp4-based transfer machinery and must be  auxotrophic.
2. **Recipient strain**
3. **Growth media** for individual strains. Donor growth media must contain **0.3mM DAP**.
4. **Non-selective plates** containing **DAP** (0.3 mM) for conjugation. These plates should be made from media that **both strains** can grow on, and contain **no antibiotic**.
5. **Selective plates** containing the same antibiotic as the marker on the plasmid you want to transfer, but **no DAP**. To plate dilutions, you need 3-4 plates per conjugation mixture.
6. 
7. 
8. Pipettes, microcentrifuge tubes, and incubator that supports growth of both donor and recipients.
9. Antibiotics for plasmid selection

	A	B
	Antibiotic	Working Concentration
	Ampicillin	100 µg/ml
	Carbenicillin	100 µg/ml
	Chloramphenicol	33 µg/ml
	Kanamycin	50 µg/ml
	Streptomycin	25 µg/ml
	Tetracycline	15 µg/ml

Protocol materials







Before start

1. Verify that the donor strain is carrying the desired plasmid to be transferred
2. Your plasmid **MUST CONTAIN an oriT sequence** in order to be transferred.
3. The donor strain (if using a different one than WM3064) must carry rp4-based transfer machinery and must be  2,6-Diaminopimelic acid **Thermo Scientific** auxotrophic in order to be eliminated during selection.
4. Ensure that appropriate growth medium is used to grow individual strains and that the donor strain growth medium contains  0.3 millimolar (mM)  2,6-Diaminopimelic acid **Thermo Scientific** .
5. Non-selective plates must contain DAP. These plates should be made from media that **both** strains can grow on, and contain no antibiotic. You can set up ~4 conjugation mixtures per plate.
6. Selective plates must contain NO DAP and the appropriate antibiotic depending on the selection marker used in the transferred payload. You may also include additional counter selectable markers to minimize E. coli background.

This protocol was developed as part of a master's thesis affiliated with the Microalgae & Microbiomes Research Group at NFH, UiT, Norway.

M2RG Website : <https://uit.no/research/micro>



Overnight culture preparation

- 1 Inoculate a loopful of donor E.coli WM3064 and recipient strains in respective culture tubes containing  5 mL of growth medium. Ensure that the donor strain is grown in presence of  0.3 millimolar (mM)  2,6-Diaminopimelic acid **Thermo Scientific** and appropriate antibiotic according to the selection marker present on the carried plasmid.
- 2 Incubate the cells  Overnight or until they have grown sufficiently. The recipient cells may take longer to grow.

Conjugation mixture preparation

5m

- 3 Transfer  1 mL of donor and recipient overnight cultures to microcentrifuge tubes.
- 4 Gently spin down the cells at  3000 rpm, 00:05:00 and wash with  500 μ L  1X PBS (Phosphate-buffered saline) .
- 5 Repeat  [go to step #4](#) .
- 6 Resuspend the cells in  500 μ L  1X PBS (Phosphate-buffered saline) .

Note

Washing removes residual antibiotic from the donor cell culture. Recipient cells can also be scraped off of plates and washed in PBS.

- 7 Prepare  100 μ L conjugation mixture by combining the desired volume ratio of donor and recipient cells in a microcentrifuge tube. Mix by pipetting twice.



Note

For fast-growing commercial strains, the recipient : donor ratio can be 1 : 1, but non-model/slow growing strains will be rapidly out-competed by the donor. In such cases, a larger recipient : donor ratio will give significantly better results. 10 : 1 or even 100 : 1 are some variations of mixtures worth trying.

Non-selective plating

- 8 Spot-plate  100 μ L of the mixture onto a **non-selective plate containing DAP and no antibiotic**. DO NOT SPREAD.

Note

Up to 4-5 conjugation mixtures can be set up on the same plate by carefully plating into quadrants on the plate.

- 9 Let the spots dry and incubate the plates  Overnight .

Note

This step may need to be optimized based on the recipient strain used. Some may require 24+ hours of incubation. Tweak the mixture ratio and conjugation time such that the donor strain does not out-compete the recipient cells entirely. The incubation temperature can also be regulated to control this.

Washing

5m

- 10 Scrape each conjugation mixture from the non-selective plates and resuspend in

 1 mL  1X PBS (Phosphate-buffered saline)

- 11 Vortex and gently spin down  3000 rpm, 00:05:00 .

5m

- 12 Resuspend in  1 mL  1X PBS (Phosphate-buffered saline)



13 Repeat  [go to step #11](#)

14 Resuspend the mixture in  100 μ L  1X PBS (Phosphate-buffered saline)

Selective plating

15 Spread-plate the entire  100 μ L or desired dilutions of washed mixtures on **selective plates containing antibiotic and no DAP.**

16 Incubate the plates according to the recipient strain requirements and pick single colonies to restreak or grow in selective media.

17 Confirm the transfer of the payload via PCR.

Protocol references

<https://barricklab.org/twiki/bin/view/Lab/ProtocolsConjugation>

Acknowledgements

This protocol was developed as part of a master's thesis affiliated with the Microalgae & Microbiomes Research Group at NFH, UiT, Norway.

M2RG Website : <https://uit.no/research/micro>

The donor strain used in this project was generously provided by the Lowe-Power Lab at UC Davis (<https://lowepowerlab.ucdavis.edu/>).