



Sep 27, 2019

Version 2

Biochemical Measures of Neuropathy - Hydrogen Peroxide Assay V.2

DOI

dx.doi.org/10.17504/protocols.io.7rghm3w



Eva Feldman¹

¹University of Michigan - Ann Arbor

Diabetic Complications Consortium

Tech. support email: rmcindoe@augusta.edu



Lili Liang

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.7rghm3w>

External link: <https://www.diacomp.org/shared/document.aspx?id=54&docType=Protocol>

Protocol Citation: Eva Feldman 2019. Biochemical Measures of Neuropathy - Hydrogen Peroxide Assay. **protocols.io**
<https://dx.doi.org/10.17504/protocols.io.7rghm3w>

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: September 27, 2019

Last Modified: September 27, 2019

Protocol Integer ID: 28168

Keywords: Biochemical Measures of Neuropathy, diabetic neuropathy, Hydrogen Peroxide , biochemical measures of neuropathy, neuropathy, oxidative stress, increased levels of glucose overload mitochondria, hydrogen peroxide assay summary, antioxidant capacity, hyperglycemia, glucose overload mitochondria, axons of the peripheral nervous system, diabetic complication, biochemical measure, metabolic change, excess glucose, peripheral nervous system, neuron, flow of excess glucose, cellular pathway, reactive oxygen species

Abstract

Summary:

Oxidative stress is highly correlated with the metabolic changes caused by hyperglycemia. Increased levels of glucose overload mitochondria and result in the production of reactive oxygen species (ROS). In addition, the flow of excess glucose through cellular pathways decreases the cell's normal ability to detoxify ROS. As a result, the neurons and axons of the peripheral nervous system contain increased levels of ROS and decreased antioxidant capacity. The following assays are used to measure these changes in rodent models of diabetic neuropathy.


Diabetic Complication:



Neuropathy

Materials

MATERIALS

 Amplex Red Hydrogen Peroxide Assay Kit **Catalog # A-12212**

 HPLC-grade water

Reagent Preparation:

Amplex Red reagent: Prepare a ~20 mM stock solution. Bring DMSO and 1 vial of Amplex Red reagent to room temp. Just prior to use dissolve the Amplex Red reagent in 200 μ L DMSO. Store stock solution at -20°C , protected from light.

Reaction Buffer (5X) (0.25M sodium phosphate, pH 7.4): Dilute 5 mL of Reaction buffer in 20 mL of deionized water.

HRP (Horseradish peroxidase): Dissolve 1 vial of HRP in 1 mL of 1X Reaction Buffer (200 U/mL). After use divide remaining stock into small aliquots and store frozen at -20°C .

20mM H_2O_2 : Dilute (check bottle for %) 17.9 μ L H_2O_2 (3.8%) in 982.1 μ L dH_2O . (Check label for exact concentration) (23 μ L 3% H_2O_2 into 977 μ L dH_2O) Use promptly.

Resorufin, sodium salt: add 1mL dH_2O to a vial of resorufin solid. (2 mM stock). Store at -20°C , protected from light.

Troubleshooting

Sample Preparation:

- 1 Prepare stock solutions above for Amplex Red Reagent, Reaction Buffer, HRP, 20 mM H₂O₂ and Resorufin, sodium salt.
- 2 Prepare H₂O₂ standard curve. Label tubes 1-5.
Add 999 µL 1X Reaction Buffer to tube 2.
Add 100 µL 1X Reaction Buffer to tubes 2 – 5.
Add 1 µL of 20 mM H₂O₂ working dilution prepared above to tube # 1 to produce concentrations of 20 µM.
Take 100 µL from tube 1 and add to tube 2.
Take 100 µL from tube 2 and add to tube 3.
Take 100 µL from tube 3 and add to tube 4.
Tube 5 is dH₂O only.
(Final concentration will be twofold lower, 0 to 10 µM.)
- 3 Dilute samples in 1X Reaction Buffer.
- 4 Pipette 100 µL of diluted standards, controls (if any) and samples into wells. (For DRG we used 25 µL)
- 5 From 20 µM stock solution of Amplex Red reagent - Prepare 400 µM dilution containing 2 U/mL HRP by adding 20 µL of Amplex Red stock solution and 100 µL of 200 U/mL HRP stock solution to 9.7 mL of 1X Reaction Buffer.
- 6 Initiate reaction by adding 100 µL from above to each well.
- 7 Place plate into Fluroskan holder and click **START**.
- 8 Take 4 readings at 15 minute intervals using 544/590 filter pairs.
- 9 Save raw data as an Excel file into the RHPx data folder. Use the naming convention RHXXXX.xls, where XXXX is the date in mmdd format.
- 10 Select Process>Organize. Choose the appropriate data to organize (usually Measure1), then click **OK**. This rearranges the data into columns.



- 11 Save organized data as an Excel file into the RHP data folder. Use the naming convention rhXXXXor.xls, where XXXX is the date in mmdd format.