

Oct 27, 2021 Version 1

Biofilm growth with starch treatment V.1

DOI

dx.doi.org/10.17504/protocols.io.eq2lypqzelx9/v1

Bjorn Bartholdy¹, a.g.henry¹

¹Leiden University

BYOC



Bjørn Peare Bartholdy

Delft University of Technology, Leiden University

OPEN  ACCESS



DOI: dx.doi.org/10.17504/protocols.io.eq2lypqzelx9/v1

Protocol Citation: Bjorn Bartholdy, a.g.henry 2021. Biofilm growth with starch treatment. **protocols.io**
<https://dx.doi.org/10.17504/protocols.io.eq2lypqzelx9/v1>

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: In development

We are still developing and optimizing this protocol

Created: May 20, 2021

Last Modified: October 27, 2021

Protocol Integer ID: 50123

Funders Acknowledgements:

European Union's Horizon 2020 research and innovation program
Grant ID: STG-677576 ("HARVEST")



Abstract

Protocol to grow mature, mineralised oral biofilm with starches.

Modified from protocol by Sissons et al.:

Sissons, C. H., Cutress, T. W., Hoffman, M. P., & Wakefield, J. S. J. (1991). A Multi-station Dental Plaque Microcosm (Artificial Mouth) for the Study of Plaque Growth, Metabolism, pH, and Mineralization: *Journal of Dental Research*. <https://doi.org/10.1177/00220345910700110301>

Guidelines

The preferred substratum for inoculation is glass or hydroxyapatite. Plastic substrata can be used but are less effective, so surface treatment of the plastic is recommended (e.g. heated HCl or acetone).

Substrata should be autoclaved prior to use, both before and during the experiment.

Materials

Solutions

Artificial saliva

Protocol



NAME

Artificial saliva

CREATED BY

Bjørn Peare Bartholdy

PREVIEW

CPMU

Protocol



NAME

CPMU

CREATED BY

Bjørn Peare Bartholdy

PREVIEW

20% (v/v) sterile glycerol in dH₂O

5% (w/v) sucrose in dH₂O

0.25% (w/v) potato starch in dH₂O

0.25% (w/v) wheat starch in dH₂O

0.50% (w/v) equal parts wheat (0.25%) + potato (0.25%) in dH₂O


Equipment


24 deepwell polypropylene plates (w. lid containing pegs suspended from the lid)


Shaking incubator



Protocol materials

 Sucrose

 Sucrose

 Sucrose

Before start

Saliva donor criteria

- Must have no/limited history of dental caries
- Must not have used antibiotics in the past 6 months
- Abstain from oral hygiene 24 hours prior to donation
- Refrain from eating and drinking (except water) 2 hours before donation

For experiments involving starches, donors should avoid eating starch-containing foods on the day of donation. To make this more bearable, saliva donation should take place in the morning before breakfast.



Saliva collection

- 1 Saliva donors rinse their mouth with water for 30 seconds.
- 2 Stimulate saliva production by chewing tasteless gum or parafilm.
- 3 Collect the saliva by spitting into 50 ml plastic centrifuge tubes.

Note

Make sure donors wear gloves to avoid contamination with non-oral bacteria.

- 4 Make a 2-fold dilution of saliva in **sterile** [M] 20 % (v/v) glycerol and vortex the solution.

Day 0: Inoculation and feeding

- 5 Before inoculation, vortex the saliva solution again.
- 6 Pipette the saliva solution into the wells, so approx. 1-2 cm of the substratum is submerged.
- 7 Place the plate in the incubator at 36°C for 4 hours for static inoculation ⌚ 04:00:00
- 8 After inoculation, transfer the samples to a new plate containing the artificial saliva, and place in a shaking incubator at ↻ 30 rpm, 36°C for ⌚ 04:00:00

4h

4h



Protocol






NAME

Artificial saliva


CREATED BY

Bjørn Peare Bartholdy



PREVIEW


8.1 Add  300 mL distilled (or deionized) dH₂O to a  1000 mL beaker, with stirring and heat  60 °C .

8.2 Add:



-  2.5 g

 Mucin from porcine stomach (Type III) **Becton Dickinson (BD) Catalog #M1778**

-  5 g  Trypticase™ Peptone **Becton Dickinson (BD) Catalog #211921**

-  10 g

 Oxoid™ Proteose Peptone **Becton Dickinson (BD) Catalog #LP0085B**


-  5 g  Bacto Yeast Extract **Becton Dickinson (BD)**

Let the reagents completely dissolve before continuing to the next step

8.3 Add:

-  2.5 g  KCl **Becton Dickinson (BD)**





-  0.35 g  NaCl **Becton Dickinson (BD)**

-  0.2 g  CaCl₂ **Becton Dickinson (BD)**


-  0.74 g

 Sodium phosphate dibasic **Becton Dickinson (BD) Catalog #7558-79-4**



-  0.54 g  NaHCO₃ **Becton Dickinson (BD)**
-  2.5 mg  Hemin **Becton Dickinson (BD)**

8.4 Add the remaining  700 mL distilled (or deionized) dH₂O

8.5 Adjust to  7 with  NaOH **Becton Dickinson (BD)** and stirring

8.6 Transfer to two 1000 ml bottles, so half of each bottle is filled.

15m

Autoclave at  121 °C ,  1 Bar for  00:15:00 minutes

Safety information


Do NOT screw bottle caps on tightly.

Loosely screw the caps on the bottles or cover the tops with foil

8.7 Once the solution has cooled, add:



-  1 mg  Menadione **Becton Dickinson (BD)**
-  0.3 g  Urea **Becton Dickinson (BD)**
-  0.17 g  L-Arginine **Becton Dickinson (BD)** Catalog #A5006

8.8 Store in fridge at ca.  4 °C

Occasionally test the pH to ensure it stays around  7

9 Transfer the samples to a plate containing a 5% (m/v)  Sucrose solution for







6m

 00:06:00 , then transfer back to the artificial saliva and leave  Overnight .

Day 1-2: Feeding

8h 12m



- 10 First thing in the morning, transfer the samples to a new plate containing a 5% (m/v)  Sucrose solution for  00:06:00 . While in the sucrose solution, add more artificial saliva to the wells on the original plate that have been partially depleted overnight. 6m
- 11 After the 6 mins. return the samples to the plate with **artificial saliva**, and cover up the sucrose plate and leave for  08:00:00 . 8h
- 12 After 8 hours, transfer the samples back to the plate with 5% (m/v)  Sucrose **Contributed by users** solution for  00:06:00 . Transfer back to the artificial saliva and leave  Overnight . Dispose of the sucrose. 6m

Day 3: Inoculation and feeding

8h

- 13 Repeat steps from saliva collection and Day 0: Inoculation and feeding. 8h

Expected result

A layer of clear plaque should be visible on the substrata

- 14 Prepare a new plate with artificial saliva. Transfer the samples from the inoculation plate to the **artificial saliva**.

Day 4: Feeding

8h

- 15 Repeat steps 10 through 12  [go to step #10](#) 8h

Note

Prepare a new plate of artificial saliva every third day throughout the experiment. Every other morning, top up the wells with artificial saliva (so ca. 1-2 mm of the substratum is submerged).



Day 5: Inoculation

16 Repeat steps 1 through 9

⇒ go to step #1

Day 6-8: Feeding

17 Repeat steps 10 through 12

⇒ go to step #10

Day 9-14: Starch treatment

8h 6m

18 Transfer the samples to a plate with the starch treatment(s) for ⌚ 00:06:00 at

6m

⌚ 60 rpm, 36°C

Safety information

If you are using multiple starch treatments within a single plate, take care to avoid cross-contamination.

19 Transfer the samples back to the artificial saliva plate for ⌚ 08:00:00 at

8h

⌚ 30 rpm, 36°C

20 Transfer the samples to a plate with the starch treatment(s) for ⌚ 00:06:00 at

6m

⌚ 60 rpm, 36°C

21 Transfer the samples back to the artificial saliva plate at ⌚ 30 rpm, 36°C and leave

⌚ Overnight

Day 15-24: Mineralisation

8h 30m

22 Transfer the samples to a plate containing the calcium phosphate monofluorophosphate urea (CPMU) solution for ⌚ 00:06:00 at ⌚ 60 rpm, 36°C

6m



Protocol



NAME

CPMU

CREATED BY

Bjørn Peare Bartholdy

PREVIEW

22.1 Add 300 mL distilled (or deionized) dH₂O to a 1000 mL beaker, with stirring and heat 60 °C

22.2 Add:

- 1.55 g Calcium Chloride **Sigma Aldrich**
- 1.44 g Sodium Phosphate monobasic **Sigma Aldrich**
- 0.72 g Sodium Fluorophosphate **Sigma Aldrich Catalog #344443**
- 0.08 g Magnesium Chloride **Sigma Aldrich Catalog #AC223210010**
- 30 g Urea **Sigma Aldrich**

22.3 Add the remaining 700 mL and keep stirring until precipitate has completely dissolved
Store in fridge at 4 °C

23 Transfer the samples back to the **artificial saliva** for 02:00:00 at 30 rpm, 36°C

2h

Put a lid on the plate with CPMU (or cover with foil) to prevent evaporation.


Repeat step 22 and 23, four more times every two hours.

24 Transfer the samples to a plate with the **starch treatment(s)** for 00:06:00 at 60 rpm, 36°C

6m



25 Transfer the samples back to the **artificial saliva plate** at  30 rpm, 36°C and leave

 Overnight

Analysis

26 Samples should be dried before sampling.

Transfer to a new plate with no liquid and leave in the incubator.

27 Once dried, samples are processed like archaeological samples.