

Apr 10, 2020 Version 1

# BioMark Single Cell Protocol (Two-Step RTSTA) V.1

DOI

[dx.doi.org/10.17504/protocols.io.wdsfa6e](https://dx.doi.org/10.17504/protocols.io.wdsfa6e)



Shaina Robbins<sup>1</sup>, Alison Moss<sup>1</sup>, Sean Nieves<sup>1</sup>

<sup>1</sup>Thomas Jefferson University

SPARC

Tech. support email: [info@neuinfo.org](mailto:info@neuinfo.org)



Sirisha Achanta

OPEN  ACCESS



**DOI:** [dx.doi.org/10.17504/protocols.io.wdsfa6e](https://dx.doi.org/10.17504/protocols.io.wdsfa6e)

**Protocol Citation:** Shaina Robbins, Alison Moss, Sean Nieves 2020. BioMark Single Cell Protocol (Two-Step RTSTA). [protocols.io https://dx.doi.org/10.17504/protocols.io.wdsfa6e](https://dx.doi.org/10.17504/protocols.io.wdsfa6e)

**License:** This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

**Protocol status:** Working

**We use this protocol and it's working**

**Created:** December 12, 2018

**Last Modified:** April 30, 2020

**Protocol Integer ID:** 18578

















**Keywords:** Biomark Fluidigm High Throughput Real Time PCR

## Abstract

This protocol utilizes Fluidigm's Biomark system which performs high-throughput real-time PCR that can assay 48 or 96 genes for 48 or 96 samples respectfully. This protocol is used for gene expression targeting samples at the single-cell scale and can be used with the 48.48 Dynamic Array integrated fluidic circuit (IFC) or the 96.96 Dynamic Array IFC. Note that this protocol assumes that single-cell samples are captured using the Arcturus Laser Capture Microdissection system.




## Materials

### MATERIALS

-  T4 Gene 32 Protein - 100 ug **New England Biolabs Catalog #M0300S**
-  Exonuclease I (E.coli) - 3,000 units **New England Biolabs Catalog #M0293S**
-  Thermocycler
-  TaqMan™ PreAmp Master Mix **Applied Biosystems (ThermoFisher Scientific) Catalog #4391128**
-  TaqMan Universal PCR Master Mix **Life Technologies Catalog #4304437**
-  TE Buffer **Invitrogen - Thermo Fisher Catalog #AM9861**
-  CapSure HS Caps **Applied Biosystems (ThermoFisher Scientific) Catalog #LCM0214**
-  CellsDirect™ One-Step qRT-PCR Kit **Invitrogen - Thermo Fisher Catalog #11753100**
-  SuperScript™ VILO™ cDNA Synthesis Kit **Invitrogen - Thermo Fisher Catalog #11754050**
-  2X SsoFast EvaGreen Supermix with Low ROX **Bio-Rad Laboratories Catalog #1725211**
-  20x DNA Binding Dye Sample Loading Reagent **Fluidigm Catalog #100-7609**
-  2X Assay Loading Reagent **Fluidigm Catalog #100-7611**
-  Control Line Fluid Kit—96.96 **Fluidigm Catalog #89000021**
-  Control Line Fluid Kit—48.48 **Fluidigm Catalog #89000020**
-  48.48 Dynamic Array™ IFC for Gene Expression **Fluidigm Catalog #BMK-M-48.48**
-  96.96 Dynamic Array™ IFC for Gene Expression **Fluidigm Catalog #BMK-M-96.96**
-  IFC Controller HX **Fluidigm**
-  BioMark HD System **Fluidigm**

## LCM Extraction


1 Prepare lysis buffer as follows:

-  0.5  $\mu\text{L}$  Lysis enhancer (CellDirect kit)
-  5  $\mu\text{L}$  Resuspension buffer (CellDirect kit)
-  5.5  $\mu\text{L}$  Total lysis buffer per sample

2 Capture single cells on CapSure HS caps.

Add  5.5  $\mu\text{L}$  lysis buffer per cap

Incubate at  75 °C on heat block for  00:15:00 cap-side down

Incubate on ice for  00:05:00 cap-side down

Incubate on ice for  00:05:00 right-side up

Spin down and transfer lysate to PCR tube or plate or freeze at  -80 °C

### Note

Sample size will vary, i.e., single cell, 10-cell pool, 100-cell pool, etc.

## RNA Dilution Series (Optional)

3 Prepare total RNA dilutions of 1 ng/ $\mu\text{L}$ , 300 pg/ $\mu\text{L}$ , 100 pg/ $\mu\text{L}$ , 30 pg/ $\mu\text{L}$ , 10 pg/ $\mu\text{L}$ , 3 pg/ $\mu\text{L}$ , 1 pg/ $\mu\text{L}$ , 300 fg/ $\mu\text{L}$  and H<sub>2</sub>O.

Add  1  $\mu\text{L}$  of RNA to  5.5  $\mu\text{L}$  lysis buffer.


Treat RNA dilution samples the same as the LCM samples.

**Note**

We assume the average RNA yield from single cell LCM is approximately 10 pg.

## Reverse Transcription

- 4 Transfer all samples to a PCR plate.







Add  1.2  $\mu\text{L}$  5X VILO Reaction Mix to each cell lysate (including dilution series).

**Note**

VILO Reaction Mix is part of the VILO cDNA synthesis kit and always runs short. Make sure sufficient reagent is available before RT.

- 5 Incubate plate at  65 °C for  00:01:30 , snap cool on ice and spin down.

- 6 Mix the following reagents and add  1  $\mu\text{L}$  to each sample (for 96.96 Biomark chip):

- 10X Superscript III Mix (VILO kit)  0.15  $\mu\text{L}$  x 110 =  16.5  $\mu\text{L}$
- T4 Gene 32 Protein  0.12  $\mu\text{L}$  x 110 =  13.2  $\mu\text{L}$
- H<sub>2</sub>O  0.73  $\mu\text{L}$  x 110 =  80.3  $\mu\text{L}$

Final sample volume including lysate =  7.5  $\mu\text{L}$



## 7 Incubate PCR plate using the following thermocycles:

- 25 °C for 00:05:00
- 50 °C for 00:30:00
- 55 °C for 00:25:00
- 60 °C for 00:05:00
- 70 °C for 00:10:00
- Hold at 4 °C

## Preparing Pre-Amp Primer Pool

### 8 For 96.96 chip, each primer pair is diluted to 50 uM

- Primer pairs (50 uM) 2 µL x 96 = 192 µL
- H<sub>2</sub>O 8 µL
- Total = 200 µL of primer pool (500 nM concentration)

## PreAmp

### 9 Mix together and add 10.5 µL to each RT product (for 96.96 chip):

- Taqman PreAmp Master Mix 9 µL x 110 = 990 µL
- 500 nM primer pool 1.5 µL x 110 = 165 µL

Final sample volume = 18 µL

## 10 Incubate PCR plate using the following thermocycles:

- 95 °C for 00:10:00

22 cycles of:

- 96 °C for 00:00:05
- 60 °C for 00:04:00

### Note

We consider 22 cycles of pre-amp optimal for LCM captured neurons. Number of cycles can be adjusted according to sample type.

## 11 Exonuclease Treatment

Mix the following and add 6 µL to each sample (for 96.96 chip):

- Exonuclease I reaction buffer 10X 0.6 µL x 110 = 66 µL
- Exonuclease I 1.2 µL x 110 = 132 µL
- H2O 4.2 µL x 110 = 462 µL


Final sample volume = 24 µL

## 12 Incubate PCR plate using the following thermocycles:

- 37 °C for 00:30:00
- 80 °C for 00:15:00



Add 51 µL of TE buffer to each reaction.

### Note



Water can be used if long term storage of amplified cDNA is not required. Samples can be stored in  -20 °C if they will not be taken to the next step right away.


## Prepare Sample Pre-Mix and Samples

13 Prepare the sample pre-mix:

- 2X SsoFastEvaGreen Supermix with low ROX  4 µL
- 20X DNA Binding Dye Sample Loading Reagent (green cap)  0.4 µL

Add  4.4 µL sample pre-mix to  3.6 µL of sample **after PreAmp** for a total of  8 µL per sample.

Vortex the Sample Mix solution for a minimum of  00:00:20 , and centrifuge for at least  00:00:30 .

Prepared reactions can be stored for short times at  4 °C until the samples are ready to be loaded into the chip.

## Prepare Assay Mix



14



### Note

The same preparation of primers can be used for the Fluidigm 48.48 Dynamic Array IFC and the Fluidigm 96.96 Dynamic Array IFC.



Combine the following:

- 2X Assay Loading Reagent (yellow cap)  4  $\mu$ L
- 20 uM primers  4  $\mu$ L

Vortex the Assay Mix for a minimum of  00:00:20 and centrifuge for at least  00:00:30 to spin down all the components.



#### Note

Vortex thoroughly and centrifuge all samples and assay solutions before pipetting into the chip inlets. Failure to do so may result in a decrease in data quality. The final concentration of each primer is 5 uM in the inlet, and 500 nM in the reaction chamber.

## Priming and Loading the Dynamic Array IFC

15

#### Note




**CAUTION!** Due to different accumulator volumes, use the appropriate control syringe for your chip type:  300  $\mu$ L for the 48.48 chip or  150  $\mu$ L for the 96.96 chip.

- Inject control line fluid into each accumulator on the chip (see diagram).
- Remove and discard the blue protective film from the bottom of the chip.
- Place the chip in the IFC Controller MX for the 48.48 chip or the IFC Controller HX for the 96.96 chip, then run the **Prime (113x)** script for the 48.48 chip or the **Prime (136x)** script for the 96.96 chip.
- When the prime script has finished, press **Eject** to remove the primed chip from the IFC Controller.



### Note

**CAUTION!** While pipetting, do not go past the first stop on the pipette. Doing so may introduce air bubbles into the inlets.

- Pipette  5  $\mu$ L of each assay and  5  $\mu$ L of each sample into their respective inlets on the chip (refer to diagram).
- **IMPORTANT:** You must pop/remove any air bubbles that are in the inlets. This must be done in under  00:30:00 .
- Return the the chip to the IFC Controller.
- Using the software, run the **Load Mix (113x)** script for the 48.48 chip or **Load Mix (136x)** script for the 96.96 chip to load the samples and assays into the chip.
- When the Load Mix script has finished, remove the loaded chip from the IFC Controller.
- Remove any dust particles or debris from the chip surface using scotch tape. you are now ready for the chip run.



96.96chip with well number and co...

## BioMark Run

- 16 1. Double-click the Data Collection Software Icon on the desktop to launch the software.
2. Click **Start a New Run**.
3. Check the status bar to verify the lamp and the camera are ready. Make sure both are green before proceeding.
4. Place the chip into the reader.
5. Click **Load**.
6. Verify chip barcode and chip type. Choose project settings (if applicable). Click **Next**.
7. Chip run file: Select **New**. Browse to a file location for data storage. Click **Next**.
8. Application, Reference, Probes: Select Application Type--**Gene Expression** for version 3.1.2 or higher software (for earlier versions, contact Fluidigm Technical Support). Select Passive References- **ROX**. Select Probe--**Single probe**. Select probe type- **EvaGreen**. Click **Next**.
9. Click **Browse** to find the thermal cycling protocol file. For BioMark HD- GE Fast 48×48 PCR+Metl v2.pcl or GE Fast 96×96 PCR+Melt v2.pcl. For BioMark- GE 48×48 PCR+Metl v2.pcl or GE 96×96 PCR+Melt v2.pcl.
10. Confirm **Auto Exposure** is selected.
11. Click **Next**.



12. Verify the chip run information.

13. Click **Start Run**.